

**National Institute of Mental Health  
Psychoactive Drug Screening Program  
(NIMH PDSP)**

**ASSAY PROTOCOL BOOK  
Version II**

March 2013

Bryan L. Roth, MD, PhD

Department of Pharmacology  
University of North Carolina at Chapel Hill

This is an updated version as of March 2013 of the previous PDSP assay protocol book, and contains detailed descriptions of experimental procedures and data analysis for radioligand binding and functional assays as performed by the NIMH-PDSP.

# Table of Contents

1. [Radioligand binding assays](#)
  - 1.1. [Drug plate preparation \(Hamilton ASM® and STAR®\)](#)
  - 1.2. [Cell culture and membrane fraction preparation for binding assays](#)
    - 1.2.1. [Calcium phosphate precipitation transfection](#)
    - 1.2.2. [Membrane preps from stable cells](#)
    - 1.2.3. [Membrane preps from transient cells](#)
    - 1.2.4. [Membrane preps from tissues](#)
      - 1.2.4.1. [Membrane preps for NMDA binding](#)
    - 1.2.5. List of cell lines for membrane pellets
  - 1.3. [Saturation binding assay](#)
  - 1.4. [Primary and secondary binding assay](#)
  - 1.5. [Radioligand binding assays for nicotinic acetylcholine receptors](#)
    - 1.5.1. [Cell culture and membrane preparation](#)
    - 1.5.2. [General procedures for nAChR radioligand binding assay](#)
    - 1.5.3. [Primary and Secondary nAChR radioligand binding assay](#)
    - 1.5.4. [Data analysis and result reporting](#)
  - 1.6. [General binding data entry, analysis, and quality control \(non nAChRs\)](#)
    - 1.6.1. [Built-in analysis tool in the PDSP Database](#)
      - 1.6.1.1. [Saturation binding result analysis and reporting](#)
      - 1.6.1.2. [Primary binding result analysis and reporting](#)

- 1.6.1.3. [Secondary binding result analysis and reporting](#)
- 1.6.2. [Data analysis using Prism v5.0](#)
  - 1.6.2.1. [Saturation binding result analysis](#)
  - 1.6.2.2. [Primary binding result analysis](#)
  - 1.6.2.3. [Secondary binding result analysis](#)
- 1.7. Table of conditions for binding assays and representative figures
  - 1.7.1. [5-HT receptors](#)
  - 1.7.2. [Muscarinic acetylcholine receptors](#)
  - 1.7.3. [Nicotinic acetylcholine receptors](#)
  - 1.7.4. [Adenosine receptors](#)
  - 1.7.5. [Adrenergic receptors](#)
  - 1.7.6. [Cannabinoid receptors](#)
  - 1.7.7. [Dopamine receptors](#)
  - 1.7.8. [Histamine receptors](#)
  - 1.7.9. [Neurotransmitter transporters \(DAT, NET, and SERT\)](#)
  - 1.7.10. [Opioid receptors](#)
  - 1.7.11. [Oxytocin and Vasopressin receptors](#)
  - 1.7.12. [Smoothed receptor](#)
  - 1.7.13. [Prostanoid receptors](#)
  - 1.7.14. [Protein Kinase C](#)
  - 1.7.15. [GABAA receptors](#)
  - 1.7.16. [Glutamate receptors \(NMDA and mGluRs\)](#)

- 1.7.17. [Sigma receptors](#)
- 1.7.18. [Angiotensin II receptors](#)
- 1.7.19. [Neurotensin receptors](#)
- 1.7.20. [Vesicular monoamine transporter \(VMAT2\)](#)
- 1.7.21. [Calcium and Sodium channels](#)
- 1.7.22. [Potassium channels \(hERG\)](#)
- 1.7.23. [Imidazoline receptor](#)
- 2. [Functional assays](#)
  - 2.1. [Drug plate preparations for functional assays](#)
  - 2.2. [General procedures for PDSP functional assays](#)
  - 2.3. [Data analysis for functional assays](#)
    - 2.3.1. [Agonist activity](#)
    - 2.3.2. [Antagonist activity](#)
    - 2.3.3. [Schild plot analysis](#)
    - 2.3.4. [Quantifying bias](#)
    - 2.3.5. [Allosteric operational analysis](#)
    - 2.3.6. List of available functional assays at NIMH-PDSP
  - 2.4. [Functional assays for G<sub>q</sub> coupled GPCRs:](#)
    - 2.4.1. [Calcium mobilization assays \(FLIPR<sup>TETRA</sup>\)](#)
    - 2.4.2. [Intracellular inositol phosphate accumulation](#)
    - 2.4.3. [Table of G<sub>q</sub> coupled GPCRs with functional assays](#)
    - 2.4.4. Representative figures

- 2.4.4.1. [Muscarinic receptors \(M1, M3, and M5\)](#)
- 2.4.4.2. [5-HT receptors \(5-HT2A, 5-HT2B, and 5-HT2C\)](#)
- 2.4.4.3. [Adrenergic receptors \( \$\alpha\$ 1A,  \$\alpha\$ 1B, and  \$\alpha\$ 1D\)](#)
- 2.4.4.4. [Bombesin receptors \(BB1 and BB2\)](#)
- 2.4.4.5. [Bradykinin receptor \(B2\)](#)
- 2.4.4.6. [CCK2 and Ghrelin receptors](#)
- 2.4.4.7. [Histamine receptors \(H1 and H2\)](#)
- 2.4.4.8. [mGluR receptors \(mGluR1 and mGluR5\)](#)
- 2.4.4.9. [Melanocortin receptors \(MC3 and MC4\)](#)
- 2.4.4.10. [Tachykinin receptors \(NK1, NK2, and NK3\)](#)
- 2.4.4.11. [Neurotensin receptor \(NTS1\)](#)
- 2.4.4.12. [Purinergic receptors \(P2Y1, P2Y2, P2Y4, P2Y6, and P2Y11\)](#)
- 2.4.4.13. [Protease-activated receptor \(PAR1\)](#)
- 2.4.4.14. [Platelet-activating factor receptor \(PTAFR\)](#)
- 2.4.4.15. [Oxytocin and Vasopressin receptors \(V1A, V1B, V2, and OT\)](#)
- 2.4.4.16. [GPR40, GPR41, GPR43](#)
- 2.4.4.17. [Representative figures for Schild analysis](#)
- 2.4.4.18. [Representative figures and tables for quantifying bias](#)
- 2.4.4.19. [Representative figures for allosteric operational analysis](#)
- 2.5. [Functional assays for  \$G\_i\$  or  \$G\_s\$  coupled GPCRs: split luciferase based biosensor assay to measure cAMP production \(GloSensor® cAMP assay\)](#)
  - 2.5.1. [Cell culture and transfection](#)

- 2.5.2. [Split luciferase based cAMP assay – Luciferin first protocol](#)
- 2.5.3. [Split luciferase based cAMP assay – Drug first protocol](#)
- 2.5.4. [Split luciferase based cAMP assay – G<sub>i</sub> coupled mGluRs](#)
- 2.5.5. [Table of G<sub>i</sub> or G<sub>s</sub> coupled GPCRs with functional assays](#)
- 2.5.6. Representative figures
  - 2.5.6.1. [Muscarinic receptors \(M2, M4, and DREADD\)](#)
  - 2.5.6.2. [5-HT receptors \(5-HT1, 5-HT4, 5-HT6, 5-HT7A\)](#)
  - 2.5.6.3. [Adrenergic receptor \( \$\beta\$ 2\)](#)
  - 2.5.6.4. [Adenosine receptor \(A2A\)](#)
  - 2.5.6.5. [Corticotropin-releasing factor receptors \(CRF1, CRF2\)](#)
  - 2.5.6.6. [Dopamine receptors \(D1, D2, D4, and D5\)](#)
  - 2.5.6.7. [Histamine receptors \(H2 and H3\)](#)
  - 2.5.6.8. [Melanocortin receptors \(MC1 and MC3\)](#)
  - 2.5.6.9. [Opioid receptors \(DOR, MOR, KOR, and NOP\)](#)
  - 2.5.6.10. [mGluRs \(mGluR2, mGluR3, mGluR4, mGluR6, mGluR8\)](#)
  - 2.5.6.11. [Somatostatin receptor \(SSTR5\)](#)
  - 2.5.6.12. [GPR68 and GPR88](#)
- 2.6. [Functional assays for G-protein independent  \$\beta\$ -arrestin recruitment](#)
  - 2.6.1.  [\$\beta\$ -arrestin recruitment \(GPCR Tango assay\)](#)
  - 2.6.2. [Tables of GPCR Tango constructs](#)
  - 2.6.3. [Representative figures](#)
- 2.7. [Parallel receptorome screening using  \$\beta\$ -arrestin recruitment assay](#)

- 2.7.1. [Cell culture](#)
- 2.7.2. [DNA plate preparation](#)
- 2.7.3. [Transfection using Calcium phosphate precipitation method](#)
- 2.7.4. [Assay procedure](#)
- 2.7.5. [Data processing and analysis](#)
- 2.7.6. [Representative figures](#)
- 2.8. [hERG functional assay: Thallium \(Tl<sup>+</sup>\) flux assay \(FluxOR\)](#)
  - 2.8.1. [Cell culture](#)
  - 2.8.2. [Assay procedure for hERG modulators](#)
  - 2.8.3. [Assay procedure for hERG trafficking modulators](#)
  - 2.8.4. [Data analysis and representative figures](#)
- 2.9. [hERG electrophysiology: automated parallel patch-clamp assay](#)
  - 2.9.1. [External and internal solutions](#)
  - 2.9.2. [Cell preparation](#)
  - 2.9.3. [Drug plate preparation](#)
  - 2.9.4. [Assay procedure](#)
  - 2.9.5. [Data analysis and representative figures](#)
- 2.10. [Neurotransmitter transporter functional assay](#)
- 2.11. [Multidrug resistance transporter -1 \(MDR-1\) modulation assay](#)
- 2.12. [Histone Deacetylase \(HDAC\) assay](#)
- 2.13. [Monoamine Oxidase \(MAO\) activity assay](#)
- 2.14. [Protein Kinase C \(PKC\) activity assay](#)

- 2.15. [Checkpoint Kinase 2 \(CHK2\) activity assay](#)
- 2.16. [Nicotinic acetylcholine receptor functional assay –  \$^{86}\text{Rb}^+\$  efflux](#)
  - 2.16.1. [Cell culture](#)
  - 2.16.2. [General assay procedure](#)
  - 2.16.3. [Primary functional assays](#)
  - 2.16.4. [Secondary functional assays](#)
  - 2.16.5. [Data analysis and representative figures](#)
- 3. [Master table of targets at NIMH-PDSP](#)
- 4. [List of publications supported by NIMH-PDSP](#)

## **Section 1: Radioligand binding assays**

### **1.1. Drug plate preparation for radioligand binding assays: Hamilton ASM® and STAR®**

The NIMH PDSP uses an **Integrated Sample Storage, Retrieval & Liquid Handling Robotic System** (see SOW) to prepare assay plates for the research staff to perform experiments. Briefly, the system consists of a refrigerated storage module that stores the PDSP samples and standard compounds; a retrieval module which selects, retrieves, thaws, uncaps and delivers the compounds; and a multi-station liquid handling platform to transfer the samples and standards to make assay plates. The system is controlled by four PC computers utilizing Hamilton's Microlab STAR Venus One and Sample Manager software suites. Considerable training is required to operate the system. Advanced training is required to program and trouble-shoot the system.

A set of steps to run the machines to create a plate is called a "method." Each method has hundreds of lines of instructions. Below is an overview of the steps required to create assay plates for primary and secondary assay plates. This is an overview of the liquid handling and major plate movement steps. Robotic movements for tip set-up, tip movements, tip recycling, plate staging, etc. are not included.

#### **A. Primary Plate Preparation:**

1. A "work list" is created by a senior staff member utilizing our PDSP Database.
2. The work list is downloaded onto the ASM computer.
3. The "run control" program is opened on the ASM server computer.
4. The assay method corresponding to the assay to be run is selected; the ASM server and store, de-capper and robotic arms perform self-checks.
5. The user is prompted to select the work list to be run.
6. The run control program verifies the work list, prompts the user to select "start" and thaw time variables or defaults.
7. Sample selection begins in the ASM Store, which sends messages to the appropriate liquid handling station.

- A. The ASM Store selects and retrieves samples using the 2-dimensional bar code on each tube.
  - B. Samples are collected in an empty rack, inside the ASM Store
  - C. When all samples are collected the rack is sent to the ASM server via an automated trolley and the server's internal single-grip arm.
  - D. Samples are thawed with forced hot air according to protocol with an option to enter non-default thaw time.
  - E. The rack is then delivered via the ASM Server single-grip arm and hand-off arm to the Capper/Decapper
  - F. Sample tubes are uncapped and sent to the appropriate STAR deck
8. The user prepares the STAR liquid handler deck with proper plates, buffer(s) and pipetting tips. (Steps 8-14 are performed in parallel to step 7.)
  9. The STAR computer run control (this is a separate instance from the ASM run control) prompts the user for tip-counting verification; the user also affirms that the protocol and plate count are correct.
  10. The STAR platform begins the plating procedure.
  11. Assay plates are distributed to the STAR deck (1 per final plate); the plate barcode is recorded and used throughout our process to verify the integrity of sample tracking.
  12. 97.5 microliters of buffer is added to each of the wells in rows A and E of each assay plate.
  13. DMSO (2.5 microliters, vehicle control) is added to A1 & E1 of each assay plate.
  14. The STAR platform pauses with user input, waits for drug delivery from the ASM (step 7E above).
  15. Samples are mixed after step 7F is completed, and the sample rack is delivered to the STAR deck.
  16. 2.5 microliters of control compound is aspirated from the sample tube and distributed into the assay plate in wells A12 and E12 using the eight single-channel pipettors.
  17. 2.5 microliters of each sample is aspirated from their sample tubes and distributed into the assay plate in wells A2 through A11 and E2 through E11 using the eight single-channel pipettors.

18. Row A is then mixed and aliquoted from A to B and A to C and A to D with a single-row selection of pipette tips by the 96-well head, each well then has 25 microliters of combined sample and buffer.
19. Step 18 is repeated for row E-H to make replicates.
20. A plate map is created corresponding to the plate bar-code, samples are added and the receptor to be assayed is assigned.
21. Steps 16-20 are repeated for each assay plate. (See Figure 1 for primary drug plate map).


**B. Secondary Plate Preparation:**

1. A “work list” is created by a senior staff member utilizing our PDSP Database.
2. The work list is down-loaded onto the ASM computer.
3. The “run control” program is opened by the user on the ASM server computer.
4. The assay “method” corresponding to the assay to be run is selected; the ASM server & store, de-capper and robotic arms perform self-checks.
5. The user is prompted to select the work list to be run.
6. The “run control” program verifies the work list, prompts the user to select “start” and the thaw time variables or defaults.
7. Sample selection and retrieval begin in the ASM Store, which sends messages to the appropriate liquid handling station.
  - A. The ASM Store selects samples using the 2-dimensional bar code on each tube
  - B. Samples are collected in an empty rack, inside the ASM Store
  - C. When all samples are collected the rack is sent to the ASM server via an automated trolley and the server’s internal single-grip arm.
  - D. Samples are thawed with forced hot air according to protocol with an option to enter non-default thaw times.
  - E. The rack is then delivered via the ASM Server single-grip arm and hand-off arm to the Capper/Decapper
  - F. Sample tubes are uncapped and sent to the appropriate STAR deck
8. The user prepares the STAR liquid handler deck with proper plates, buffer(s) and pipetting tips. (Steps 8-14 are performed in parallel to step 7.)

9. The STAR computer run control program (this is a separate instance from the ASM run control program) prompts the user for tip-counting verification; the user also affirms that the protocol and plate count are correct.
10. The STAR platform begins the plating procedure.
11. Assay plates are distributed to the STAR deck (1 dilution plate and 3 plates for replicates per set); the plate barcode is recorded and used throughout our process to verify the integrity of sample tracking.
12. 180 microliters of buffer is added to all wells of each assay plate with the 96CORE Head.
13. 15 microliters of buffer is added to column 12, 40 microliters is aspirated from column 11.
14. The STAR platform pauses and waits for drug delivery from the ASM (step 7F above).
15. Samples are mixed.
16. 5 microliters of sample (or control) is aspirated from the appropriate sample tube and distributed into the first dilution plate in column 12 using the eight single-channel pipettors.
17. 40 microliters of buffer is removed from column 11 of the dilution plate.
18. After mixing, 60 microliters is aspirated from column 12 and dispensed into column 11 of the dilution plate.
19. A serial dilution is performed; 20 microliters is aspirated from/dispensed to column 12 to 10, column 10 to 8, column 8 to 6, column 6 to 4, column 4 to 2, column 11 to 9, column 9 to 7, column 7 to 5 and, column 5 to 3, resulting in a half-log dilution series across the plate.
20. The entire plate is mixed and then three replicates containing 25 microliters of drug solution per well are made.
21. Steps 16-20 are repeated for each additional set of plates.
22. Plate maps are created corresponding to the plate barcode, the samples added and the receptor assigned for all plates on the deck. (See Figure 2 for secondary drug plate map).

Managing the workflow of the automated system is extremely important. First, the overall flow, that is, what assays to select, when to make assay plates and, assigning assays to a research technician must be performed by a senior staff member. Second, the workflow for the system

also requires careful planning and scheduling. Operating the system can be performed by a research technician. The operation of the system *requires* careful attention to detail and significant planning in order to achieve optimal performance. Preparing the liquid handling deck for a method, being attentive to tip usage, having all the reagents and consumables ready prior to start, as well as the sequence of primary versus secondary methods to minimize transition time between runs, are all required have the system operate for >35 hours per week (its current utilization).



Input Barcodes:

- BAGS33

### Primary Binding Plates:


Barcode: BAGS33

Receptor: H3-0 Primary BAGS33 03


Date: 3-14-2013 09-13-18

|   | 1     | 2     | 3     | 4     | 5     | 6     | 7     | 8     | 9     | 10    | 11    | 12        |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-----------|
| A | Empty | 27054 | 22578 | 22580 | 22581 | 22582 | 22583 | 22584 | 22585 | 26656 | 26657 | Histamine |
| B | Empty | 27054 | 22578 | 22580 | 22581 | 22582 | 22583 | 22584 | 22585 | 26656 | 26657 | Histamine |
| C | Empty | 27054 | 22578 | 22580 | 22581 | 22582 | 22583 | 22584 | 22585 | 26656 | 26657 | Histamine |
| D | Empty | 27054 | 22578 | 22580 | 22581 | 22582 | 22583 | 22584 | 22585 | 26656 | 26657 | Histamine |
| E | Empty | 26658 | 26772 | 26787 | 26788 | 26789 | 26790 | 26791 | 26793 | 26794 | 26796 | Histamine |
| F | Empty | 26658 | 26772 | 26787 | 26788 | 26789 | 26790 | 26791 | 26793 | 26794 | 26796 | Histamine |
| G | Empty | 26658 | 26772 | 26787 | 26788 | 26789 | 26790 | 26791 | 26793 | 26794 | 26796 | Histamine |
| H | Empty | 26658 | 26772 | 26787 | 26788 | 26789 | 26790 | 26791 | 26793 | 26794 | 26796 | Histamine |

**Figure 1.** 96-well drug plate setup diagram (upper panel) and an actual barcoded 96-well drug plate (lower panel) for primary radioligand binding assays. The drug plate is designed to contain 20 compounds (compounds “a” to “t” in Columns 2 - 11, each in quadruplicate) and one reference compound (Column 12), 25  $\mu$ l per well at 5x of the final concentration of 10  $\mu$ M. Column 1 contains buffer only, and therefore is designated as “total binding” (0% inhibition). Column 12 contains a reference compound (Histamine in this case) and is designated as nonspecific binding (100% inhibition).



**Figure 2a.** 96-well drug plate map for secondary radioligand binding assays. The drug plate is designed to contain 7 PDSP compounds (Rows A to G) and 1 reference compound (Row H) as indicated at right of the plate, one compound per row and 25  $\mu$ l per well at 5x of the final concentrations as indicated above the plate. Well H12 is designated for nonspecific binding and column 1 is designated for total binding for the plate. Each secondary assay contains 3 sets of identical 96-well plates, therefore all compounds are tested in triplicate.



**Figure 2b.** 96-well drug plate map for secondary radioligand binding assays. The drug plate is designed to contain 2 PDSP compounds (Rows C to H) and 1 reference compound (Rows A and B) as indicated at right of the plate, therefore the PDSP compounds are assayed in triplicate and the reference compound in duplicate. Each well contains 25  $\mu$ l at 5x of final concentrations as indicated above the plate.

## 1.2. Cell culture and membrane fraction preparation

The majority of assays performed by the NIMH PDSP staff employ transfected (primarily stably transfected) cell lines expressing mainly human recombinant receptors, monoamine transporters, or ion channels. A detailed description of the cell lines, the culture and sub-culture conditions, and the appropriate media are listed in [Table 1](#).

### 1.2.1. Calcium phosphate precipitation transfection.

**Reference:** Jordan M, Schallhorn A, Wurm FM (1996). Transfecting mammalian cells: optimization of critical parameters affecting calcium-phosphate precipitate formation. *Nucleic Acids Res* **24**(4): 596-601.

**2x HBS:** 140 mM NaCl, 1.5 mM Na<sub>2</sub>HPO<sub>4</sub>, 50 mM HEPES, pH 7.05, room temperature

The Calcium phosphate precipitation method was adapted from Jordan et al (1996). In brief, HEK293 T cells are subcultured in 15-cm dishes at a density of 8 – 10 million cells per dish and incubated overnight. For each 15 cm dish, 18 µg DNA, 100 µl of 2.5 M CaCl<sub>2</sub>, 100 µl of TE (1 mM Tris HCl, 0.1 mM EDTA, pH 7.60) are diluted into a final volume of 1 ml water. The mixture is then quickly added to an equal volume of 2x HBS and added dropwise to cells. Cells are then incubated overnight and for another day or two in fresh growth medium. For transient transfections with 5-HT receptors, we use medium containing dialyzed fetal bovine serum (FBS) overnight before harvesting. Cells are usually harvested 48 – 72 hours after transfection.

**1.2.2. Membrane preparations from stably transfected cells.** To make membrane fractions from stably transfected cell lines, cells are subcultured in 15-cm dishes and grown to 90% confluency. For 5-HT receptors, cells are incubated overnight in either serum-free medium or medium containing 1% dialyzed FBS before harvesting. The next day, the cells are rinsed with PBS, scraped off into 50 ml conical tubes, and pelleted by centrifugation (1000 x g, 10 min at 4°C). The cell pellet is resuspended in chilled (4°C) lysis buffer (50 mM Tris HCl buffer, pH 7.4) and triturated gently for hypotonic lysis. The suspension is then centrifuged at 21,000 x g for 20 min at 4°C to obtain a crude membrane fraction pellet. The fresh membrane pellet is then resuspended in 3 volumes of cold lysis buffer, and is immediately subjected to the Bradford

protein assay to determine protein concentration, followed by a saturation binding assay (see following section for detail) to determine receptor expression level and the affinity of a selected radioligand. Based on the receptor expression level and the  $K_d$  value, the fresh membrane suspension is stored at  $-80^\circ\text{C}$  freezer in small aliquots, such that one aliquot is sufficient for one 96-well plate to have at least 500 cpm/well when assayed at  $0.5 - 1.0 \times K_d$  value of the appropriate radioligand.

**1.2.3. Membrane preparations from transiently transfected cells.** To make membrane fractions from transiently transfected cells (usually HEK T cells), cells are transfected as described in **Section 1.2.1**. Transfected cells are usually harvested 48 – 72 hours after transfection, and processed as in **Section 1.2.2** to make membrane pellets.

**1.2.4. Membrane preparations from tissues.** To make membrane fractions from tissues, crude membrane fractions are prepared from rodent (typically rat or guinea pig) brain or kidney (purchased from PelFreeze Biologicals). The following is the general procedure used for membrane preparations unless otherwise stated. Frozen tissue (maintained at  $-80^\circ\text{C}$ ) is thawed on ice, homogenized on ice in 10 volumes of cold lysis buffer (50 mM Tris HCl, pH 7.4, containing protease inhibitor cocktail from Roche) using a Polytron homogenizer (6 pulses and 10 seconds per pulse). The homogenate is centrifuged at  $1,000 \times g$  for 10 min at  $4^\circ\text{C}$  to obtain supernatant. The supernatant is then centrifuged at  $40,000 \times g$  for 20 min at  $4^\circ\text{C}$ , and the resulting supernatant is decanted and replaced with the same lysis buffer. Two or three additional rounds of homogenization-centrifugation are performed to ensure thorough homogenization and also to wash out endogenous ligands (particularly important for GABA binding assays). The final pellet is resuspended in the same buffer and homogenized one last time. The fresh suspension is subjected to protein concentration measurement and saturation binding assay, and is then aliquoted for storage at  $-80^\circ\text{C}$  for future use. Aliquots are also made such that each aliquot is sufficient for one binding assay in one 96-well plate to have at least 500 cpm/well when assayed with a  $0.5 - 1.0 \times K_d$  level of the appropriate hot ligand.

**1.2.4.1. Membrane preparations from rat brain for NMDA radioligand binding assays.** The protocol is adapted from Chiu et al., 1999 (Alcohol 17:215-221). In brief, rat brains are thawed

on ice in SHE buffer (300 mM Sucrose, 10 mM HEPES, 2 mM EDTA, pH 7.3) and homogenized on ice in 20 ml ice-cold SHE buffer per gram of wet tissue using a glass homogenizer (at least 6 strokes). The homogenate is centrifuged at 1,000 x g for 20 min at 4°C to collect supernatant. The supernatant is then centrifuged at 16,000 x g for 1h at 4°C. The pellet is suspended in the same cold SHE buffer and stored in aliquots (10 mg per aliquot after protein concentration determination) at -80°C until use. The protein concentration is determined using the Bradford method. Immediately prior to assay, the aliquoted pellet is resuspended in 1 ml HE buffer (20 mM HEPES, 1 mM EDTA, pH 7.0) and centrifuged (16,000 x g) briefly in a benchtop centrifuge at 4°C to eliminate sucrose. The pellet is then resuspended in 0.5 ml HE buffer and incubated at 37°C for 30 min. The suspension is then centrifuged again at 13,000 x g for 10 min at 4°C. The pellet is washed four times by resuspending it in 1 ml HE buffer. The final pellet is resuspended to about 2 mg/ml concentration after protein concentration determination, and is then ready for binding assays at 100 µg per well, using the HE buffer supplemented with 100 µM Glutamate Sodium and 100 µM Glycine.

**Table 1.** List of cell lines and targets used by the PDSP to make membrane pellets for binding assays. All clones are stable lines, whereas transiently transfected cells are marked with “\*”. Clones are human unless noted. Detailed information regarding culture media is listed **Section 1.2.5.**

| Receptor                             | Note | Parental cells    | Media (see detail below the table) |
|--------------------------------------|------|-------------------|------------------------------------|
| <b>Serotonin (5HT)</b>               |      |                   |                                    |
| 5-HT1A                               |      | stable CHO        | 500 G418                           |
| 5-HT1B                               |      | stable HEK        | 500 G418                           |
| 5-HT1D                               | *    | HEKT              | COS/HEK                            |
| 5-HT1E                               |      | stable HEK        | 500 G418                           |
| 5-HT2A (rat)                         |      | stable 3T3        | 500 G418                           |
| 5-HT2A                               | *    | HEKT              | COS/HEK                            |
| 5-HT2B                               |      | stable HEK        | 2 µg/ml Puromycin                  |
| 5-HT2C                               |      | Flp-IN HEK        | DMEM 100 µg/ml Hygromycin B        |
| 5-HT3                                | *    | HEKT              | COS/HEK                            |
| 5-HT5A                               |      | Flp-In CHO        | DMEM/F-12 200 µg/ml Hygromycin B   |
| 5-HT6                                |      | stable HEK        | 500 G418                           |
| 5-HT7A                               |      | stable HEK        | 2 µg/ml Puromycin                  |
| <b>Dopamine</b>                      |      |                   |                                    |
| D1                                   | *    | HEKT              | COS/HEK                            |
| D2                                   |      | stable fibroblast | COS/HEK                            |
| D2L                                  |      | stable CHO        | F-12/10%FBS 400G418                |
| D3 (rat)                             | *    | HEKT              | COS/HEK                            |
| D3                                   | *    | HEKT              | COS/HEK                            |
| D4                                   |      | stable            | DMEM/F12 10% CS Fe+                |
| D5                                   | *    | HEKT              | COS/HEK                            |
| <b>Opioid</b>                        |      |                   |                                    |
| Mu, MOR                              |      | stable HEK        | 200 G418                           |
| Delta, DOR                           |      | stable HEK        | 200 G418                           |
| Kappa, KOR (rat)                     |      | stable HEK        | 500 G418                           |
| Kappa, KOR                           |      | stable HEK        | 500 G418                           |
| Nociceptin, NOP                      | *    | HEKT              | COS/HEK                            |
| <b>Neurotransmitter Transporters</b> |      |                   |                                    |
| SERT                                 |      | stable HEK        | 500 G418                           |
| NET                                  |      | stable HEK        | hNET (250 G418)                    |
| DAT                                  |      | stable HEK        | hDAT (350 G418)                    |
|                                      |      |                   |                                    |
|                                      |      |                   |                                    |

| Receptor                        | Note | Parental cells | Media (see detail below the table) |
|---------------------------------|------|----------------|------------------------------------|
| <b>Vasopressin and Oxytocin</b> |      |                |                                    |
| V1A                             |      | stable CHO     | V1A & OT media                     |
| V2                              |      | stable CHO     | V2 & V1B media                     |
| V1B                             |      | stable CHO     | V2 & V1B media                     |
| OT                              |      | stable CHO     | V1A & OT media                     |
| <b>Prostaglandin</b>            |      |                |                                    |
| EP-3                            | *    | HEKT           | COS/HEK                            |
| EP-4                            | *    | HEKT           | COS/HEK                            |
| <b>Adrenergic</b>               |      |                |                                    |
| alpha 1A                        |      | stable         | 500 G418                           |
| alpha 1B                        | *    | HEKT           |                                    |
| alpha 1D                        |      | stable         | 500 G418                           |
| alpha 2A                        |      | stable MDCK    | 500 G418                           |
| alpha 2B                        | *    | HEKT           | COS/HEK                            |
| alpha 2C                        |      | stable MDCK    | 500 G418                           |
| beta 1                          |      | CHO Flp-In     | DMEM/F12 200 µg/ml Hygromycin B    |
| beta 2                          |      | HEK Flp-In     | DMEM 100 µg/ml Hygromycin B        |
| beta 3                          |      | HEK Flp-In     | DMEM 100 µg/ml Hygromycin B        |
| <b>Muscarinic acetylcholine</b> |      |                |                                    |
| M1                              |      | stable CHO     | 500 G418                           |
| M2                              |      | stable CHO     | 500 G418                           |
| M3                              |      | stable CHO     | 500 G418                           |
| M3D                             |      | CHO Flp-In     | DMEM/F12 100 µg/ml Hygromycin B    |
| M4                              |      | stable CHO     | 10% FBS F12                        |
| M5                              |      | stable CHO     | 500 G418                           |
| <b>Nicotinic acetylcholine</b>  |      |                |                                    |
| α2β3                            |      | HEK            | 500 G418                           |
| α2β4                            |      | HEK            | 500 G418                           |
| α3β2                            |      | HEK            | 500 G418                           |
| α3β4                            |      | HEK            | 500 G418                           |
| α4β2                            |      | HEK            | 500 G418                           |
| α4β4                            |      | HEK            | 500 G418                           |
| α7                              |      | HEK            | 500 G418                           |
| <b>Histamine</b>                |      |                |                                    |
| H1                              |      | stable HEK     | 500 G418                           |
| H2 (in progress)                |      | stable HEK     | 500 G418                           |
| H3                              |      | HEK Flp-In     | DMEM 100 µg/ml Hygromycin B        |
| H4 (in progress)                |      |                | 500 G418                           |
|                                 |      |                |                                    |
|                                 |      |                |                                    |

| Receptor                     | Note | Parental cells | Media (see detail below the table) |
|------------------------------|------|----------------|------------------------------------|
| <b>Cannabinoid</b>           |      |                |                                    |
| CB1 (in progress)            |      | HEK            | 500 G418                           |
| CB1                          |      | HEK Flp-In     | DMEM 100 µg/ml Hygromycin B        |
| CB2                          |      | HEK Flp-In     | DMEM 100 µg/ml Hygromycin B        |
| <b>Adenosine</b>             |      |                |                                    |
| A1                           | *    | HEKT           | COS/HEK                            |
| A2A                          | *    | HEKT           | COS/HEK                            |
| A2A                          |      | HEK            | 500 G418                           |
| A2B                          | *    | HEKT           | COS/HEK                            |
| A3                           | *    | HEKT           | COS/HEK                            |
| <b>Melanocortin</b>          |      |                |                                    |
| MC-1                         | *    | HEKT           | COS/HEK                            |
| MC-2                         | *    | HEKT           | COS/HEK                            |
| MC-3                         | *    | HEKT           | COS/HEK                            |
| MC-4                         | *    | HEKT           | COS/HEK                            |
| MC-5                         | *    | HEKT           | COS/HEK                            |
| <b>Purinergic P2Y</b>        |      |                |                                    |
| P2Y1                         |      | Astrocyte line | 500 G418                           |
| P2Y2                         |      | Astrocyte line | 500 G418                           |
| P2Y4                         |      | Astrocyte line | 500 G418                           |
| P2Y6                         |      | Astrocyte line | 500 G418                           |
| P2Y11                        |      | Astrocyte line | 500 G418                           |
| P2Y12                        |      | Astrocyte line | 500 G418                           |
| <b>Trace Amine</b>           |      |                |                                    |
| TA-1                         | *    | HEKT           | COS/HEK                            |
| TA-3                         | *    | HEKT           | COS/HEK                            |
| TA-4                         | *    | HEKT           | COS/HEK                            |
| TA-5                         | *    | HEKT           | COS/HEK                            |
| <b>Lysophospholide (LPA)</b> |      |                |                                    |
| LPA1                         | *    | HEKT           | COS/HEK                            |
| LPA2                         | *    | HEKT           | COS/HEK                            |
| LPA3                         | *    | HEKT           | COS/HEK                            |
| <b>Tachykinin (NK)</b>       |      |                |                                    |
| NK1                          |      | HEK            | 500 G418                           |
| NK2                          |      | HEK            | 500 G418                           |
| NK3                          |      | HEK            | 500 G418                           |
| <b>mGluRs</b>                |      |                |                                    |
| mGluR1 (in progress)         | *    | HEKT           | 500 G418                           |
| mGluR2 (in progress)         | *    | HEKT           | 500 G418                           |
| mGluR3 (in progress)         | *    | HEKT           | 500 G418                           |
| mGluR4 (in progress)         | *    | HEKT           | 500 G418                           |

| Receptor                    | Note | Parental cells              | Media (see detail below the table) |
|-----------------------------|------|-----------------------------|------------------------------------|
| mGluR5                      |      | CHO                         | 2 µg/ml Puromycin                  |
| mGluR5 (in progress)        | *    | HEKT                        | 500 G418                           |
| mGluR6 (in progress)        | *    | HEKT                        | 500 G418                           |
| mGluR7 (in progress)        | *    | HEKT                        | 500 G418                           |
| mGluR8 (in progress)        | *    | HEKT                        | 500 G418                           |
| Others                      |      |                             |                                    |
| Ghrelin                     |      | HEK Flp-In                  | DMEM 100 µg/ml Hygromycin B        |
| PAR1                        |      | Lung Fibroblast,<br>PAR1 KO | 500 G418                           |
| SMO                         | *    | HEKT                        | COS/HEK                            |
| SMO (in progress)           |      | HEK                         | 500 G418                           |
| CCK2                        |      | CHO                         | 500 G418                           |
| Orexin-2                    | *    | HEKT                        | COS/HEK                            |
| GPR58                       | *    | HEKT                        | COS/HEK                            |
| GPR61                       | *    | HEKT                        | COS/HEK                            |
| GPR62                       | *    | HEKT                        | COS/HEK                            |
| GPR40                       | *    | HEKT                        | COS/HEK                            |
| GPR41                       | *    | HEKT                        | COS/HEK                            |
| GPR43                       | *    | HEKT                        | COS/HEK                            |
| Il-1 imidazoline            | *    | HEKT                        | COS/HEK                            |
| HERG-K <sup>+</sup> Channel |      | HEK                         | 500ug/ml G418                      |
| PC12                        |      |                             | COS/HEK                            |

### 1.2.5. Media Recipes

Details of media composition for cells are as follows. DMEM is purchased as 1X high D-glucose (4,500 mg/L) with L-glutamine and supplemented upon receipt with penicillin/streptomycin (100 U/ml).

**a. COS/HEK Medium (10% serum) – 1 L: also used for transient transfection**

1L DMEM, 100 ml FBS, 10 ml Pen/Strep

**b. Standard G418 Media (500 µg/ml; 10% serum) – 1L**

1 L DMEM, 500 mg G418, 100 ml FBS, 10 ml Pen/Strep

**c. Standard Puromycin Medium (2 µg/ml Puromycin:10% Serum)-1L**

1 L DMEM, 2 aliquots (1mg/ml Puromycin in PBS), 100 ml FBS, 10 ml Pen/Strep

**d. Standard Puromycin Medium +G418 Selection Medium**

1L DMEM, 2 aliquots (1mg/ml Puromycin in PBS), 500 mg geneticin, 100 ml FBS, 10 ml Pen/Strep

**e. Dialyzed Medium (1% dialyzed serum) – 1 L**

1 L DMEM, 10 ml Dialyzed FBS, 10 ml Pen/Strep

**f. hNET Medium (250 µg/mL G418; 10% serum) – 1 L**

1 L DMEM, 250 mg geneticin, 100 ml FBS, 10 ml Pen/Strep

**g. hDAT Medium (350 µg/ml G418; 10% serum) - 1L**

1L DMEM, 350mg geneticin, 100ml FBS, 10ml Pen/Strep

**h. V1A and OT Medium (400 µg/ml G418, 10% FBS, 15 mM HEPES) – 1L**

1L Hams F12, 400 mg Geneticin, 100 ml Calf serum, 15 ml HEPES (1M stock), 10 ml Pen/Strep

**i. V2 & V1B Medium (150 µg/mL Zeocin; 10% Serum; 15 mM HEPES) – 500 mL**

1L Hams F12, 1.5 ml Zeocin (100mg/mL), 100ml Calf Serum, 15 ml HEPES (1M), 10 ml Pen/Strep

**j. MOR/DOR selection medium (200 µg/ml G418, 10% serum) – 1L**

1 L DMEM, 200 mg G418, 100 ml FBS, 10 ml Pen/Strep

**k. Alpha 1A and 1D selection medium (500 µg/ml G418, 10% serum) – 1L**

1 L DMEM, 500 mg G418, 100 ml FBS, 10 ml Pen/Strep

**l. M4 medium (10% serum) -1L**

1L Hams F12, 100 ml FBS, 10 ml Pen/Strep

**m. D4 medium (10% serum) – 1L**

1L DMEM/F12, 15mM HEPES; with pyridoxine HCl, 100 ml Donor Calf Serum with Iron, 10 ml Pen/Strep

**n. FLP-In CHO medium (10% serum) – 1L**

1L DMEM/F12 50/50, 10% FBS, 200 µg/ml Hygromycin B, 10 ml Pen/Strep

**o. FLP-In HEK medium (10% serum)– 1L**

1L DMEM, 10% FBS, 100 µg/ml Hygromycin B, 10 ml Pen/Strep

**p. PC-12 medium (10% serum) – 1L**

1L DMEM, 5% Horse Serum, 5% FBS, 10 ml Pen/Strep

**q. HTLA medium – 1L**

1L DMEM, 100 ml FBS, 100 µg/ml Hygromycin B, 5 µg/ml Puromycin, 10 ml Pen/Strep

**r. Inositol/Inositol-free medium**


500 ml BME (LONZA), 5 ml Pen/Strep, 5 ml L-Glutamine

**s. McCoy medium for U2OS – 1L**

1L McCoy, 100 ml dialyzed FBS, 50 µg/ml Hygromycin B, 200 µg/ml Zeocin, 100 µg/ml G418, 0.1 mM non-essential amino acids, 25 mM HEPES, 1 mM Sodium Pyruvate, 10 ml Pen/Strep

### 1.3. Saturation binding assays

Saturation binding assays are usually performed immediately after the membrane fraction is obtained and protein concentration is determined (see above section for membrane preparations) to measure receptor expression level ( $B_{max}$ ) and binding affinity ( $K_d$ ) of a selected radioligand. **Tables 2-24** list the identity of hot ligands, reference compounds, and buffers for each family of targets. Saturation binding assays are carried out in 96-well plates in a final volume of 125  $\mu$ l per well. In brief, 25  $\mu$ l of radioligand is added to each well of a 96-well plate according to the setup as in [Figure 3](#); followed by addition of 25  $\mu$ l binding buffer (for total binding) or 25  $\mu$ l reference compound at final of 10  $\mu$ M (for nonspecific binding). The reaction starts upon addition of 75  $\mu$ l of fresh membrane protein (typically 25 to 50  $\mu$ g per well) and the reaction is usually incubated in the dark at room temperature for 90 min. The reaction is stopped by vacuum filtration onto cold 0.3% polyethyleneimine (PEI) soaked 96-well filter mats using a 96-well Filtermate harvester, followed by three washes with cold wash buffers (for detail, see **Tables 2-24**). Scintillation cocktail is then melted onto the microwave-dried filters on a hot plate and radioactivity is counted in a Microbeta counter.



**Figure 3.** Radioligand saturation binding plate set-up in a 96-well plate format. Total and nonspecific binding are determined in the absence and presence of 10  $\mu$ M of appropriate reference compound, respectively, in the indicated final concentrations of hot ligand (nM). Final radioligand concentrations are shown as an example, and different concentration ranges are used for different receptors depending on their  $K_d$  values for selected radioligand. The  $K_d$  value is usually in the middle row of the plate. Amount of membrane protein ( $\mu$ g/well) is adjusted to a point that bound radioactivity is less than 10% of total added radioactivity.

#### 1.4. Primary and secondary radioligand binding assays

Compounds are typically subjected to primary radioligand binding assays at targets selected by investigators and approved by the PDSP Director, Dr Bryan Roth. In the primary binding assays, compounds are usually tested at single concentrations (10  $\mu$ M) in quadruplicate in 96-well plates (see **Figure 1** for drug plate map). Compounds that show a minimum of 50% inhibition at 10  $\mu$ M

are tagged for secondary radioligand binding assays to determine equilibrium binding affinity at specific targets. In the secondary binding assays, selected compounds are usually tested at 11 concentrations (0.1, 0.3, 1, 3, 10, 30, 100, 300 nM, 1, 3, 10  $\mu$ M) and in triplicate (3 sets of 96-well plates, see **Figure 2a** for drug plate map). An alternative drug plate set up is shown in **Figure 2b** which contains 2 PDSP compounds and 1 reference compound. The drug map in **Figure 2a** uses resources more efficiently than the one in **Figure 2b**; drug plates such as those in **Figure 2b** are usually made manually for urgent projects with a small number of high-priority compounds to be tested immediately upon request.

Both primary and secondary radioligand binding assays are carried out in a final volume of 125  $\mu$ l per well in appropriate binding buffer (for details, see **Tables 2-24**). The hot ligand concentration is usually at a concentration close to the  $K_d$  (unless otherwise indicated). Total binding and nonspecific binding are determined in the absence and presence of 10  $\mu$ M appropriate reference compound, respectively. In brief, plates are usually incubated at room temperature and in the dark for 90 min (unless otherwise indicated). Reactions are stopped by vacuum filtration onto 0.3% polyethyleneimine (PEI) soaked 96-well filter mats using a 96-well Filtermate harvester, followed by three washes with cold wash buffer (for detail, see **Tables 2-24**). Scintillation cocktail is then melted onto the microwave-dried filters on a hot plate and radioactivity is counted in a Microbeta counter. A general procedure for binding assays is shown below (details for radioligand binding assays for nAChRs are described in a separate section).

1. Retrieve barcoded primary or secondary drug plates from the Hamilton STAR<sup>®</sup> or cold room
2. Check out drug plates and mark them “in progress” in the PDSP database
3. Prepare appropriate binding and wash buffers
4. Check out membrane pellets
5. Check pellet box for  $K_d$  value and concentration to be used for the radioligand
6. Prepare 2.5x of final concentration of radioligand working solution
7. Count 50  $\mu$ l of radioligand working solution to confirm concentration and activity
8. Add 50  $\mu$ l radioligand to each well

9. Add 50  $\mu$ l membrane suspension to each well
10. Mix by gentle and brief shaking
11. Incubate the plates in the dark for desired period of time (usually 60 – 90 min at RT)
12. Soak filters in cold 0.3% PEI
13. Stop the reaction by vacuum filtration and washing
14. Microwave-dry the filters
15. Melt scintillation cocktail on top of filters
16. Wrap up filters in plastic wrap
17. Count radioactivity
18. Download results, process and upload to the PDSP database
19. Preview and submit results
20. Report progress in the PDSP database (completion or redo or repeat) after receiving approving email

**1.5. Radioligand binding assays for nicotinic acetylcholine receptors (nAChRs).** Radioligand binding assays with nAChRs follow slightly different protocols from those outlined in the above sections and are detailed in the following sections.

#### **Main references for radioligand binding assays with nAChRs**

- Xiao Y, Meyer EL, Thompson JM, Surin A, Wroblewski J, Kellar KJ (1998). Rat  $\alpha$ 3/ $\beta$ 4 subtype of neuronal nicotinic acetylcholine receptor stably expressed in a transfected cell line: pharmacology of ligand binding and function. *Mol Pharmacol* 54(2): 322-333.
- Xiao Y, Fang H, Musachio JL, Wei ZL, Chellappan SK, Kozikowski AP, Kellar KJ (2006). Sazetidine-A, a novel ligand that desensitizes  $\alpha$ 4 $\beta$ 2 nicotinic acetylcholine receptors without activating them. *Mol Pharmacol* 70(4): 1454-1460.

**1.5.1. Cell culture.** The six stable cell lines expressing human  $\alpha$ 2 $\beta$ 2,  $\alpha$ 2 $\beta$ 4,  $\alpha$ 3 $\beta$ 2,  $\alpha$ 3 $\beta$ 4,  $\alpha$ 4 $\beta$ 2, or  $\alpha$ 4 $\beta$ 4 nAChR subtypes were established by stably co-transfecting HEK293 cells with a combination of one  $\alpha$  nicotinic receptor subtype gene and one  $\beta$  subunit gene. The cell line

expressing  $\alpha 7$  nAChRs was established by stably transfecting HEK293 cells with the rat  $\alpha 7$  nAChR subunit. Cells are grown in minimum essential medium (MEM) supplemented with 10% FBS, 100 units/ml penicillin G, 100  $\mu\text{g}/\text{ml}$  streptomycin, and selective antibiotics at  $37^\circ\text{C}$  with 5%  $\text{CO}_2$  in a humidified incubator.

**1.5.2. Nicotinic receptor radioligand binding assays.** Radioligand binding assays of nAChRs to [ $^3\text{H}$ ]-epibatidine use the above stably expressed nAChRs, or rat forebrain tissue. In brief, cells stably expressing nAChRs are harvested in 50 mM Tris HCl (pH 7.4), washed, homogenized with a Brinkmann polytron homogenizer and centrifuged at 36,000 x g. The resulting washed membranes are then incubated with [ $^3\text{H}$ ]-epibatidine for 4h at room temperature in a final volume of 0.5 ml. Nonspecific binding is assessed in parallel incubations in the presence of 300  $\mu\text{M}$  nicotine. Bound and free ligands are separated by vacuum filtration through Whatman GF/C filters treated with 0.5% polyethylenimine. The filter-retained radioactivity is measured by liquid scintillation counting. Specific binding is defined as the difference between total binding and nonspecific binding. [Examples are listed in Table 4 and Figure 12.](#)

**1.5.3. Primary and secondary nAChR binding assays.** Primary nAChR binding assays are performed with 100 pM [ $^3\text{H}$ ]-epibatidine and a single concentration of PDSP compound (10  $\mu\text{M}$ ) in quadruplicate. Results are expressed as percentage inhibition of [ $^3\text{H}$ ]-epibatidine specific binding. Compounds with a minimum of 25% inhibition are subjected to secondary binding assays to determine binding affinity. In the secondary binding assays, 0.5 nM [ $^3\text{H}$ ]-epibatidine and 10 concentrations of test PDSP compounds in singlets are tested to generate a competition binding curve. Results are analyzed in Prism 5.0 by nonlinear least-squares regression to obtain  $K_i$  values. Nicotine at 300  $\mu\text{M}$  is included in all assays to define nonspecific binding, and a nicotine concentration-response curve is included as a positive control in all the secondary binding assays.

**1.5.4. Nicotinic acetylcholine receptor binding data entry and reporting.** Binding assay results with nAChRs are analyzed as indicated in **Section 1.6.2**, and the  $K_i$  values are then manually

entered into the PDSP database. Figures and raw data for these assays are available to appropriate investigators (who submitted the compounds) upon request.

### **1.6. General binding data entry, analysis, and quality control for targets other than nAChRs.**

Except radioligand binding results with nAChRs, whose analysis is outlined in **Section 1.5**, binding assays for targets other than nAChRs are analyzed and reported as outlined in the following sections. Raw cpm data from the Microbeta counters are uploaded into the PDSP database and analyzed online using built-in analysis tools in the PDSP database. If necessary, binding assay results are also analyzed in Microsoft Excel (for primary binding results) or Prism v5.0 (for saturation binding and secondary binding results). The following sections provide detailed procedures for analysis and reporting of binding results, using both built-in tools in the PDSP database, or third party software (Microsoft Excel and GraphPad Prism).

**1.6.1. Results analysis using built-in tools in the PDSP database.** The PDSP database is coded to analyze results from saturation binding and competition binding assays (primary and secondary binding assays). The code for analyzing saturation and competition binding results is written in PHP using external library JQuery UI. The non-linear regression calculation is computed by the **R** statistical programming language using an Automatic Differentiation Model Builder (ADMB) to support robust non-linear regression calculation. To use the built-in tools in the PDSP database to analyze binding assay results, technicians upload raw counting results in Excel spreadsheets from Microbeta counters directly to the PDSP database through corresponding reporting systems designed for saturation binding, primary binding, or secondary binding assays. The online reporting system and analysis tools are designed to optimize and streamline both the data processing and quality control.

**1.6.1.1. Analysis of saturation binding results.** Saturation binding assays are carried out whenever a new batch of membrane preparations is made. The protein amount per well is adjusted such that less than 10% total added radiolactivity is bound. Upon completion of saturation binding assays, technicians can upload results in Excel spreadsheets directly to the

membrane pellet section of the PDSP database. During the uploading and reporting process, technicians must provide (1) receptor identity and species information; (2) stable line or transient transfection or animal tissue; (3) amount of protein per well used; (4) radioligand identity and concentrations. An automatic email is sent to the designated PDSP administrator when saturation binding results are submitted and results are placed in a list pending review and approval as indicated in **Figure 4**. The PDSP database is coded to analyze saturation binding results and generate a hyperbolic curve using non-linear least squares curve-fitting algorithm to determine  $B_{\max}$  and  $K_d$  values. The PDSP administrator reviews and analyzes the binding results online with options to (1) reject the whole set of results; (2) exclude apparent outliers among replicates; (3) request a repeat; or (4) approve the data set. **Figure 5** shows a screen captured during review of a saturation binding assay. When the data set is approved, the PDSP administrator decides on (1) amount of protein per well and (2) radioligand concentration to be used for competition binding assays. Upon approval, the saturation binding curve is immediately available in the “approved membrane pellet” section with unique timestamp associated with each batch of membrane preparation. Technicians make aliquots accordingly so that each aliquot is sufficient for one plate of a binding assay at the decided concentration of radioligand. Aliquots are stored in 1.5 ml Eppendorf tubes at  $-80^{\circ}\text{C}$  in assorted boxes until use at specified concentration of radioligand, with one aliquot being used per assay plate.

**1.6.1.2. Analysis of primary binding assay results.** Upon completion of primary binding assays, technicians upload raw results in Excel spreadsheets directly to the PDSP database via the primary binding assay reporting menu. When uploading primary binding results, technicians must identify plates by barcodes and match the barcodes generated by the Hamilton STAR liquid handling system when the drug plates were made. The designated PDSP administrator double checks the plate identity and sample order before making the results available online for reviewing. The PDSP database is coded to calculate percentage inhibition in each assay plate with total binding (with buffer) as 0% inhibition and nonspecific binding (in the presence of 10  $\mu\text{M}$  of a reference compound) as 100% inhibition. It also calculates average values from 4

replicates for each sample. The PDSP database automatically marks those entries with 50% or higher inhibition from secondary binding assays and also highlights those with variances greater than 20% among the quadruplicate determinations for further inspection. In addition, those compounds with negative inhibition (more binding in the presence of test compound) by over 20% are also highlighted for inspection. For the results of any compound, there are four options available: “Redo” deletes the data from database and puts the compound back on the primary schedule list; “Repeat” accepts the results and puts the compound back on the primary schedule list for a repeat assay; “Force Secondary” accepts the primary results and adds the compound on the secondary schedule list for secondary assays, even if the primary percent inhibition is less than 50%; “Approve” accepts the results, and the database puts the compound on the secondary schedule list if primary inhibition is over 50%. Either “Repeat”, or “Force Secondary”, or “Approve” immediately makes primary binding results available online to investigators. A representative screen capture (**Figure 6**) shows a reviewing screen for primary binding assay results.

**1.6.1.3. Analysis of secondary binding assay results.** Upon completion of secondary binding assays, technicians upload raw results in Excel spreadsheets directly to the PDSP database via a secondary binding assay reporting menu. In contrast to the reporting of primary binding results, the PDSP database analysis tools can accept two types of 96-well drug plate layouts for this purpose, as indicated in [Figures 2a and 2b](#), with the plates being made either manually or by the Hamilton STAR<sup>®</sup>. When uploading and submitting secondary binding results, technicians must match drug plate maps and barcodes with those generated by Hamilton STAR<sup>®</sup> liquid handling system when the drug plates were made. They also must choose buffers, radioligand identity, and reference compound from pull-down menus and enter total radioligand activity (dpm), specific activity, and the  $K_d$  value (from a saturation binding assay, see above) for the radioligand; the later values will be used by the PDSP database to calculate the radioligand concentration to be used for the assay, and convert  $IC_{50}$  to  $K_i$  values. A representative screen capture (**Figure 7**) is shown below as it appears while technicians are uploading secondary binding results for

submission. Technicians can view and examine the preliminary curves on a plate or individual compound basis before submission.

After secondary binding results are submitted, the PDSP database pools results from 3 replicates of 96-well plates as a single file, assigns it a serial number with extra information such as receptor identity and PDSP compound numbers, and places the results in a pending list for review (**Figure 4**). When PDSP administrators review and analyze the results for approval by clicking on any set of results in the pending list, the PDSP database parses the uploaded file, and generates an editable table along with a preliminary graph for analysis purpose (**Figure 8**). Clicking on individual PDSP compound leads to a new screen (**Figure 9**) showing a curve of the selected compound along with the reference curve. The raw data for both the selected compound and reference compound are shown in tables at left of the screen. The PDSP database is also coded to run the **Robust regression and Outlier removal (ROUT)** algorithm to identify and highlight potential outliers. PDSP administrators decide whether to exclude or include the potential outliers when reviewing each data set. The raw data tables are also editable and allow for exclusion of points that are known to be errors, for example, if a mistake is noticed in pipetting radioligand or membrane by a technician.

As an internal control, each set of secondary binding assays contains a reference concentration-response curve, also done in triplicate. The reference  $K_i$  value is automatically compared with its historical average value and/or the value in our  $K_i$  database. If the reference  $K_i$  is more than 3 fold away from its historical average value, the entire data set is flagged for further inspection. Senior PDSP scientists work with technicians to discover the underlying causes for such discrepancies, and arrive at solutions. For example, the secondary binding assay may be repeated, or the old stock of reference compound is replaced with a fresh one, if the same reference compound has low affinity in multiple data sets.

Each curve has its  $K_i$  value listed at the right side of the screen and can be marked as either “Redo” or “Accept”. If a curve is marked as “Accept”, a new graph is generated for the PDSP

compound together a reference curve and status of the assay is flagged as “completed”; the result is immediately available to investigators online. Otherwise, “Redo” puts the compound back on to the secondary schedule list. When an assay is approved, the PDSP database checks whether that compound has been previously tested at the same receptor. If so, the database compares the results of the assays, and generates a warning email in the case that it finds that one  $K_i$  value is at least 4x different (0.6 log unit) from another one. All of the secondary assays with the same compound are put onto the “Experiments under review” list. PDSP administrators review the results, and decide how to proceed with the compounds. Options include (1) one more repeat with the same stock; (2) one more repeat with a fresh stock; (3) approval with no further action, since lowest  $K_i$  value could be 4x away from highest  $K_i$  value, but still be less than 3x away from mean value. If the PDSP database has multiple accepted entries for the same compound, the database calculates an average from all previous repeats. If a compound shows allosteric potentiation in radioligand binding assays, the result is usually extracted and analyzed separately using Prism v5.0 ([See representative figures for 5-HT1E secondary binding](#)).

### Saturation Binding Analysis Tool

Pellets pending approval:

| PelletID            | Receptor | Radioligand   | Date Created | Technician |
|---------------------|----------|---------------|--------------|------------|
| <a href="#">673</a> | 5-HT1A   | [3H]WAY100635 | 2013-03-19   | tom        |

PDSP Home [Ki DB Home](#) [Tracer Database](#) [Roth Lab](#) [U N C](#) [Binding Assay](#)  
[Home](#) [Log Out](#) Add... Edit... View...

### View Pending Experiments

#### Primary Binding Analysis Tool


| Experiment Name                   | Experiment Type           |
|-----------------------------------|---------------------------|
| <a href="#">TC312135HT2B</a>      | primaryBindingExperiments |
| <a href="#">TC312135HT2C</a>      | primaryBindingExperiments |
| <a href="#">TC31213ALPHA2A2</a>   | primaryBindingExperiments |
| <a href="#">TC31213ALPHA2C1</a>   | primaryBindingExperiments |
| <a href="#">TC31213D1</a>         | primaryBindingExperiments |
| <a href="#">TC31213D2</a>         | primaryBindingExperiments |
| <a href="#">TC31213D3</a>         | primaryBindingExperiments |
| <a href="#">TC31213D4</a>         | primaryBindingExperiments |
| <a href="#">TC31213D5</a>         | primaryBindingExperiments |
| <a href="#">899-Beta1ASMa tjm</a> | primaryBindingExperiments |

#### Secondary Binding Analysis Tool

Plates pending approval:

| PlateID              | Receptor                 | Compounds                                       | Technician |
|----------------------|--------------------------|---|------------|
| <a href="#">2336</a> | BZP Rat Brain Site 90001 |   | rwh        |
| <a href="#">2871</a> | 5-ht5a                   | 1, 2, 3, 4, 5, 6, 90001                         | sal        |
| <a href="#">4345</a> | DAT                      | 27103, 27266, 27267, 27101, 27364, 27365        | jfw        |
| <a href="#">4896</a> | DAT                      | 26610, 26611, 26614, 26615, 26616, 26617, 26618 | jfw        |

**Figure 4.** Representative pending lists for Saturation binding (Top), Primary binding (Middle), and Secondary binding (bottom) results. When binding assay results are reported and uploaded into the PDSP database, they are placed in separate pending lists.



**Figure 5.** A representative screen capture showing a reviewing screen of saturation binding results. Both total and nonspecific binding (cpm/well) are shown in the table on the left side of the screen. Specific binding is fitted to a hyperbolic curve to determine  $B_{max}$  (cpm/well) and  $K_d$  (nM); both values are listed at the right side of the curve. The pellet information is listed above the figure, showing that this is a 5-HT1A membrane pellet, made from a stable line, assayed at 50 µg/well protein, on March 19, 2013.

**View Pending Experiment**

Experiment: TC312135HT2C  
 Receptor: 5-HT2C  
 Date Completed: 2013-03-12

| Compound | Total 1 | Total 2 | Mean Specific | Mean % Inhibition | Total 3 | Total 4 | Mean Specific | Mean % Inhibition | Mean % Inhibition (N=4) | Redo                  | Repeat                | Force Secondary       | Approve               |
|----------|---------|---------|---------------|-------------------|---------|---------|---------------|-------------------|-------------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| None     | 2341    | 2128    | 2187          | 0                 | 2961    | 2100    | 2485          | 0.00              | 0                       |                       |                       |                       |                       |
| 27475    | 881     | 757     | 771.5         | 84.723385         | 1085    | 877     | 915.5         | 82.86             | 63.791702               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27476    | 1653    | 1982    | 1770          | 19.067215         | 2370    | 2057    | 2148          | 12.86             | 15.963627               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27478    | 1315    | 1156    | 1188          | 45.679012         | 1384    | 2008    | 1630.5        | 33.85             | 39.766483               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27782    | 303     | 464     | 336           | 84.636488         | 494     | 525     | 444           | 81.99             | 83.312158               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27796    | 2049    | 2303    | 2128.5        | 2.6748971         | 2285    | 2397    | 2275.5        | 7.69              | 5.1812619               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27797    | 1515    | 1894    | 1657          | 24.234110         | 2517    | 751     | 1568.5        | 36.37             | 30.301639               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27798    | 1407    | 1798    | 1555          | 28.898033         | 1913    | 2029    | 1905.5        | 22.70             | 25.797901               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27799    | 3434    | 2463    | 2901          | 32.64746          | 2149    | 2786    | 2402          | 2.56              | 15.04584                | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27800    | 428     | 537     | 435           | 80.109739         | 658     | 770     | 648.5         | 73.69             | 76.900711               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 27801    | 1728    | 1686    | 1664.5        | 23.891175         | 2126    | 920     | 1457.5        | 40.87             | 32.381693               | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

Figure 6. A representative screen capture showing a reviewing screen of primary binding assay results. The barcoded plate is for a 5-HT<sub>2C</sub> binding assay carried out on March 12, 2013. Total binding values are shown in the top row and PDSP compounds are shown in subsequent rows. The nonspecific binding values are at the bottom of the page and do not appear in the Figure. Highlighted in yellow are the samples with more than 50% inhibition and these samples are automatically tagged for secondary binding assays. Highlighted in red are the ones with more 20% variations among 4 replicates. There are four options for each compound: “Redo”, “Repeat”, “Force Secondary”, and “Approve”, as indicated at the right side of the screen capture.





Figure 7. A representative screen capture showing the final stage in an actual secondary binding result submission. On the left, one complete set of raw data are pooled from three 96-well plates (Plate 7, Plate 8, and Plate 9), and organized from Rows A to G (in triplicate) from top to bottom with PDSP compound # showing above the table. On the upper right, detailed information for the binding assay is shown, including target receptor identity (5-HT<sub>2A</sub> receptor in this case), K<sub>d</sub> value (1.6 nM) for the radioligand [<sup>3</sup>H]-ketanserin, total radioactivity (17939 dpm) added to each well, specific activity of the radioligand (50.3 Ci/mmol), and final volume of the assay (125 µl/well) in standard binding buffer. On the lower right, detailed information on the tested compounds (PDSP numbered compounds and reference compound with corresponding concentrations) are shown. The concentrations for each compound are set at default values, but can be changed individually by sliding the corresponding bar button, if necessary.




Re-Calculate

| Name      | Quality |
|-----------|---------|
| Clozapine |         |

|       | Reco                  | Accept                |
|-------|-----------------------|-----------------------|
| 27741 | <input type="radio"/> | <input type="radio"/> |
| 27746 | <input type="radio"/> | <input type="radio"/> |
| 27782 | <input type="radio"/> | <input type="radio"/> |
| 27748 | <input type="radio"/> | <input type="radio"/> |
| 27750 | <input type="radio"/> | <input type="radio"/> |
| 27752 | <input type="radio"/> | <input type="radio"/> |
| 27753 | <input type="radio"/> | <input type="radio"/> |

Finalize Plate

**Figure 8.** A screen capture showing a representative overview of a set of secondary binding results. On the left, counts (cpm/well) for the reference compounds are shown in a table. In the middle, up to 8 curves are shown on the same graph with a pre-determined color scheme with PDSP numbers above the graph (one on each tab) and corresponding  $pK_i$  and  $K_i$  values on the right side.



**Figure 9.** A screen capture showing a representative reviewing process of a secondary binding assay. On the right, counts (cpm/well) for both a reference compound and a selected PDSP compound are shown in corresponding tables. One potential outlier has been highlighted (in dark gray with a checkmark “✓”) by the PDSP database using a robust outlier identification algorithm. In the middle, the selected compound is plotted side by side with the reference compound, and corresponding pK<sub>i</sub> and K<sub>i</sub> values are listed on the right side. At the lower left corner, each compound can be marked either “Redo” or “Accept”, and the long box can accommodate a brief text note.

**1.6.2. Analysis of binding assay results using Prism v5.0:** Binding assay results are also analyzed using Prism v5.0, if necessary, using Prism's built-in functions for corresponding assay types. The following sections provide detailed procedures for analysis of binding data.

**1.6.2.1. Analysis of saturation binding assays.** For saturation binding results, total binding and nonspecific binding results are analyzed in Prism v5.0 by fitting results to the following equations to determine  $B_{max}$  and  $K_d$  values. To do this, total binding in cpm is entered in column A and nonspecific binding in cpm is entered in column B, corresponding to concentrations of free radioligand.

$$\begin{aligned} \text{Nonspecific binding} &= NS * X + \text{Background} \\ \text{Total binding} &= \text{Nonspecific binding} + \frac{B_{max} * X}{(X + K_d)} \end{aligned}$$

in which "**Nonspecific binding**" and "**Total binding**" are measured radioactivity (cpm per well) in the absence and presence of 10  $\mu$ M of a reference compound, respectively, at the corresponding concentration [**X**] in nM.  $B_{max}$  is the receptor expression level in cpm/well or fmol/mg membrane protein;  $K_d$  is the equilibrium binding affinity, corresponding to the radioligand concentration at which point 50% of  $B_{max}$  is bound to the radioligand; **NS** and **Background** are two fitting values for nonspecific binding and background counts.

**1.6.2.2. Analysis of primary binding assays.** For primary binding results, non-specific binding in the presence of 10  $\mu$ M of an appropriate reference compound is set as 100% inhibition; total binding in the absence of test compound or reference compound is set as 0% inhibition. The radioactivity in the presence of test compound is calculated with the following equation and expressed as a percent inhibition. The normalization process is carried out in Prism or Excel.

$$\% \text{ Inhibition} = 100 - \frac{\text{Sample cpm} - \text{Nonspecific cpm}}{\text{Total cpm} - \text{nonspecific cpm}} \times 100$$

**1.6.2.3. Analysis of secondary binding assays.** For secondary binding results, counts (cpm/well) are pooled and fitted to a three-parameter logistic function for competition binding in Prism v 5.0 to determine  $IC_{50}$  values,

$$Y = Bottom + \frac{(Top - Bottom)}{1 + 10^{X - LogIC_{50}}}$$

in which **Y** is total binding in the presence of corresponding concentration of testing drug (**X**); **Top** and **Bottom** are the total and nonspecific binding in the absence and presence of 10  $\mu$ M of the reference compound;  $IC_{50}$  is the concentration at which 50% observed binding is inhibited and is converted to  $K_i$  according to Cheng-Prusoff equation,

$$K_i = \frac{IC_{50}}{1 + \frac{L}{K_d}} \text{ or } LogK_i = LogIC_{50} - Log\left(1 + \frac{L}{K_d}\right)$$


in which **L** is the radioligand concentration used in the competition binding assay;  $K_d$  is the radioligand equilibrium binding affinity determined in the above saturation binding assays. In the curve-fitting analysis, top and bottom values are shared among all binding curves from the same plate if necessary.



## 1.7. Tables of binding assay conditions and representative figures



**Table 2.** 5-HT receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value. The  $K_d$  values listed in this table are the mean  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. BB for binding buffer; WB for wash buffer. Historical reference  $K_i$  values from the last 6 months are also listed.



| 5-HT receptors   |                              |  |                    |                      |
|--|------------------------------|--|--------------------|----------------------|
| Standard BB (SBB): 50 mM Tris HCl, 10 mM MgCl <sub>2</sub> , 0.1 mM EDTA, pH 7.4, RT |                              |  |                    |                      |
| Standard WB (SWB): 50 mM Tris HCl, pH 7.4, cold                                      |                              |  |                    |                      |
| Target   | Radioligand                  | $K_d$ or [ <sup>3</sup> H] for binding in nM (N) | Reference Compound | Reference $K_i$ (nM) |
| 5-HT1A   | [ <sup>3</sup> H]Way100635   | 0.43 $\pm$ 0.06 (21)                             | 8-OH-DPAT          | 0.60 $\pm$ 0.10      |
| 5-HT1B   | [ <sup>3</sup> H]5-CT        | 1.91 $\pm$ 0.29 (8)                              | Ergotamine         | 4.75 $\pm$ 0.51      |
| 5-HT1D   | [ <sup>3</sup> H]5-CT        | 1.47 $\pm$ 0.24 (13)                             | Ergotamine         | 5.04 $\pm$ 0.47      |
| 5-HT1E   | [ <sup>3</sup> H]5-HT        | 3.36 $\pm$ 0.46 (8)                              | 5-HT               | 11.8 $\pm$ 1.2       |
| 5-HT1F   |                              |  |                    |                      |
| 5-HT2A   | [ <sup>3</sup> H]Ketanserin  | 1.57 $\pm$ 0.26 (19)                             | Clozapine          | 9.50 $\pm$ 1.03      |
| 5-HT2B   | [ <sup>3</sup> H]LSD         | 1.04 $\pm$ 0.15 (19)                             | SB206553           | 18.2 $\pm$ 2.3       |
| 5-HT2C   | [ <sup>3</sup> H]Mesulergine | 2.92 $\pm$ 0.83 (13)                             | Ritanserin         | 2.12 $\pm$ 0.22      |
| 5-HT3  | [ <sup>3</sup> H]GR65630     | 1.0 for binding                                  | Zacopride          | 0.65 $\pm$ 0.08      |
| 5-HT4  | [ <sup>3</sup> H]GR113808    | 0.2 – 0.5 for binding                            | SD205557           | 36.5 $\pm$ 8.5       |
| 5-HT5A   | [ <sup>3</sup> H]LSD         | 2.32 $\pm$ 0.14 (36)                             | Ergotamine         | 36.5 $\pm$ 3.9       |
| 5-HT6  | [ <sup>3</sup> H]LSD         | 3.97 $\pm$ 0.46 (10)                             | Chlorpromazine     | 10.5 $\pm$ 1.6       |
| 5-HT7A   | [ <sup>3</sup> H]LSD         | 5.30 $\pm$ 0.70 (2)                              | Chlorpromazine     | 12.9 $\pm$ 0.9       |



Figure 10. Representative competition binding curves with 5-HT receptors






















**Table 3.** Muscarinic acetylcholine receptors, radioligands and corresponding  $K_d$  values, reference compounds, and buffers for primary and secondary radioligand binding assays. For muscarinic receptors, we usually use 0.5 – 1.0 nM  $^3\text{H}$ -QNB or  $^3\text{H}$ -NMS for inhibition binding assays. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. BB for binding buffer; WB for wash buffer. Historical reference  $K_i$  values from last the 6 months are also included for quality control.

| Muscarinic acetylcholine receptors (mAChRs)  |                     |                      |            |                      |
|--|---------------------|----------------------|------------|----------------------|
| Muscarinic BB (MBB) #1: 50 mM Tris HCl, pH 7.7, RT + SWB   |                     |                      |            |                      |
| Muscarinic BB (MBB) #2: 25 mM Sodium Phosphate, 5 mM MgCl <sub>2</sub> , pH 7.4, RT (cold for washing) |                     |                      |            |                      |
| Muscarinic wash buffer #2: Same as muscarinic binding buffer #2, cold                                  |                     |                      |            |                      |
| Use 0.5 – 1.0 nM hot ligand for binding assays   |                     |                      |            |                      |
| Target   | Radioligand         | $K_d$ in nM (N)      | References | Reference $K_i$ (nM) |
| M1   | [ $^3\text{H}$ ]QNB | 1.02 $\pm$ 0.17 (9)  | Atropine   | 1.80 $\pm$ 0.12      |
| M1   | [ $^3\text{H}$ ]NMS | 0.90 $\pm$ 0.11 (5)  | Atropine   |                      |
| M2   | [ $^3\text{H}$ ]QNB | 0.34 $\pm$ 0.06 (28) | Atropine   | 1.52 $\pm$ 0.42      |
| M2   | [ $^3\text{H}$ ]NMS | 0.34 $\pm$ 0.03 (5)  | Atropine   |                      |
| M3   | [ $^3\text{H}$ ]QNB | 0.57 $\pm$ 0.05 (17) | Atropine   | 0.86 $\pm$ 0.12      |
| M3   | [ $^3\text{H}$ ]NMS | 0.36 $\pm$ 0.10 (7)  | Atropine   |                      |
| M4   | [ $^3\text{H}$ ]QNB | 0.27 $\pm$ 0.05 (28) | Atropine   | 0.86 $\pm$ 0.13      |
| M4   | [ $^3\text{H}$ ]NMS | 0.19 $\pm$ 0.07 (7)  | Atropine   |                      |
| M5   | [ $^3\text{H}$ ]QNB | 0.67 $\pm$ 0.13 (21) | Atropine   | 0.88 $\pm$ 0.08      |
| M5   | [ $^3\text{H}$ ]NMS | 0.45 $\pm$ 0.08 (12) | Atropine   |                      |

**Figure 11.** Representative competition binding curves with muscarinic acetylcholine receptors











**Table 4.** Radioligand binding assays for nicotinic acetylcholine receptors (nAChRs). Primary binding assays use 100 pM of [<sup>3</sup>H]-epibatidine and secondary binding assays use 0.5 nM [<sup>3</sup>H]-epibatidine. Historical reference K<sub>i</sub> values from the last 6 months are also included for quality control.



| Nicotinic acetylcholine receptors (nAChRs)                                  |                              |                               |           |                               |
|---|------------------------------|-------------------------------|-----------|-------------------------------|
| Nicotinic acetylcholine receptor binding buffer: 50 mM Tris HCl, pH 7.4, RT |                              |                               |           |                               |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold                          |                              |                               |           |                               |
| Target  | Radioligand                  | [ <sup>3</sup> H] for binding | Reference | Reference K <sub>i</sub> (nM) |
| α2β2  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 10.4 ± 2.2                    |
| α2β4  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 113.0 ± 7.4                   |
| α3β2  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 62.8 ± 5.4                    |
| α3β4  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 466.2 ± 6.4                   |
| α4β2  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 10.5 ± 1.3                    |
| α4β2*   | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 10.9 ± 1.0                    |
| α4β4  | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 35.8 ± 2.0                    |
| α7, rat   | [ <sup>3</sup> H]Epibatidine | 0.5 nM                        | Nicotine  | 690                           |

\* Rat forebrain

**Figure 12.** Representative competition binding curves with nicotinic acetylcholine receptors








**Table 5.** Adenosine receptors, radioligand and corresponding concentration, reference compound, and buffer for primary and secondary radioligand binding assays.



| Adenosine receptors   |                        |                                     |            |                               |
|---|------------------------|-------------------------------------|------------|-------------------------------|
| Adenosine Binding Buffer: 50 mM Tris HCl, 1U/ml adenosine deaminase, pH 7.4, RT |                        |                                     |            |                               |
| Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold                              |                        |                                     |            |                               |
| Target  | Radioligand            | [ <sup>3</sup> H] for binding in nM | References | Reference K <sub>i</sub> (nM) |
| A1  | [ <sup>3</sup> H] NECA | 3 nM                                | NECA       |                               |
| A2A   | [ <sup>3</sup> H] NECA | 3 nM                                | NECA       |                               |
| A2B   | [ <sup>3</sup> H] NECA | 3 nM                                | NECA       |                               |
| A3  | [ <sup>3</sup> H] NECA | 3 nM                                | NECA       |                               |



**Table 6.** Adrenergic receptors, radioligands and corresponding  $K_d$  values, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentrations of radioligand used for competition binding assays are usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also listed.



| Adrenergic receptors   |                              |  |               |                      |
|--|------------------------------|--|---------------|----------------------|
| $\alpha$ 1 Binding Buffer: 20 mM Tris HCl, 145 mM NaCl, pH 7.4, RT             |                              |  |               |                      |
| $\alpha$ 2 Binding Buffer: 50 mM Tris HCl, 5 mM MgCl <sub>2</sub> , pH 7.7, RT |                              |  |               |                      |
| $\beta$ Binding Buffer: 50 mM Tris HCl, 3 mM MgCl <sub>2</sub> , pH 7.7, RT    |                              |  |               |                      |
| Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold                             |                              |  |               |                      |
| Target   | Radioligand                  | $K_d$ or [ <sup>3</sup> H] for binding in nM (N) | References    | Reference $K_i$ (nM) |
| $\alpha$ 1A  | [ <sup>3</sup> H]Prazosin    | 0.64 $\pm$ 0.14 (8)                              | Prazosin      | 0.95 $\pm$ 0.08      |
| $\alpha$ 1B  | [ <sup>3</sup> H]Prazosin    | 0.5 – 1.0 nM for binding                         | Prazosin      | 1.28 $\pm$ 0.13      |
| $\alpha$ 1D  | [ <sup>3</sup> H]Prazosin    | 0.89 $\pm$ 0.14 (7)                              | Prazosin      | 1.04 $\pm$ 0.10      |
| $\alpha$ 2A  | [ <sup>3</sup> H]Rauwolscine | 0.5 - 1.0 nM for binding                         | Oxymetazoline | 5.15 $\pm$ 0.21      |
| $\alpha$ 2B  | [ <sup>3</sup> H]Rauwolscine | 1.36 $\pm$ 0.50 (7)                              | Prazosin      | 5.86 $\pm$ 0.52      |
| $\alpha$ 2C  | [ <sup>3</sup> H]Rauwolscine | 0.73 $\pm$ 0.24 (6)                              | Oxymetazoline | 36.9 $\pm$ 3.3       |
| $\beta$ 1  | [ <sup>3</sup> H]CGP12177    | 1.0 nM for binding                               | Alprenolol    | 1.86 $\pm$ 0.21      |
| $\beta$ 1  | [ <sup>125</sup> I]Pindolol  | 0.05 nM for binding                              | Alprenolol    |                      |
| $\beta$ 2  | [ <sup>3</sup> H]CGP12177    | 1.03 $\pm$ 0.47 (6)                              | Alprenolol    | 2.03 $\pm$ 0.28      |
| $\beta$ 3  | [ <sup>3</sup> H]CGP12177    | 13.8 $\pm$ 9.1 (5)<br>10 nM for binding          | Propranolol   | 49.7 $\pm$ 14.1      |
| $\beta$ 3  | [ <sup>125</sup> I]Pindolol  | 0.5 nM for binding                               | Propranolol   |                      |


Figure 13. Representative competition binding curves with adrenergic receptors












**Table 7.** Cannabinoid receptors, radioligand and corresponding concentration, reference compound, and buffers for primary and secondary radioligand binding assays. Historical reference  $K_i$  values from the last 6 months are also included.

| Cannabinoid receptors   |                          |                                     |                       |                          |
|---|--------------------------|-------------------------------------|-----------------------|--------------------------|
| Cannabinoid Binding Buffer: 50 mM Tris HCl, 5 mM MgCl <sub>2</sub> , 1 mM EDTA, 1 mg/ml BSA, pH 7.4, RT |                          |                                     |                       |                          |
| Cannabinoid Wash Buffer: cannabinoid binding buffer + 1 mg/ml BSA, pH 7.4, cold                         |                          |                                     |                       |                          |
| Target  | Radioligand              | [ <sup>3</sup> H] for binding in nM | References            | Reference $K_i$ (nM)     |
| CB1 (rat brain)   | [ <sup>3</sup> H]CP55940 | 1 nM                                | WIN55212-3            | 18.5 ± 3.0               |
| CB2   | [ <sup>3</sup> H]CP55940 | 2 nM                                | WIN55212-3<br>CP55940 | 9.4 ± 2.6<br>4.95 ± 0.55 |



Figure 14. Representative competitive binding curves for cannabinoid receptors.






**Table 8.** Dopamine receptors, radioligands and corresponding  $K_d$  values, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

| Dopamine receptors   |                                    |                                   |                |                      |
|--|------------------------------------|-----------------------------------|----------------|----------------------|
| Dopamine Binding Buffer: 50 mM HEPES, 50 mM NaCl, 5 mM MgCl <sub>2</sub> , 0.5 mM EDTA, pH 7.4, RT |                                    |                                   |                |                      |
| Standard Wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                                    |                                   |                |                      |
| Target   | Radioligand                        | [ <sup>3</sup> H] $K_d$ in nM (N) | References     | Reference $K_i$ (nM) |
| D1   | [ <sup>3</sup> H]SCH23390          | 0.89 $\pm$ 0.12 (16)              | SKF38393       | 3.97 $\pm$ 0.37      |
| D2   | [ <sup>3</sup> H]N-methylspiperone | 0.67 $\pm$ 0.11 (21)              | Haloperidol    | 7.20 $\pm$ 0.75      |
| D3   | [ <sup>3</sup> H]N-methylspiperone | 0.66 $\pm$ 0.24 (7)               | Chlorpromazine | 10.30 $\pm$ 1.08     |
| D4   | [ <sup>3</sup> H]N-methylspiperone | 0.84 $\pm$ 0.09 (9)               | Chlorpromazine | 37.37 $\pm$ 4.21     |
| D5   | [ <sup>3</sup> H]SCH23390          | 2.08 $\pm$ 0.27 (18)              | SKF38393       | 3.08 $\pm$ 0.18      |

Figure 15. Representative competitive binding curves with Dopamine receptors.











**Table 9.** Histamine receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

| Histamine receptors   |  |   |                  |                      |
|---|--|---|------------------|----------------------|
| Histamine Binding Buffer: 50 mM Tris HCl, 0.5 mM EDTA, pH 7.4, RT |  |   |                  |                      |
| Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold                |  |   |                  |                      |
| Filter: GF/B  |  |   |                  |                      |
| Target  | Radioligand                                | $K_d$ or [ $^3\text{H}$ ] for binding in nM (N) | References       | Reference $K_i$ (nM) |
| H1  | [ $^3\text{H}$ ]Pyrilamine                 | $1.01 \pm 0.13$ (21)                            | Chlorpheniramine | $3.22 \pm 0.33$      |
| H2  | [ $^3\text{H}$ ]Cimetidine                 | 3 nM for binding                                | ORG-5222         | $3.05 \pm 0.47$      |
| H2  | [ $^3\text{H}$ ]Tiotidine                  | 2 nM for binding                                | Cimetidine       |                      |
| H3  | [ $^3\text{H}$ ] $\alpha$ -methylhistamine | $0.74 \pm 0.11$ (19)                            | Histamine        | $6.43 \pm 1.60$      |
| H4  | [ $^3\text{H}$ ]Histamine                  | $1.20 \pm 0.10$ (2)                             | Clozapine        |                      |

Figure 16. Representative competitive binding curves with Histamine receptors.








**Table 10.** Neurotransmitter transporters, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included. BB for binding buffer; WB for wash buffer.

| Neurotransmitter transporters  |                             |  |               |                      |
|--|-----------------------------|--|---------------|----------------------|
| Transporter BB: 10 mM HEPES, 135 mM NaCl, 5 mM KCl, 0.8 mM MgCl <sub>2</sub> , 1 mM ETGA, pH 7.4, RT |                             |  |               |                      |
| Transporter WB: Transporter binding buffer, pH 7.4, cold   |                             |  |               |                      |
| Target   | Radioligand                 | $K_d$ Or [ <sup>3</sup> H] for binding in nM (N) | References    | Reference $K_i$ (nM) |
| DAT  | [ <sup>3</sup> H]WIN35428   | 11.1 $\pm$ 1.4 (10)                              | GBR12909      | 3.04 $\pm$ 0.21      |
| NET  | [ <sup>3</sup> H]Nisoxetine | 4.9 $\pm$ 0.7 (4)                                | Desipramine   | 3.16 $\pm$ 0.28      |
| SERT   | [ <sup>3</sup> H]Citalopram | 3.8 $\pm$ 2.1 (4)                                | Amitriptyline | 4.46 $\pm$ 0.55      |

Figure 17. Representative competitive binding curves with neurotransmitter transporters.








**Table 11.** Opioid receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentrations of radioligand used for competition binding assay are usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

| Opioid receptors   |                         |  |                   |                                     |
|--|-------------------------|--|-------------------|-------------------------------------|
| Standard binding buffer: 50 mM Tris HCl, 10 mM MgCl <sub>2</sub> , 0.1 mM EDTA, pH 7.4, RT |                         |  |                   |                                     |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, 4 °C to 8 °C                                 |                         |  |                   |                                     |
| Target   | Radioligand             | $K_d$ or [ <sup>3</sup> H] for binding in nM (N) | References        | Reference $K_i$ (nM)                |
| DOR  | [ <sup>3</sup> H]DADLE  | 1.85 $\pm$ 0.15 (2)                              | Naltrindole       | 0.81 $\pm$ 0.08                     |
| KOR  | [ <sup>3</sup> H]U69593 | 1.07 $\pm$ 0.10 (21)                             | Salvinorin A      | 1.93 $\pm$ 0.45                     |
| MOR  | [ <sup>3</sup> H]DAMGO  | 1.73 $\pm$ 0.14 (6)                              | DAMGO, Morphine   | 2.62 $\pm$ 0.22                     |
| NOP  | [ <sup>3</sup> H]N/OFQ  | 0.74 $\pm$ 0.22 (4)                              | JDTiC<br>SB612111 | 12.05 $\pm$ 1.47<br>6.58 $\pm$ 1.42 |

Figure 18. Representative competitive binding curves with opioid receptors.





**Table 12.** Oxytocin and Vasopressin receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentrations of radioligand used for competition binding assays is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

Oxytocin and Vasopressin receptors


Oxytocin binding buffer: 50 mM HEPES, 10 mM  $MgCl_2$ , pH 7.4



Vasopressin Binding Buffer: 20 mM Tris HCl, 100 mM NaCl, 10 mM  $MgCl_2$ , 0.1 mg/ml bacitracin, 1 mg/ml BSA, pH 7.4, RT

Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold

| Target   | Radioligand          | $K_d$ in nM (N)      | References  | Reference $K_i$ (nM) |
|----------|----------------------|----------------------|-------------|----------------------|
| Oxytocin | [ $^3H$ ]Oxytocin    | $3.10 \pm 0.50$ (30) | Oxytocin    | $6.07 \pm 0.78$      |
| V1a      | [ $^3H$ ]Vasopressin | $2.12 \pm 0.98$ (32) | Vasopressin | $1.40 \pm 0.17$      |
| V1b      | [ $^3H$ ]Vasopressin | $1.33 \pm 0.15$ (9)  | Vasopressin | $1.53 \pm 0.33$      |
| V2       | [ $^3H$ ]Vasopressin | $3.03 \pm 0.48$ (10) | Vasopressin | $4.60 \pm 0.57$      |

Figure 19. Representative competitive binding curves with Oxytocin and Vasopressin receptors.






**Table 13.** Smoothened receptor, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. Historical reference  $K_i$  values from last the 6 months are also included.

| Smoothened receptor   |                              |                         |                  |                              |
|---|------------------------------|-------------------------|------------------|------------------------------|
| SMO binding buffer: 50 mM HEPES, 3 mM MgCl <sub>2</sub> , Protease inhibitors, 0.1mg/ml BSA, pH 7.2, RT |                              |                         |                  |                              |
| SMO wash buffer: PBS, pH 7.2, cold  |                              |                         |                  |                              |
| Target  | Radioligand                  | [ <sup>3</sup> H] in nM | References       | Reference $K_i$ (nM)         |
| SMO   | [ <sup>3</sup> H]Cyclopamine | 3.94 ± 0.67 (3)         | SAG<br>LY2490680 | 12.19 ± 5.01<br>18.06 ± 2.86 |

**Figure 20.** Representative competitive binding curves with SMO receptors.



**Table 14.** Prostanoid receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays.

| Prostanoid receptors   |                       |                         |            |
|--|-----------------------|-------------------------|------------|
| Prostanoid binding buffer: 25 mM Tris HCl, 10 mM MgCl <sub>2</sub> , 1 mM EDTA, pH 7.4, RT |                       |                         |            |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                       |                         |            |
| Target   | Radioligand           | [ <sup>3</sup> H] in nM | References |
| EP3  | [ <sup>3</sup> H]PGE2 | 10                      | EP2        |
| EP4  | [ <sup>3</sup> H]PGE2 | 10                      | EP2        |

**Table 15.** PKC subunits, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays.

| PKC  |                       |                         |            |
|--|-----------------------|-------------------------|------------|
| PKC binding buffer: 50 mM Tris HCl, 1 mM CaCl <sub>2</sub> , 4 mg/ml BSA, 100 µg/ml phosphatidylserine, pH 7.4, RT |                       |                         |            |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                       |                         |            |
| Target   | Radioligand           | [ <sup>3</sup> H] in nM | References |
| PKC $\alpha$   | [ <sup>3</sup> H]PDBU | 3                       | PDBU       |
| PKC $\beta$  | [ <sup>3</sup> H]PDBU | 3                       | PDBU       |
| PKC $\gamma$   | [ <sup>3</sup> H]PDBU | 3                       | PDBU       |
| PKC $\delta$   | [ <sup>3</sup> H]PDBU | 3                       | PDBU       |
| PKC $\epsilon$   | [ <sup>3</sup> H]PDBU | 3                       | PDBU       |


**Table 16.** GABAA receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. Historical reference  $K_i$  values from the last 6 months are also included.


| GABA receptors   |                                |                         |                        |                      |
|--|--------------------------------|-------------------------|------------------------|----------------------|
| GABA/PBR binding buffer: 50 mM Tris Acetate, pH 7.4, RT                                    |                                |                         |                        |                      |
| Benzodiazepine (BZP) binding buffer: 50 mM Tris HCl, 2.5 mM CaCl <sub>2</sub> , pH 7.4, RT |                                |                         |                        |                      |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                                |                         |                        |                      |
| Target   | Radioligand                    | [ <sup>3</sup> H] in nM | References             | Reference $K_i$ (nM) |
| GABA/PBR<br>(rat brain)  | [ <sup>3</sup> H]PK11195       | 1                       | PK11195<br>Ro5-4864    | 27.6 ± 2.3           |
| GABAA<br>(rat brain)   | [ <sup>3</sup> H]Muscimol      | 5.0                     | GABA                   | 241 ± 26             |
| GABAA/BZP<br>(rat brain)   | [ <sup>3</sup> H]Flunitrazepam | 0.5                     | Diazepam<br>Clonazepam | 1.50 ± 0.08          |
| α1β2γ2   | [3H]Flunitrazepam              | 0.5                     | Diazepam               | 1.2                  |
| α2β2γ2   | [3H]Flunitrazepam              | 0.5                     | Diazepam               | 0.8                  |
| α3β2γ2   | [3H]Flunitrazepam              | 0.5                     | Diazepam               | 1.4                  |
| α5β2γ2   | [3H]Flunitrazepam              | 0.5                     | Diazepam               | 0.6                  |

Notes:


Protocol was adapted from Nadler et al., Brain Res Dev Brain Res 97(2): 216-225 (1996)

Figure 21. Representative competitive binding curves with GABA receptors.







**Alpha1 of GABAA Receptor**  
**[<sup>3</sup>H]Flunitrazepam (0.5 nM)**  
**50 mM Tris-acetate**




Alpha2 of GABAA Receptor  
[<sup>3</sup>H]Flunitrazepam (0.5 nM)  
50 mM Tris-acetate



Alpha3 of GABAA  
[<sup>3</sup>H]Flunitrazepam (0.5 nM)  
50 mM Tris-acetate



Alpha5 GABAA Receptor  
[<sup>3</sup>H]Flunitrazepam (0.5 nM)  
50 mM Tris-acetate





**Table 17.** Glutamate receptors radioligands, corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. Membrane pellets for NMDA binding are made with rat brain cortex with a protocol modified from Chiu et al., 1999.

| Glutamate receptors   |                        |                         |            |                               |
|---|------------------------|-------------------------|------------|-------------------------------|
| NMDA binding buffer: 20 mM HEPES, 1 mM EDTA, 100 $\mu$ M L-Glu(Na), 100 $\mu$ M Glycine, pH 7.0, RT |                        |                         |            |                               |
| mGluR binding buffer: 50 mM Tris HCl, 10 mM MgCl <sub>2</sub> , 0.1 mM EDTA, pH 7.4, RT             |                        |                         |            |                               |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold  |                        |                         |            |                               |
| Target  | Radioligand            | [ <sup>3</sup> H] in nM | References | Reference K <sub>i</sub> (nM) |
| NMDA (rat brain)  | [ <sup>3</sup> H]MK801 | 1                       | MK801      | 4.17 $\pm$ 0.74               |
| mGluR1  |                        |                         |            |                               |
| mGluR2  |                        |                         |            |                               |
| mGluR3  |                        |                         |            |                               |
| mGluR4  |                        |                         |            |                               |
| mGluR5  | [ <sup>3</sup> H]MPEP  | 0.9                     | MTEP       | 27                            |
| mGluR6  |                        |                         |            |                               |
| mGluR7  |                        |                         |            |                               |
| mGluR8  |                        |                         |            |                               |

**Notes:**

- 1) The mGluR binding protocol is adapted from Kozikowski et al., J Med Chem 33(6):1561-1571 (1990)
- 2) The NMDA binding protocol is modified from Chiu et al., Alcohol 17:215-221 (1999).
- 3) The mGluR1, 2, 3, 4, 6, 7, 8 binding assays are under development


Figure 22. Representative competitive binding curves for NMDA and mGluR5 receptors.



**Table 18.** Sigma receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are the average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

| Sigma receptors  |                              |                       |             |                      |
|--|------------------------------|-----------------------|-------------|----------------------|
| Sigma Binding Buffer: 50 mM Tris HCl, pH 8.0, RT<br>Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold |                              |                       |             |                      |
| Target   | Radioligand                  | $K_d$ in nM (N)       | References  | Reference $K_i$ (nM) |
| Sigma1<br>(Guinea pig)   | [ <sup>3</sup> H]Pentazocine | 4.45 $\pm$ 0.55 (2)   | Haloperidol | 3.29 $\pm$ 0.50      |
| Sigma2<br>(PC12)   | [ <sup>3</sup> H]DTG         | 10.08 $\pm$ 1.38 (10) | Haloperidol | 24.7 $\pm$ 1.8       |

Figure 23. Representative competitive binding curves with Sigma receptors.




**Table 19.** Angiotensin II receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed.

| Angiotensin receptors  |                                 |                         |             |                               |
|--|---------------------------------|-------------------------|-------------|-------------------------------|
| Angiotensin Binding Buffer: 50 mM Tris HCl, 150 mM NaCl, 5 mM MgCl <sub>2</sub> , 0.5 mg/ml BSA, 100 mM Bacitracin, pH 7.4, RT |                                 |                         |             |                               |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                                 |                         |             |                               |
| Target   | Radioligand                     | [ <sup>3</sup> H] in nM | References  | Reference K <sub>i</sub> (nM) |
| AT1  | [ <sup>3</sup> H]Angiotensin II | 0.1                     | Candesartan |                               |
| AT2  | [ <sup>3</sup> H]Angiotensin II | 0.1                     | PD123319    |                               |

**Table 20.** Neurotensin receptors, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed. The  $K_d$  values listed in this table are an average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2011 to 2012.

| Neurotensin receptors                                |                             |   |             |     |
|--|-----------------------------|---|-------------|-----|
| Neurotensin binding buffer: 50 mM Tris HCl, 0.2% BSA |                             |   |             |     |
| Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold   |                             |   |             |     |
| Target   | Radioligand                 | $K_d$ or [ $^3\text{H}$ ] in nM for binding (N) | References  |     |
| NTS1   | [ $^3\text{H}$ ]Neurotensin | $6.33 \pm 2.03$ (3)                             | Neurotensin | 1.8 |
| NTS2   | [ $^3\text{H}$ ]Neurotensin | 2 nM  | Neurotensin |     |


Figure 24. Representative binding curve for NTS1 receptor.



**Table 21.** VMAT2 transporter, radioligand and corresponding concentration, reference compound, and buffers for primary and secondary radioligand binding assays.

| VMAT2 transporter  |                                 |   |            |                               |
|--|---------------------------------|---|------------|-------------------------------|
| VMAT2 binding buffer: 50 mM Tris HCl, 0.2% BSA<br>Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold |                                 |   |            |                               |
| Target   | Radioligand                     | K <sub>d</sub> or [ <sup>3</sup> H] in nM for binding (N) | References | Reference K <sub>i</sub> (nM) |
| VMAT2  | [ <sup>3</sup> H]Tetrabenzazine | 1.5   | Reserpine  | 280                           |
|  |                                 |   |            |                               |


Figure 25. Representative binding curve for VMAT2 transporters.



**Table 22.** Calcium and sodium channels, radioligands and corresponding concentrations, reference compounds, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value, or as listed. The  $K_d$  values listed in this table are an average  $\pm$  SEM from a minimum of 3 saturation binding assays before 2012.

| Calcium and Sodium channels   |                                |  |             |                      |
|---|--------------------------------|--|-------------|----------------------|
| Calcium channel Binding Buffer: 50 mM Tris HCl, 50 mM NaCl, 1 mM CaCl <sub>2</sub> , pH 7.4, RT<br>Sodium channel Binding Buffer: 50 mM HEPES, 130 mM Choline Cl, 5.4 mM KCl, 0.8 mM MgSO <sub>4</sub> , 5.5 mM Glucose, 1 $\mu$ M tetrodotoxin, 1 mg/ml BSA, 30 $\mu$ g/well scorpion verom, pH 7.4, 37°C.<br>Standard Wash Buffer: 50 mM Tris HCl, pH 7.4, cold |                                |  |             |                      |
| Target  | Radioligand                    | $K_d$ or [ <sup>3</sup> H] in nM for binding (N) | References  | Reference $K_i$ (nM) |
| Ca <sup>2+</sup> channel  | [ <sup>3</sup> H]Nitrendipine  | 4.77 $\pm$ 1.75 (3)                              | Nifendipine |                      |
| Na <sup>+</sup> channel   | [ <sup>3</sup> H]Batrachotoxin |  | Veratridine |                      |


Figure 26. Representative figure for  $\text{Ca}^{2+}$  binding assay.



**Table 23.** HERG potassium channel, radioligand and corresponding concentration, reference compound, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value, or as listed. The  $K_d$  values listed in this table are an average  $\pm$  SEM from a minimum of 3 saturation binding assays from 2009 to 2012. Historical reference  $K_i$  values from the last 6 months are also included.

| hERG channel   |                             |                      |                         |                                 |
|--|-----------------------------|----------------------|-------------------------|---------------------------------|
| hERG binding buffer: 10 mM HEPES, 135 mM NaCl, 5 mM KCl, 0.8 mM MgCl <sub>2</sub> , 1 mM EGTA, 1 mg/ml BSA, pH 7.4, RT |                             |                      |                         |                                 |
| hERG wash buffer: hERG binding buffer, pH 7.4, cold  |                             |                      |                         |                                 |
| Target   | Radioligand                 | $K_d$ in nM (N)      | References              | Reference $K_i$ (nM)            |
| hERG   | [ <sup>3</sup> H]Dofetilide | 2.93 $\pm$ 0.55 (11) | Dofetilide<br>Cisapride | 4.1 $\pm$ 0.4<br>38.9 $\pm$ 3.9 |

Figure 27. Representative binding curve with hERG potassium channel.




**Table 24.** Imidazoline receptor, radioligand and corresponding concentration, reference compound, and buffers for primary and secondary radioligand binding assays. The concentration of radioligand used for competition binding assay is usually at or near the  $K_d$  value or as listed.


| Imidazoline  |                            |                                     |             |
|--|----------------------------|-------------------------------------|-------------|
| Imidazoline binding buffer: 5 mM Tris HCl, 5 mM HEPES, 0.5 mM EGTA, 0.5 mM EDTA, 0.5 mM MgCl <sub>2</sub> , pH 8.0, RT |                            |                                     |             |
| Standard wash buffer: 50 mM Tris HCl, pH 7.4, cold   |                            |                                     |             |
| Target   | Radioligand                | [ <sup>3</sup> H] in nM for binding | References  |
| Imidazoline 1  | [ <sup>3</sup> H]Clonidine | 0.1 nM                              | Naphazoline |
|  |                            |                                     |             |

## Section 2: Functional assays


**2.1. Drug plate preparations for functional assays:** Drug plates for functional assays are either manually prepared or made by the STAR robotic system. The following drug plate maps are the ones most commonly used for functional assays.



**Figure 28.** 384-well drug plate map #1 for primary screening: Singlet format (first addition). A total of 320 different drugs (test compounds) can be plated in singlet format in a 384-well drug plate. Columns 1, 2, 23, and 24 are usually used for positive and negative controls



**Figure 29.** 384-well drug plate map #2 for primary screening: quadrant format (first addition). A total of 80 different drugs can be plated in quadrant format (q1, q2, q3, and q4) in a 384-well drug plate. Columns 1, 2, 23, and 24 are usually used for positive and negative controls.




**Figure 30.** 384-well drug plate map #3: second or third addition in primary screening assays (second or third addition). Except for Column 1 (which serves as a negative control with assay buffer for all additions), all the other wells receive EC<sub>20</sub> or EC<sub>80</sub> of a reference agonist to determine allosteric modulator or antagonist activity.

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | [X] (M) |       |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|---------|-------|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    | 0       |       |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-12 |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-11 |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-11 |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-10 |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-10 |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-09 |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-09 |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-08 |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-08 |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-07 |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-07 |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-06 |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-06 |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-05 |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-05 |

|   |        |        |        |        |        |        |        |        |
|---|--------|--------|--------|--------|--------|--------|--------|--------|
| 8 | Drug A | Drug B | Drug C | Drug D | Drug E | Drug F | Drug G | Drug H |
|---|--------|--------|--------|--------|--------|--------|--------|--------|

|   |         |         |         |         |         |         |
|---|---------|---------|---------|---------|---------|---------|
| 6 | Drug A' | Drug B' | Drug C' | Drug D' | Drug E' | Drug F' |
|---|---------|---------|---------|---------|---------|---------|

**Figure 31.** 384-well drug plate map #4: agonist or antagonist concentration-responses for secondary screening (first addition): Eight serial drug dilutions are made in triplicate (Drug A to H) or six drug serial dilutions in quadruplicate (Drug A' to F') (indicated below the plate template) from high to low concentrations (final concentrations are indicated to the right of the plate template). Either format can be used; one of the compounds is a positive control with a known agonist and/or antagonist, such as acetylcholine and/or atropine for muscarinic receptors.



**Figure 32.** 384-well drug plate map #5: EC<sub>50</sub> or EC<sub>80</sub> of a reference agonist as second addition for antagonist activity (second addition). The even numbered wells in row A serve as negative controls with assay buffer only, while the other wells on the plate contain EC<sub>50</sub> – EC<sub>80</sub> concentrations of a reference agonist.

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |

| X | Buffer | 30 nM | 100 nM | 300 nM | 1 $\mu$ M | 3 $\mu$ M | 10 $\mu$ M | 30 $\mu$ M |
|---|--------|-------|--------|--------|-----------|-----------|------------|------------|
|---|--------|-------|--------|--------|-----------|-----------|------------|------------|

**Figure 33.** 384-well drug plate map #6 for Schild plot analysis (first addition): Antagonist or allosteric modulator is made in 7 concentrations as indicated at the bottom of the drug plate template. Each concentration (including buffer control) has three columns.

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | [X] (M) |       |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|---------|-------|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    | 0       |       |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-12 |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-11 |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-11 |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-10 |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-10 |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-09 |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-09 |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-08 |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-08 |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-07 |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-07 |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-06 |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-06 |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 1E-05 |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |         | 3E-05 |

**Figure 34.** 384-well drug plate map #7 for Schild plot analysis (second addition). Reference agonist is made in the following serial dilutions (as indicated at the right of the drug plate template). Each concentration (including buffer control) has one complete row.

## 2.2. General procedures for PDSP functional assays:

In general, PDSP functional assays are carried out in two steps, primary and secondary assays (unless otherwise stated). In primary screening assays, compounds are tested in triplicate or quadruplicate at a final concentration of 10  $\mu\text{M}$  or 10  $\mu\text{g/ml}$  (for crude extracts), or at specific concentrations upon request for agonist and antagonist activities (see 384-well drug plate maps #1, #2, and #3 for detailed setup). Results are normalized and transformed to percentage values. For agonists, the reference agonist activity at 10  $\mu\text{M}$  is set as 100% and basal activity with buffer as 0%. **For orphan GPCRs without a known agonist as a reference, activity is expressed as percentage value of basal (with buffer).** For antagonists, the basal activity with buffer only is set as 100% inhibition and the activity of the  $\text{EC}_{80}$  of the reference agonist as 0% inhibition. For allosteric potentiators, the activity of the  $\text{EC}_{20}$  of a reference agonist is set as 0% potentiation. Compounds with a minimum of 30% agonist activity, or 50% antagonist activity, or 30% potentiation above control are selected for secondary screening using concentration-response assays.

In secondary functional assays, potential agonist hits are subjected to full concentration-response studies to determine their efficacy and potency (see 384-well drug plate map #4). Potential antagonist hits are subjected to full concentration-response studies to determine  $\text{IC}_{50}$  values against  $\text{EC}_{50}$  to  $\text{EC}_{80}$  concentrations of a reference agonist (see 384-well drug plate map #5 for detailed setup). Potential antagonist hits and positive allosteric modulators are further characterized, if necessary, by Schild plot analysis to obtain  $\text{pK}_B$  and  $\text{pA}_2$  values (see 384-well drug plate maps #6 and #7 for detailed setup). Schild plot analysis is designed to obtain full agonist concentration-response curves in the absence and presence of increasing concentrations of potential antagonists or allosteric modulators. For those antagonists that shift agonist concentration-response curves to the right and also reduce  $\text{E}_{\text{max}}$ , a modified Lew and Angus analysis is applied to obtain  $\text{pK}_B$  and  $\text{pA}_2$  values. For allosteric modulators, secondary concentration-response curves are also analyzed using the allosteric operational model to obtain modulation parameters (such as  $\alpha$  and  $\beta$ ). For studies of ligand functional selectivity and bias among multiple signaling pathways, concentration-response results are analyzed using the Black

and Leff operational model to quantify bias. Detailed procedures for these secondary assays are given in subsequent data analysis sections.

Secondary functional assays provide the following parameters for the PDSP database: maximum activation or inhibition after normalization (in percentage values), concentration range used in the assay (lowest and highest concentrations), Hill slope of the concentration-response curve, and corresponding potency.

### 2.3. Data analysis for functional assays

Unless otherwise stated, all functional results are analyzed using GraphPad Prism v5.0 using its built-in functions.

**2.3.1: Agonist activity:** Functional assay results are plotted against concentrations and analyzed with the following 4-parameter built-in agonist dose-response function in GraphPad Prism v5.0.

$$Response = Bottom + \frac{(Top - Bottom)}{1 + 10^{(LogEC_{50} - X)n}}$$

In which, **Top** and **Bottom** are maximum response ( $E_{max}$ ) and basal activity, respectively; **X** is the agonist concentration and **n** is the Hill slope;  $EC_{50}$  is the concentration that generates a 50% level of activity.

**2.3.2: Antagonist activity:** To determine antagonist activity, agonist responses are measured at a fixed  $EC_{80}$  concentration of the reference agonist in the presence of a series of dilutions of antagonist. Results are fitted with the following inhibitory concentration-response function to determine the  $IC_{50}$  value.

$$Response = Bottom + \frac{(Top - Bottom)}{1 + 10^{(LogIC_{50} - X)n}}$$

The  $IC_{50}$  is then converted to  $K_i$  using the Cheng-Prusoff equation:

$$K_i = \frac{IC_{50}}{1 + \frac{L}{EC_{50}}}$$

in which  $K_i$  is the ligand binding affinity determined from an antagonist concentration-response assay;  $IC_{50}$  is the antagonist concentration at which point 50% inhibition is reached;  $L$  is the reference agonist concentration used in the assay (usually the  $EC_{50}$  to  $EC_{80}$  concentration of the reference agonist);  $EC_{50}$  is the predetermined potency of the reference agonist.

**2.3.3. Schild plot analysis (see [Figure 36](#) for examples):** The PDSP has adopted a nonlinear regression analysis with modified Lew and Angus method to estimate antagonist potency  $pA_2$  values. For Schild plot analysis, agonist concentration-response studies are designed and performed in the absence and presence of 7 concentrations of antagonist in a 384-well plate (see **Figures 33 and 34**). Each agonist concentration-response curve is based on triplicate values. Eight agonist concentration-response curves share the same X-axis and each set of measured responses is arranged in a Y column under antagonist concentrations (M) as a title. Results are analyzed using the agonist concentration-response function as above to determine Bottom and Top values for each concentration-response curve. To normalize using Prism's built-in normalization function, the data set is transformed into percentage values with the shared **bottom** as 0% and **Top** of the reference agonist under control conditions as 100%. An equiactive agonist concentration (such as an agonist concentration to generate 20 – 50% activity) is selected for all or most curves if possible, and corresponding agonist concentration-response curves in the absence and presence of increasing concentrations of antagonist are obtained. Prism 5 has a built-in feature to provide equiactive agonist concentrations if the radio button for "Interpolate unknowns from standard curve" is selected. The feature is at the bottom of the Analyze>Fit tab. The equiactive agonist concentrations are plotted in  $-\log$  format (i.e.  $pEC_{20}$  or  $pEC_{50}$ ) against corresponding antagonist (B) concentrations. The results are fitted using following equation to obtain  $\text{Log}K_B$  and  $n$  values.

$$Y = -\text{Log}(X^n + 10^{\text{Log}K_B}) - \text{Log}C$$

In which Y is the equiactive agonist concentration in  $-\log$  format (such as  $pEC_{25}$  or other equal-active agonist concentration values) at corresponding antagonist concentration (X); n is the Schild slope;  $K_B$  is the apparent binding affinity of testing antagonist. The  $pA_2$  value could be calculated from  $pK_B/n$ . If  $n = 1$  or not significantly different from 1,  $pA_2 = pK_B$ , the testing antagonist is competitive against agonist; otherwise, i.e., if  $n \neq 1$ , the antagonist may not be competitive against agonist.

**2.3.4. Quantifying bias and functional selectivity analysis (see [Figures 39 and 40](#), [Tables 26 and 27](#) for examples).** To quantify functional selectivity and to calculate a bias factor, we analyze functional concentration-response results using the Black and Leff operation model as outlined by Dr. Terry Kenakin et al., 2012 (ACS Chemical Neuroscience, 3:193-203, 2012).

$$Response = Basal + \frac{(E_m - Basal)[A]^n \tau^n}{[A]^n \tau^n + ([A] + K_A)^n}$$

Functional dose-response results are analyzed using the above equation in GraphPad Prism V5.0, in which  $E_m$  is the maximal possible response of the system and **basal** is the response in the absence of test drugs (i.e. buffer only). The  $E_m$  and **Basal** are usually 100 and 0 if responses are normalized to a percentage of the activity of a reference compound (100%).  $K_A$  is the equilibrium dissociation constant of the agonist (**A**),  $\tau$  is the operational agonist efficacy of the agonist (**A**) and is defined as  $R_T/K_E$  (where  $R_T$  is the receptor density and  $K_E$  is the intrinsic efficacy of the agonist (**A**) in a particular signaling pathway),  $n$  is the transducer slope for the function between agonist occupancy and measured responses. The fitting parameters  $E_m$  and  $n$  are cell-specific and are shared by all agonists that are being tested for the same pathway.  $\log(\tau/K_A)$  is defined as “transduction coefficient” for a particular agonist at a measured signaling pathway.

To quantify ligand bias (see the following equations and steps), drug activity is measured in two or more pathways for a group of ligands (e.g., agonists) to determine their corresponding  $\log(\tau/K_A)$  values. For each pathway,  $\Delta\log(\tau/K_A)$  values are calculated by subtracting the  $\log(\tau/K_A)$  value of the reference agonist (usually endogenous agonist). The same reference agonist should

be used for the different pathways. For each agonist,  $\Delta\Delta\log(\tau/K_A)$  values are calculated by subtracting  $\Delta\log(\tau/K_A)$  for pathway I from the corresponding  $\Delta\log(\tau/K_A)$  for pathway II. Bias is quantified as  $10^{\Delta\Delta\log(\tau/K_A)}$ . If the value of  $10^{\Delta\Delta\log(\tau/K_A)}$  is 1, the agonist is unbiased; if the value is larger than 1, it is biased towards pathway I and if it is less than 1, it is biased towards pathway II.

**Step 1:** Concentration-response curves are generated for 2 or more pathways to determine the 'transduction coefficient'  $\text{Log}(\tau/K_A)$  value for each ligand at each pathway;

**Step 2:** The transduction coefficient difference  $\Delta\text{Log}(\tau/K_A)$  between different ligands in the same pathway is calculated by subtracting transduction coefficient of the test sample from that of the reference compound, usually the endogenous ligand, such as 5-HT for serotonin receptors;

$$\Delta\text{Log}\left(\frac{\tau}{K_A}\right) = \text{Log}\left(\frac{\tau}{K_A}\right) \text{ of sample} - \text{Log}\left(\frac{\tau}{K_A}\right) \text{ of reference}$$

**Step 3:** The difference between transduction coefficients  $\Delta\Delta\text{Log}(\tau/K_A)$  is calculated between different pathways for the same test ligand by subtracting the  $\Delta\text{Log}(\tau/K_A)$  for one pathway from that of the other pathway;

$$\Delta\Delta\text{Log}\left(\frac{\tau}{K_A}\right) = \Delta\text{Log}\left(\frac{\tau_1}{K_{A1}}\right) - \Delta\text{Log}\left(\frac{\tau_2}{K_{A2}}\right)$$

**Step 4:** The bias factor is calculated as below

$$\text{Bias} = 10^{\Delta\Delta\log\left(\frac{\tau}{K_A}\right)}$$

**Step 5:** to determine if this bias value is statistically significant, a statistical analysis is done with  $\Delta\Delta\text{Log}(\tau/K_A)$  values. In addition to the comprehensive methods to estimate SEM values outlined by Dr. Terry Kenakin (2012), an alternative (and simpler) method to estimate SEM values is given below, as suggested by Dr. Arthur Christopoulos and his colleagues (Gregory et al., 2012, JBC 287:37066-37077).

To calculate SEM for  $\Delta\log(\tau/K_A)$ :

$$SEM = \sqrt{\text{Sample sem}^2 + \text{Reference sem}^2}$$

To calculate SEM for  $\Delta\Delta\log(\tau/K_A)$ :

$$SEM = \sqrt{\text{pathway I sem}^2 + \text{pathway II sem}^2}$$

**2.3.5. Allosteric operational analysis (see [Figure 41](#) for an example).** For ligands with allosteric modulator activity, we apply the allosteric operational model to analyze functional results. Functional assays are carried out in the same way as for Schild plot analysis (**Section 2.3.3**), in which orthosteric agonist concentration-response curves are measured in the absence and presence of increasing concentrations of a potential allosteric modulator. Results are then analyzed with the following equation as described by Leach et al., 2007 (Trends in Pharmacological Sciences 28: 382-389).

$$Effect = Basal + \frac{(E_{Max} - Basal) (\tau_A[A](K_B + \alpha\beta[B]) + \tau_A[B]K_A)^n}{([A]K_B + K_A K_B + K_A[B] + \alpha[A][B])^n + (\tau_A[A](K_B + \alpha\beta[B]) + \tau_B[B]K_A)^n}$$

In which **Effect** is the measured functional readout in the presence of an orthosteric agonist [**A**] and an allosteric modulator [**B**]; **E<sub>max</sub>** and **Basal** are the maximal possible system activity and basal activity under control conditions; **K<sub>A</sub>** and **K<sub>B</sub>** are the equilibrium binding affinities for the orthosteric agonist (**A**) and the allosteric modulator (**B**), respectively; **α** and **β** are the allosteric effects on ligand binding (mutual effect between **A** and **B**) and agonist efficacy (with **α** = 1 for neutral cooperativity; **α** >1 for positive cooperativity; **α** <1 for negative cooperativity); **τ<sub>A</sub>** and **τ<sub>B</sub>** are the capacity of the orthosteric agonist (**A**) and the allosteric ligand (**B**) to exhibit agonism, respectively; **n** is a fitting slope factor.

While doing non-linear least-squares regression curve-fitting in Prism v5.0, the  $E_{\max}$  should be a shared value for all assays carried out under the same conditions, and should be set to a fixed value which is the maximal system efficacy. The  $K_A$  should be the equilibrium binding affinity and can be estimated by radioligand binding assay. If the test allosteric modulator itself has no agonist activity, then  $\tau_B = 0$ .

### 2.3.6. List of functional assays carried out routinely by the PDSP:

- Calcium mobilization assay with FLIPR<sup>TETRA</sup> for  $G_q$  coupled GPCRs.
- Intracellular inositol phosphate accumulation assay.
- Split luciferase based biosensor cAMP<sup>®</sup> assay for  $G_i$  or  $G_s$  coupled GPCRs.
- GPCR Tango assay for G-protein independent  $\beta$ -arrestin translocation.
- Thallium flux (FluxOR) assay for hERG potassium channel function.
- PatchXpress automated patch clamp assay for hERG potassium channels.
- Neurotransmitter transporter assays for DAT, NET, and SERT transporters.
- Multidrug resistance transporter-1 (MDR-1) assay
- Enzyme activity assays
  - HDAC assay
  - MAO-A and -B assays
  - PKC assay
  - CHK2 assay
- Nicotinic acetylcholine receptors (nAChRs) activity assay ( $^{86}\text{Rb}^+$  efflux)

## 2.4. Functional assays for G<sub>q</sub> coupled GPCRs

### 2.4.1. Calcium mobilization assays (with FLIPR<sup>TETRA</sup>)

**Main equipment:** FLIPR<sup>TETRA</sup> from Molecular Devices (Sunnyvale, CA)

**Main reagents:** Fluo-4 Direct<sup>®</sup> from Invitrogen (Carlsbad, CA)

**FLIPR Drug buffer:** 20 mM HEPES, 1x HBSS, 2.5 mM Probenecid, pH 7.40, room temperature

**2.4.1.1. Cell culture:** Cells, either stably expressing target receptors or transiently transfected with target receptor DNA (see below for transfection protocol) overnight, are plated into Poly-L-Lysine (PLL) coated 384-well black clear bottom cell culture plates with DMEM supplemented with 1% dialyzed FBS (dFBS) and at density of 15 – 20,000 cells in 40 µl per well. The plates are cultured for minimum of 6 hour or overnight before assays.

**2.4.1.2. Calcium Precipitation transfection:** The PDSP usually uses HEK293 T cells for transient transfections. HEK293 T cells are subcultured into either 10-cm dishes (3 million cells per dish) or 15-cm dishes (8 million cells per dish) and are incubated overnight. For each 10-cm dish of HEK293 T cells, 10 µg receptor of DNA construct are mixed in 440 µl distilled water and with 60 µl of 2.5 M CaCl<sub>2</sub>, then add dropwise into 500 µl 2x HBS solution (50 mM HEPES, 280 mM NaCl, 10 mM KCl, 1.5 mM Na<sub>2</sub>HPO<sub>4</sub>, pH 7.00) while shaking. The mixture is incubated at room temperature for 10 min, and then added to cells dropwise, which are then incubated overnight. For transfections in 15-cm dishes, reagents and DNA amounts are increased by 2.5 fold per dish.

**2.4.1.3. Calcium mobilization assays with FLIPR<sup>TETRA</sup>.** On the day of the assay, medium is removed and cells are loaded with 20 µl/well of 1x Fluo-4 Direct Calcium dye (prepared in FLIPR drug buffer). Plates are incubated for 60 min at 37°C, followed by a 10 min incubation at room temperature in the dark, and then loaded in the FLIPR. The FLIPR is programmed to take 10 readings (1 read per second) first as a baseline before addition of 10 µl of 3x drug solutions. The fluorescence intensity is recorded for 2 minutes after drug addition (first addition, see 384-well drug plate maps #1 and #2 for detailed setup) for agonist activity. To measure antagonist activity, drug stocks are prepared at 4x of the final concentration, and added as above for

potential effects on basal levels for 2 minutes first, followed by a 5 to 10 min incubation before addition of 10  $\mu$ l of 4x of a reference agonist at a final concentration equivalent to the  $EC_{80}$  (see 384-well drug plate map #3) to measure remaining agonist activity. To detect positive allosteric modulator activity, the second addition is at 4x of reference agonist to give a final concentration equivalent to the  $EC_{20}$  (usually the endogenous agonist, see 384-well drug plate map #3) to determine effects on agonist activity. The  $EC_{20}$  or  $EC_{80}$  concentrations are predetermined separately using the same batch of transfected cells. Before the second and/or third addition, the FLIPR is programmed to wash tips first with 10% EtOH, and then with distilled water while recording the fluorescence intensity. Incubation times between each addition can be adjusted accordingly to accommodate any special preincubation requirements.

**2.4.1.4. Schild plot analysis:** To further examine antagonist activity at a particular receptor, we also carry out agonist concentration-response curves in the absence and presence of 7 concentrations of selected PDSP compounds. In brief, cell plates are prepared in the same way as for primary or secondary functional assays using PLL coated 384-well black clear bottom cell culture plates with DMEM + 1% dFBS for a minimum of 6 hours or, more usually, overnight. Selected antagonist drug plates are prepared as shown above (384-well drug map #6) and are used for the first addition to determine if the test drug has any effect on basal levels. A reference agonist drug plate is prepared as shown above (384-well drug map #7) and is used for the second addition to measure effects of antagonists on the remaining agonist activity.

**2.4.1.5. FLIPR data processing:** In Calcium mobilization assays done with the FLIPR<sup>TETRA</sup>, every well of a 384-well plate has its own basal level, which is defined as the average value of the 10 readings before corresponding drug addition. We measure maximal fluorescence intensity readings (RFU, Relative Fluorescence Unit) within a minute after drug addition as agonist activity, and values are exported as “fold of basal” using FLIPR’s ScreenWorks® built-in batch export function. Normalization, when needed, is carried out using Prism’s built-in normalization function, in which basal levels are transformed to 0% and the  $E_{max}$  of the reference agonist in the same assay plate is transformed to 100%.

## 2.4.2. Intracellular inositol phosphate accumulation for Gq coupled GPCRs

**Receptors:** mGluR<sub>1</sub> and mGluR<sub>5</sub>

**Assay buffer:** LiCl-containing Locke's buffer (156 mM NaCl, 5.6 mM KCl, 3.6 mM NaHCO<sub>3</sub>, 1 mM MgCl<sub>2</sub>, 1.3 mM CaCl<sub>2</sub>, 5.6 mM glucose, and 20 mM HEPES, 20 mM LiCl, pH 7.4).

The protocol is adapted from: Emery AC, DiRaddo JO, Miller E, Hathaway HA, Pshenichin S, Takoudjou GR, Grajkowska E, Yasuda RP, Wolfe BB, Wroblewski JT. Ligand bias at metabotropic glutamate 1a receptors: molecular determinants that distinguish b-arrestin-mediated from G protein-mediated signaling. *Mol Pharmacol* 82:291-301 (2012).

**2.4.2.1. General experimental procedure:** Compounds are tested using cell lines stably expressing G<sub>q</sub> protein-coupled mGluR<sub>1</sub> and mGluR<sub>5</sub> receptors in functional assays. Cells are seeded and cultured in 96-well plates to confluency (3 to 4 days), and are incubated overnight in glutamine-free culture medium supplemented with 0.625 Ci/well of *myo*-[<sup>3</sup>H]-inositol (NEN) to label cell membrane phosphoinositides. After two washes with 0.1 ml of Locke's buffer, incubations with receptor ligands are carried out for 45 min at 37°C in Locke's buffer containing 20 mM LiCl to block inositol phosphate degradation. The reaction is terminated by aspiration of media and [<sup>3</sup>H]inositol phosphates are extracted in 60 µl of 10 mM formic acid for 30 min. Samples (40 µl) are transferred to opaque-welled plates and incubated with 60 µl of poly lysine-coated yttrium scintillation proximity assay (SPA) beads (GE Healthcare) at room temperature for 1 hour with vigorous shaking. After an additional 8 hours of incubation with SPA beads, inositol phosphates are detected by scintillation counting. All of these studies are performed in the absence and in the presence of agonists appropriate for the different mGluRs, and used at concentrations equivalent to their EC<sub>50</sub> values: 15 µM and 5 µM glutamate is used for mGluR1 and mGluR5, respectively.

**2.4.2.2. Primary screening assays – Single concentration and data analysis:** Each new compound is tested on selected receptors at a single concentration (10 µM) for activity as an agonist or an antagonist. Testing for antagonism is performed in presence of the EC<sub>50</sub>

concentration of a typical agonist (as described above). Each compound is tested in duplicate in two separate experiments performed on different cell passages. In addition to the tested compounds, each 96-well plate contains points for determination of basal activity, maximal agonist stimulation, agonist EC<sub>50</sub> concentrations (i.e., concentration-response isotherm), and the IC<sub>50</sub> concentration of a known antagonist for purposes of positive control and for activity calculations. The reported results for each compound are calculated for agonists as the percent of maximal activity (as obtained with maximal agonist concentrations) and for antagonist as the percent of inhibition of receptor activity (in presence of an EC<sub>50</sub> concentration of the agonist). Results are expressed as means ± SEM from four replicates.

#### **2.4.2.3. Secondary screening assays: Concentration-response curves and data analysis.**

Compounds determined to be active as agonists or antagonists may be tested for their potency in concentration-response experiments. Six-point concentration-responses curves are performed in duplicate twice on two separate passages of cells (sometimes a third curve may be needed if in the first experiment the range of concentrations used is outside of the active range). For antagonists, these curves are performed in the presence of the EC<sub>50</sub> concentration of the agonist. For each compound, the results from four replicates are averaged and then either EC<sub>50</sub> or IC<sub>50</sub> values are calculated by non-linear regression using the 4-parameter logistic equation as indicated in **Section 2.3**. Results are reported as EC<sub>50</sub> or IC<sub>50</sub> values for each tested compound (and receptor) and also include the EC<sub>50</sub> or IC<sub>50</sub> values of a known agonist or antagonist for comparison.


**2.4.2.4.** The following table lists GPCRs for which the PDSP has functional assays to measure G<sub>q</sub> protein activation via Calcium mobilization using the FLIPR<sup>TETRA</sup> or intracellular inositol phosphate accumulation. Representative results were analyzed using Prism v5.0.



**Table 25.** List of GPCRs for which the PDSP has assays for G<sub>q</sub> protein mediated Calcium mobilization or inositol phosphate (IP) accumulation and their pharmacological parameters. The PDSP will also develop Calcium mobilization assays for other GPCRs upon request and approval. Assays were carried out according to above procedures and results were analyzed in Prism. Representative results are from single assays done in either triplicate or quadruplicate.



| Receptor    | Cell line       | Reference                       | Emax (fold) | pEC50 | Hill slope |
|-------------|-----------------|---------------------------------|-------------|-------|------------|
| M1          | CHO             | Acetylcholine                   | 3.4         | 8.68  | 1.59       |
| M3          | CHO             | Acetylcholine                   | 2.8         | 8.91  | 0.95       |
| M5          | CHO             | Acetylcholine                   | 3.4         | 8.25  | 1.03       |
| M1D         | CHO             | CNO                             | 1.5         | 7.68  | 0.72       |
| M3D         | Flp-In CHO      | CNO                             | 2.7         | 7.45  | 1.01       |
| M5D         | CHO             | CNO                             | 2.1         | 7.44  | 0.81       |
| 5-HT2A      | Flp-In HEK      | 5-HT                            | 2.4         | 8.82  | 1.19       |
| 5-HT2B      | Flp-In HEK      | 5-HT                            | 4.8         | 8.85  | 1.25       |
| 5-HT2C INI  | Flp-In HEK      | 5-HT                            | 3.3         | 10.11 | 1.92       |
| 5-HT2C VNV  | Flp-In HEK      | 5-HT                            | 4.6         | 9.17  | 2.32       |
| 5-HT2C VSV  | Flp-In HEK      | 5-HT                            | 6.0         | 9.06  | 2.12       |
| α1A         | Fibroblast, rat | Norepinephrine                  | 3.0         | 8.60  | 1.34       |
| α1B         | Fibroblast, rat | Norepinephrine                  | 1.6         | 8.72  | 0.92       |
| α1D         | Fibroblast, rat | Norepinephrine                  | 4.1         | 8.61  | 1.14       |
| AT1A        | HEK             | Angiotensin II                  | 1.8         | 9.41  | 0.69       |
| BB1         | 3T3 Balb-C      | Bombesin                        | 5.6         | 7.82  | 1.36       |
| BB1         | 3T3 Balb-C      | BIM 187                         | 5.8         | 7.21  | 1.73       |
| BB1         | 3T3 Balb-C      | Neuromedin B                    | 5.5         | 8.40  | 1.40       |
| BB2         | 3T3 Balb-C      | Bombesin                        | 4.4         | 8.76  | 0.67       |
| BB2         | 3T3 Balb-C      | BIM 187                         | 4.7         | 8.49  | 0.95       |
| BB2         | 3T3 Balb-C      | Neuromedin B                    | 4.6         | 6.58  | 1.29       |
| BB2         | 3T3 Balb-C      | Gastrin-releasing peptide (GRP) | 4.9         | 8.49  | 0.70       |
| BB3         | 3T3 Balb-C      | PDSP compound                   | 4.2         | 5.87  | 4.03       |
| BB2 (mouse) | 3T3 Balb-C      | Bombesin                        | 5.3         | 9.06  | 0.82       |
| BB2 (mouse) | 3T3 Balb-C      | BIM 187                         | 5.6         | 8.74  | 1.23       |
| BB2 (mouse) | 3T3 Balb-C      | Neuromedin B                    | 5.4         | 6.74  | 2.61       |
| BB2 (mouse) | 3T3 Balb-C      | Gastrin-releasing peptide (GRP) | 5.5         | 8.80  | 0.91       |
| B2          | HEK T           | Bradykinin                      | 1.7         | 9.57  | 0.88       |
| CCK2        | HEK T           | Gastrin                         | 3.1         | 9.83  | 1.29       |
| Ghrelin     | HEK             | L-692,585                       | 1.7         | 9.05  | 0.80       |
| H1          | HEK             | Histamine                       | 2.4         | 7.77  | 1.22       |



| Receptor    | Cell line   | Reference      | Emax (fold) | pEC50 | Hill slope |
|-------------|-------------|----------------|-------------|-------|------------|
| H2          | HEK T       | Histamine      | 1.5         | 6.24  | 1.42       |
| mGluR1 (IP) | CHO         | Glutamate      | 6.4         | 4.66  | 1.20       |
| mGluR5 (IP) | CHO         | Glutamate      | 7.1         | 5.15  | 0.94       |
| MC3         |             | $\alpha$ -MSH  | 3.2         | 8.37  | 1.00       |
| MC4         |             | $\alpha$ -MSH  | 6.2         | 6.92  | 1.00       |
| NK1         | HEK         | Substance P    | 2.7         | 7.65  | 0.77       |
| NK2         | HEK         | Neurokinin A   | 4.6         | 7.53  | 0.72       |
| NK2         | HEK         | Neurokinin B   | 4.4         | 6.90  | 0.82       |
| NK3         | HEK         | Neurokinin B   | 4.0         | 8.41  | 0.83       |
| NTS1        | HEK         | Neurotensin    | 3.5         | 8.00  | 0.77       |
| P2Y1        | 1321N1      | ADP            | 3.1         | 6.64  | 1.33       |
| P2Y2        | 1321N1      | ATP            | 2.1         | 6.87  | 1.32       |
| P2Y2        | 1321N1      | UTP            | 1.9         | 8.03  | 1.48       |
| P2Y4        | 1321N1      | UTP            | 2.5         | 7.65  | 1.59       |
| P2Y6        | 1321N1      | UDP            | 3.6         | 7.78  | 1.17       |
| P2Y11       | 1321N1      | ATP            | 3.6         | 5.70  | 1.78       |
| PAR1        | KOLF, mouse | Thrombin       | 3.9         | 6.86  | 1.17       |
| PTAFR       | HEK         | PTAF           | 1.9         | 8.06  | 0.65       |
| V1a         | CHO         | Vasopressin    | 3.7         | 8.07  | 0.81       |
| V1b         | CHO         | Vasopressin    | 3.3         | 8.20  | 0.99       |
| V2          | CHO         | Vasopressin    | 2.7         | 7.42  | 1.57       |
| Oxytocin    | CHO         | Oxytocin       | 3.1         | 7.29  | 1.21       |
| GPR40       | HEK T       | Linolenic acid | 12          | 5.05  |            |
| GPR40       | HEK T       | Lauric acid    | 10          | 5.15  |            |
| GPR41       | HEK T       | Propionate     | 4.2         | 7.18  |            |
| GPR43       | HEK T       | Propionate     |             | 6.43  |            |
|             |             |                |             |       |            |



**Figure 35.** Representative dose-response curves of Calcium mobilization assays. Some curves are presented after normalization with reference agonist activity as 100% and basal as 0%.
















































### mGlu1 Receptor




### mGlu5 Receptor







### MC3 Receptor



Hank's Balanced Salt Solution,  
20 mM HEPES, 2.5 mM probenecid







**MC4 Receptor**  
**Hank's Balanced Salt Solution,**  
**20 mM HEPES, 2.5 mM probenecid**


























**GPR40 Receptor**  
Hank's Balanced Salt Solution,  
20 mM HEPES, 2.5 mM probenecid




**GPR41 Receptor**  
Hank's Balanced Salt Solution,  
20 mM HEPES, 2.5 mM probenecid





**GPR43 Receptor**  
**Hank's Balanced Salt Solution,**  
**20 mM HEPES, 2.5 mM probenecid**




**Figure 36.** Schild analysis with a PDSP compound at 5-HT<sub>2A</sub> receptors.

**Figure 37.** Schild analysis of a PDSP compound at 5-HT<sub>2B</sub> receptors



**Figure 38.** Schild analysis with a PDSP compound at 5-HT<sub>2C</sub> receptors.

**Figure 39.** Representative figures for bias analysis. LSD and Ergotamine (ERG) agonist activity in  $G_i$  (for 5-HT<sub>1B</sub>) or  $G_q$  (for 5-HT<sub>2B</sub>) or  $\beta$ -arrestin signaling pathways were determined at 5-HT<sub>1B</sub> and 5-HT<sub>2B</sub> receptors as outlined in the functional assay section, and results were normalized to corresponding 5-HT activity and analyzed in Prism using the Black and Leff operational model to determine the transduction coefficient,  $\text{Log}(\tau/K_A)$ , as listed in **Tables 26 and 27**.



**Table 26.** Transduction coefficients,  $\text{Log}(\tau/K_A)$ , and bias factor calculation for indicated pathways of agonists at 5-HT<sub>1B</sub>. See representative concentration-response curves in **Figure 39**. Bias factor =  $10^{\Delta\Delta\text{Log}(\tau/K_A)}$ .

| Ligands     | Bias calculation at 5-HT <sub>1B</sub> receptors |                              |                           |                              |                                    |             |
|-------------|--|------------------------------|---------------------------|------------------------------|------------------------------------|-------------|
|             | G <sub>i</sub> pathway                           | $\Delta\text{Log}(\tau/K_A)$ | $\beta$ -arrestin pathway | $\Delta\text{Log}(\tau/K_A)$ | $\Delta\Delta\text{Log}(\tau/K_A)$ | Bias Factor |
| 5-HT        | 9.54 ± 0.08                                      | 0                            | 7.47 ± 0.13               | 0                            | 0                                  | 1.0         |
| LSD         | 9.53 ± 0.09                                      | -0.01                        | 8.08 ± 0.29               | 0.61                         | 0.62                               | 4.2         |
| Ergotamine  | 8.53 ± 0.18                                      | -1.01                        | 7.82 ± 0.20               | 0.35                         | 1.36                               | 22.9        |
| DHE         | 8.89 ± 0.21                                      | -0.65                        | 7.87 ± 0.19               | 0.40                         | 1.05                               | 11.2        |
| MTE         | 9.41 ± 0.12                                      | -0.13                        | 7.72 ± 0.15               | 0.25                         | 0.38                               | 2.4         |
| Pergolide   | 7.82 ± 0.25                                      | -1.72                        | 6.41 ± 0.11               | -1.06                        | 0.66                               | 4.6         |
| Cabergoline | 6.58 ± 0.08                                      | -2.96                        | 6.32 ± 0.14               | -1.15                        | 1.81                               | 64.6        |
| Ro 60-0175  | 6.37 ± 0.35                                      | -3.17                        | 5.32 ± 0.53               | -2.15                        | 1.02                               | 10.5        |
| Sumatriptan | 8.16 ± 0.19                                      | -1.38                        | 6.94 ± 0.11               | -0.53                        | 0.85                               | 7.1         |
| Donitriptan | 9.51 ± 0.37                                      | -0.03                        | 8.13 ± 0.18               | 0.66                         | 0.69                               | 4.9         |
| Eletriptan  | 8.13 ± 0.29                                      | -1.41                        | 7.07 ± 0.29               | -0.40                        | 1.01                               | 10.2        |
| Rizatriptan | 7.80 ± 0.43                                      | -1.74                        | 6.63 ± 0.21               | -0.84                        | 0.90                               | 7.9         |


DHE = Dihydroergotamine; MTE = Methylergonovine


**Table 27.** Transduction coefficients,  $\text{Log}(\tau/K_A)$ , and bias factor calculation for indicated pathways of agonists at 5-HT<sub>2B</sub>. See representative concentration-response curves in **Figure 39**.

| Ligands     | Bias calculation at 5-HT <sub>2B</sub> receptors |                              |                           |                              |                                    |             |
|-------------|--|------------------------------|---------------------------|------------------------------|------------------------------------|-------------|
|             | G <sub>q</sub> pathway                           | $\Delta\text{Log}(\tau/K_A)$ | $\beta$ -arrestin pathway | $\Delta\text{Log}(\tau/K_A)$ | $\Delta\Delta\text{Log}(\tau/K_A)$ | Bias Factor |
| 5-HT        | 9.61 ± 0.05                                      | 0                            | 8.30 ± 0.09               | 0                            | 0                                  | 1.0         |
| LSD         | 7.63 ± 0.15                                      | -1.98                        | 8.38 ± 0.09               | 0.08                         | 2.06                               | 114.8       |
| Ergotamine  | 5.95 ± 0.10                                      | -3.66                        | 7.05 ± 0.16               | -1.25                        | 2.41                               | 257.0       |
| DHE         | 6.03 ± 0.07                                      | -3.58                        | 7.13 ± 0.33               | -1.17                        | 2.41                               | 257.0       |
| MTE         | 7.93 ± 0.23                                      | -1.68                        | 8.05 ± 0.20               | -0.25                        | 1.43                               | 26.9        |
| Pergolide   | 7.50 ± 0.22                                      | -2.11                        | 7.96 ± 0.08               | -0.34                        | 1.77                               | 58.9        |
| Cabergoline | 6.92 ± 0.18                                      | -2.69                        | 7.76 ± 0.22               | -0.54                        | 2.15                               | 141.3       |
| Ro 60-0175  | 9.04 ± 0.15                                      | -0.57                        | 9.00 ± 0.10               | 0.7                          | 1.27                               | 18.6        |
| Sumatriptan | 5.24 ± 0.12                                      | -4.37                        | 5.17 ± 0.89               | -3.13                        | 1.24                               | 17.4        |
| Donitriptan | 6.65 ± 0.12                                      | -2.96                        | 5.92 ± 0.04               | -2.38                        | 0.58                               | 3.8         |
| Eletriptan  | 6.17 ± 0.05                                      | -3.44                        | 5.41 ± 0.20               | -2.89                        | 0.55                               | 3.5         |
| Rizatriptan | 5.86 ± 0.15                                      | -3.75                        | 5.03 ± 0.03               | -3.27                        | 0.48                               | 3.0         |

DHE = Dihydroergotamine; MTE = Methylergonovine

**Figure 40.** Comparison of bias factors of various ligands at 5-HT<sub>1B</sub> and 5-HT<sub>2B</sub> receptors. Values are taken from Tables 26 and 27 and Figure 33. \* indicates a significant difference ( $p < 0.0001$  by two-way ANOVA, selected drugs at two pathways). Results (Figure 39, Tables 26 and 27, and Figure 40) are from a recent *Science* paper (DOI: [10.1126/science.1232808](https://doi.org/10.1126/science.1232808)).





**Figure 41.** Representative curves for the allosteric operational model. GPR68-mediated cAMP production was determined in the absence and presence of increasing concentrations of PDSP compound 102022. Results were analyzed using the allosteric operational model and best-fit values are listed in the following table. In the curve-fitting,  $\text{Log}(K_A)$  is set at -6.50 (equivalent to the  $\text{pEC}_{50}$  of protons in the absence of 102022), while the  $E_{\text{max}}$  is constrained at 800 (which is the maximum activity of the system), and  $\tau_B$  is constrained at “0” since 102022 has no agonist activity by itself.

| [102022] M       | 0         | 3.00E-7   | 1.00E-6   | 3.00E-6   | 1.00E-5   | 3.00E-5   | Global (shared) |
|------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------------|
| Best-fit values  |           |           |           |           |           |           |                 |
| $\text{Log}K_A$  | (= -6.50) | (= -6.50) | (= -6.50) | (= -6.50) | (= -6.50) | (= -6.50) |                 |
| $\text{Log}K_B$  | -6.14     | -6.14     | -6.14     | -6.14     | -6.14     | -6.14     | -6.14           |
| Basal            | 16.32     | 16.32     | 16.32     | 16.32     | 16.32     | 16.32     | 16.32           |
| $E_{\text{max}}$ | (=800)    | (=800)    | (=800)    | (=800)    | (=800)    | (=800)    |                 |
| $\tau_A$         | 0.678     | 0.678     | 0.678     | 0.678     | 0.678     | 0.678     | 0.678           |
| $\alpha$         | 0.877     | 0.877     | 0.877     | 0.877     | 0.877     | 0.877     | 0.877           |
| $\beta$          | 4.055     | 4.055     | 4.055     | 4.055     | 4.055     | 4.055     | 4.055           |
| B                | = 0.0     | 3.00E-07  | 1.00E-06  | 3.00E-06  | 1.00E-05  | 3.00E-05  |                 |
| $\tau_B$         | = 0.0     | = 0.0     | = 0.0     | = 0.0     | = 0.0     | = 0.0     |                 |
| n                | 1.443     | 1.993     | 2.441     | 2.799     | 3.116     | 3.451     |                 |

## 2.5. Functional assays for G<sub>i</sub> or G<sub>s</sub> coupled GPCRs - Split luciferase biosensor

### cAMP assays

**Main equipment:** luminescence counter

**Reagents:** GloSensor cAMP construct from Promega and Luciferin

**Assay buffer:** 20 mM HEPES, 1x HBSS, pH 7.40

**Note:** See section 2.5.5 protocols for G<sub>i</sub> coupled mGluRs

**2.5.1. Cell culture and transfection:** To determine G<sub>i</sub> or G<sub>s</sub> GPCR mediated cAMP production, the PDSP uses Promega's split luciferase based GloSensor cAMP biosensor technology. With cells stably expressing target receptors, we transfect with the GloSensor cAMP DNA construct overnight; otherwise, we co-transfect HEK293 T cells with target receptor DNA and GloSensor cAMP DNA construct overnight. For detailed calcium phosphate transfection protocol, see [Section 2.4.1.2](#). To prepare plates for assays, cells are seeded into PLL coated 384-well white clear bottom cell culture plates with DMEM supplemented with 1% dFBS at a density of 15-20,000 cells in 40 µl medium per well. The plates can be used for assays after 6 hours or overnight.

**2.5.2. Split luciferase biosensor cAMP assay – Luciferin first protocol:** On the day of the assay, wells are emptied of culture medium and loaded for 90 min at 37°C with 20 µl of 4 mM Luciferin prepared in assay buffer. All the following steps are carried out at room temperature. To measure agonist activity at G<sub>s</sub> coupled receptors, 10 µl of 3x drug solutions is added per well and plates are counted for chemiluminescence after 15 minutes. To measure antagonist activity at G<sub>s</sub> coupled receptors, cells are preincubated with test compounds for 15 minutes before addition of an EC<sub>80</sub> concentration of a reference agonist, and chemiluminescence is counted after 15 minutes. To measure agonist activity at G<sub>i</sub> coupled receptors, 10 µl 3x test drug solution is added for 15 minutes before addition of 5 µl of isoproterenol at a final concentration of 200 nM, followed by counting of the plate for chemiluminescence after 15 minutes. To measure antagonist activity at G<sub>i</sub> coupled receptors, cells are preincubated with test drug for 15 minutes

before an EC<sub>80</sub> concentration of a reference agonist is added for another 15 minutes, then, 5 µl of isoproterenol at a final concentration of 200 nM is added, and counting is done after 15 minutes. In these experiments, isoproterenol is used to activate the endogenous Gs protein via the endogenous β<sub>2</sub> adrenergic receptors expressed in HEK293 T cells. Different receptors and different cell lines may need different preincubation times; when these are done, especially before larger scale screening assays, preliminary assays are performed to optimize count windows and incubation times.

**2.5.3. Split luciferase biosensor cAMP assay – Drug first protocol:** On the day of the assay, wells are emptied of culture medium and receive 20 µl/well assay buffer, followed by addition of 10 µl of 3x drug solutions for 15 minutes at room temperature. To measure agonist activity for Gs coupled receptors, 10 µl of 4 mM Luciferin prepared in assay buffer is added and counting is done after 15 minutes. To measure agonist activity for Gi coupled receptors, 10 µl of 4 mM Luciferin supplemented with isoproterenol at a final concentration of 200 nM is added, and counting is done after 15 minutes. To measure antagonist activity at Gi coupled receptors, cells are preincubated with drugs for 15 minutes before addition of an EC<sub>80</sub> concentration of a reference agonist for another 15 minutes; this is followed by addition of 10 µl of 4 mM Luciferin supplemented with isoproterenol at a final concentration of 200 nM, and counting is done after 15 minutes. Different receptors and different cell lines may need different preincubation times; preliminary assays are performed to optimize incubation times and count windows before large scale screening assays.

**2.5.4. Split luciferase biosensor cAMP assay protocol for G<sub>i</sub>-coupled mGluRs: mGluR<sub>2</sub>, mGluR<sub>3</sub>, mGluR<sub>4</sub>, mGluR<sub>6</sub>, and mGluR<sub>8</sub>.**

Assay Buffer 1 (mGluR<sub>2</sub>, mGluR<sub>4</sub>, mGluR<sub>6</sub>, mGluR<sub>8</sub>): Modified Locke's buffer containing 125 mM NaCl, 31 mM NaGluconate, 5.6 mM KCl, 3.6 mM NaHCO<sub>3</sub>, 1 mM MgCl<sub>2</sub>, 1.3 mM CaCl<sub>2</sub>, 5.6 mM glucose, and 20 mM HEPES, pH 7.4)

Assay Buffer 2 (mGluR<sub>3</sub>): Modified Locke's buffer containing 156 mM Na Gluconate, 5.6 mM KCl, 3.6 mM NaHCO<sub>3</sub>, 1 mM MgCl<sub>2</sub>, 1.3 mM CaCl<sub>2</sub>, 5.6 mM glucose, and 20 mM HEPES, pH 7.4)

**2.5.4.1. General experimental procedure:** The second type of glutamate receptor functional assays consists of measurements of cAMP production and is used with receptors negatively coupled to adenylyl cyclase (Gi-coupled receptors). These are mGluR<sub>2</sub>, mGluR<sub>3</sub>, mGluR<sub>4</sub>, mGluR<sub>6</sub>, and mGluR<sub>8</sub> receptors expressed in CHO cells also stably expressing the GloSensor cAMP biosensor. The activity of agonists at these receptors is determined by measurement of their ability to decrease forskolin-induced elevation of cyclic AMP formation. Cells cultured in white-walled, clear-bottom, 96-well culture plates are preincubated for 1 hour at room temperature in 100  $\mu$ l of the aforementioned Locke's medium containing 6% w/v D-Luciferin, as a substrate for the GloSensor enzyme. After 1 hour, the basal RLU's are measured 5 times, every 2 min. 1  $\mu$ M forskolin and 6% w/v D-Luciferin is added without or with mGluR agonists, and the incubation is continued for 16 min. After the incubation, relative luminescence units (RLUs) are measured. All of these studies are performed in the absence and in the presence of agonists appropriate for the different mGluRs, and used at concentrations equivalent to their EC<sub>50</sub> values. Hence, 3  $\mu$ M glutamate is used for mGluR<sub>2</sub>, mGluR<sub>3</sub>, 4  $\mu$ M glutamate is used for mGluR<sub>4</sub>, 2  $\mu$ M glutamate is used for mGluR<sub>6</sub>, and 0.7  $\mu$ M glutamate is used for mGluR<sub>8</sub>.

**2.5.4.2. Primary assays - Single concentration assays.** Each new compound is tested on selected receptors at a single concentration (10  $\mu$ M) for activity as an agonist or an antagonist. Testing for antagonism is performed in presence of the EC<sub>50</sub> concentration of a typical agonist (as described above). Each compound is tested in duplicate in two separate experiments performed on different cell passages. In addition to the tested compounds, each 96-well plate contains points for the determination of basal activity, maximal agonist stimulation, agonist EC<sub>50</sub> concentrations (i.e., concentration-response isotherm), and the IC<sub>50</sub> concentration of a known antagonist for purposes of positive control and for activity calculations. The reported results for each compound are calculated for agonists as the percent of maximal activity (as obtained with maximal agonist concentrations), and for antagonists as the percent of inhibition of receptor activity (in presence of an EC<sub>50</sub> concentration of the agonist). Results are expressed as means  $\pm$  SEM from four replicates.

**2.5.4.3. Secondary assays - Concentration-response assays.** Compounds determined to be active as agonists or antagonists may be tested for their potency in concentration-response experiments. Eight-point concentration-response curves are performed in duplicate twice on two separate passages of cells (sometimes a third curve may be needed if in the first experiment the range of concentrations used is outside of the active range). For antagonists, these curves are performed in the presence of the EC<sub>50</sub> concentration of the agonist. For each compound, the results from four replicates are averaged and then either the EC<sub>50</sub> or IC<sub>50</sub> values are calculated by non-linear regression using the 4-parameter logistic equation. Results are reported as EC<sub>50</sub> or IC<sub>50</sub> values for each tested compound (and receptor) and also include the EC<sub>50</sub> or IC<sub>50</sub> values of a known agonist or antagonist for comparison purposes.

**2.5.5. Data processing and analysis:** The luminescence counter records chemiluminescence in relative luminescence units (RLU) and saves files in Excel spreadsheets for easy processing. Results in RLU are plotted and analyzed in GraphPad Prism v5.0 as outlined in **Section 2.3**.

**2.5.6. Table and Figures.** List of GPCRs for which the PDSP provides functional assays to determine Gi or Gs activity and their corresponding concentration-response curves. The PDSP will also provide functional assays for other Gi or Gs coupled GPCRs upon request.

**Table 28.** List of GPCRs for which the PDSP provides cAMP measurement using GloSensor cAMP technology for Gi or Gs coupled receptors, and representative figures.

| Receptor | Cell line   | Reference     | Emax (fold) <sup>#</sup> | pEC50 | Hill slope |
|----------|-------------|---------------|--------------------------|-------|------------|
| M2       | HEK T       | Acetylcholine | 2.2                      | 7.63  | 0.72       |
| M2D      | HEK T       | CNO           | 3.2                      | 6.27  | 0.86       |
| M4       | HEK T       | Acetylcholine | 2.6                      | 8.16  | 1.07       |
| M4D      | HEK T       | CNO           | 3.6                      | 7.07  | 1.78       |
| GsDREADD | HEK T       | CNO           | 9.9                      | 7.87  | 1.53       |
| 5-HT1A   | HEK T       | 5-HT          | 5.8                      | 8.71  | 1.02       |
| 5-HT1B   | HEK T       | 5-HT          | 5.1                      | 9.71  | 0.94       |
| 5-HT1D   | HEK T       | 5-HT          | 2.0                      | 9.44  | 1.77       |
| 5-HT1E   | HEK T       | 5-HT          | 3.2                      | 9.12  | 0.87       |
| 5-HT1F   |             |               |                          |       |            |
| 5-HT4    | HEK T       | 5-HT          | 2.3                      | 10.18 | 0.98       |
| 5-HT5A   |             |               |                          |       |            |
| 5-HT6    | HEK T       | 5-HT          | 24.5                     | 7.06  | 1.25       |
| 5-HT7A   | HEK T       | 5-HT          | 60.1                     | 8.60  | 1.58       |
| A2A      | HEK T       | NECA          | 2.4                      | 9.58  | 0.74       |
| A2A      | HEK T       | CCPA          | 1.9                      | 7.80  | 1.42       |
| A2A      | HEK T       | CGS21680      | 2.2                      | 9.67  | 0.72       |
| CRF-1    | HEK T       | CRF           | 100                      | 8.57  | 0.76       |
| CRF-2    | HEK T       | CRF           | 122                      | 7.26  | 0.95       |
| D1       | HEK T       | Dopamine      | 148.1                    | 8.18  | 0.94       |
| D2       | HEK T       | Dopamine      | 2.8                      | 8.99  | 1.01       |
| D3       |             |               |                          |       |            |
| D4       | HEK T       | Dopamine      | 2.0                      | 9.19  | 1.67       |
| D5       | HEK T       | Dopamine      | 37.6                     | 9.24  | 0.80       |
| H2       | HEK T       | Histamine     | 40.0                     | 8.07  | 1.08       |
| H3       | HEK T       | Histamine     | 4.6                      | 8.90  | 0.70       |
| β2       | HEK T       | Isoproterenol | 54.3                     | 10.71 | 0.94       |
| DOR      | HEK T       | DADLE         | 2.7                      | 7.36  | 0.94       |
| KOR      | HEK T       | Salvinorin A  | 2.4                      | 9.43  | 1.33       |
| MOR      | HEK T       | DAMGO         | 3.2                      | 8.94  | 0.75       |
| NOP      | HEK T       | Nociceptin    | 3.6                      | 8.64  | 1.06       |
| MC1      | HEK T       | α-MSH         | 31.5                     | 7.47  | 1.01       |
| MC3      | HEK T       | α-MSH         | 188.7                    | 7.31  | 1.05       |
| mGluR2   | CHO         | Glutamate     | 75.2%*                   | 5.62  | 1.30       |
| mGluR3   | CHO         | Glutamate     | 83.2%*                   | 5.54  | 1.60       |
| mGluR4   | CHO         | Glutamate     | 91.0%*                   | 5.42  | 1.50       |
| mGluR6   | CHO         | Glutamate     | 60.2%*                   | 5.92  | 1.50       |
| mGluR7   | In progress |               |                          |       |            |


| Receptor | Cell line | Reference      | Emax (fold) <sup>#</sup> | pEC50 | Hill slope |
|----------|-----------|----------------|--------------------------|-------|------------|
| mGluR8   | CHO       | Glutamate      | 52.6%*                   | 6.18  | 1.70       |
| SSTR5    | HEK T     | Somatostatin   | 7.8                      | 8.63  | 1.66       |
| GPR88    | HEK T     | PDSP compd     | 6.1                      | 6.00  | 0.74       |
| GPR68    | HEK T     | H <sup>+</sup> | 23.2                     | 6.50  | 1.64       |
|          |           |                |                          |       |            |
|          |           |                |                          |       |            |



#### Notes



\*: percentage inhibition of Forskolin-activated cAMP production



#: Emax in fold for Gi pathway represents the ratio of basal vs maximal activity; Emax in fold for Gs pathway represents the ratio of maximal vs basal activity.



**Figure 42.** Representative concentration-response curves for the GloSensor cAMP assay. The assays were conducted according to the above procedure and were analyzed using GraphPad Prism 5.0.





























































### mGlu4 Receptor



### mGlu6 Receptor







## 2.6. GPCR Tango assays: G-protein independent $\beta$ -arrestin recruitment

**Main equipment:** Luminescence counter

**Main reagents:** BrightGlo<sup>®</sup> reagents from Promega


**Tango assay buffer:** 20 mM HEPES, 1x HBSS, pH 7.40


**2.6.1. Tango construct design and cell culture:** To measure GPCR mediated  $\beta$ -arrestin translocation activity, we adopted the Tango assay system originally developed by Richard Axel and his colleagues [Barnea et al., PNAS (2008) 105:64-69]. GPCR tango constructs (**Figure 43**, below) were codon-optimized for better expression in mammalian cell lines, synthesized by Blue Heron Biotech (Bothell, WA), and finally confirmed by sequencing. The HTLA cells (an HEK293 cell line stably expressing a tTA-dependent luciferase reporter and a  $\beta$ -arrestin2-TEV fusion gene) were gifted from Richard Axel's lab, and are maintained in DMEM supplemented with 10% FBS and 2  $\mu$ g/ml Puromycin and 100  $\mu$ g/ml Hygromycin B. An N-terminal Flag epitope tag was designed into the GPCR tango constructs for confirmation of surface expression (**Figure 44**) and comparison of expression levels.

**2.6.2. GPCR  $\beta$ -arrestin recruitment assay:** HTLA cells are transfected with GPCR tango constructs overnight (see [Section 2.4.1.2](#) for detailed transfection protocol) and are plated in Poly-L-Lys (PLL) coated 384-well white clear-bottom cell culture plates at a density of 15,000 cells in 40  $\mu$ l per well of DMEM with 1% dFBS. The cells are incubated for at least 6 hours (usually overnight) for them to recover before receiving drug stimulation. Drug stimulation solutions are prepared in filtered Tango assay buffer at 5x and added to cells (10  $\mu$ l per well) for overnight. To measure antagonist activity, drug solutions are made at 6x of the final concentration and are preincubated with cells for 30 min before addition of 10  $\mu$ l of and EC<sub>80</sub> concentration of a reference agonist. The EC<sub>80</sub> concentration is determined in separate preliminary dose-response assays. On the day of measurement, medium and drug solutions are removed and 20  $\mu$ l per well of BrightGlo reagent (diluted by 20-fold with Tango assay buffer) is added. The plate is

incubated for 20 minutes at room temperature in the dark before being counted on a luminescence counter.

**Figure 43.** Design of Tango constructs (**upper panel**) and GPCR Tango assay principle (**lower panel**). GPCR Tango constructs are designed as indicated in the top panel. Each construct contains the following elements in this order: (1) signal/FLAG tag; (2) receptor gene of interest; (3) Vasopressin 2 receptor C-tail; (4) TEV protease cleavage site; (5) Tta transcription factor. The GPCR Tango assay is carried out in transiently transfected HEK293 T cells genetically modified to express  $\beta$ -arrestin fused with a TEV protease and a Tta mediated Luciferase reporter gene. Activation of the transfected GPCR leads to  $\beta$ -arrestin translocation, which guides TEV protease to cleave the Tta transcription factor from the GPCR tail. The free Tta transcription factors then is transported to the cell nucleus, where it activates the luciferase reporter gene. The principle diagram is adapted from Barnea et al., (2008) *The genetic design of signaling cascades to record receptor activation. PNAS 105:4-69.*





**Figure 44.** Confirmation of surface expression of FLAG-tagged non-orphan non-olfactory GPCRs used for the  $\beta$ -arrestin recruitment assay. Cells were pre-fixed with paraformaldehyde (PFA) for 30 min, blocked, probed using rabbit anti-FLAG antibody (primary) incubated for 1h at room temperature and then overnight at 4°C. On day 2, the plate was incubated with Alexa Fluor 594-conjugated goat anti-rabbit antibody (secondary) and Hoechst 33342 dye for 1 h. After thorough washing, the plate was post-fixed with PFA and stored at 4°C in the dark. Images were taken using the Pathway High Throughput Bioimager. Blue crosses indicated constructs that were not expressed on the surface of cells; 146 of 152 receptors (96%) were expressed on the cell surface.

**2.6.3. Data processing and analysis:** The luminescence counter records chemiluminescence in relative luminescence units (RLU) and saves files in Excel spreadsheets for easy processing. Results in RLU are plotted against concentrations and analyzed in GraphPad Prism v5.0 as outlined in **Section 2.3**.

**Table 29.** List of GPCR Tango constructs and validation. All DNA constructs are human clones, designed according to Barnea et al., (2008), codon optimized, synthesized by Blue Heron Biotech, and confirmed by sequencing. The  $\beta$ -arrestin recruitment assays using GPCR Tango constructs were conducted according to the above  $\beta$ -arrestin recruitment assay procedures and results were analyzed in GraphPad Prism 5.0.

\* “Yes” for those assays are available in the PDSP; “No” for those not yet available; “In Progress” for those being synthesized or validated; “TBO” for those being optimized and ready soon.

| Target Gene | IUPHAR Name | Reference Ligand(s) | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page       | Status*     |
|-------------|-------------|---------------------|-------------------------|-------------------|----------------------|-------------|
| HTR1A       | 5-HT1A      | 5-HT                | 15                      | 6.35              | <a href="#">P198</a> | Yes         |
| HTR1A       | 5-HT1A      | LSD                 | 13                      | 7.21              | <a href="#">P198</a> | Yes         |
| HTR1B       | 5-HT1B      | 5-HT                | 3.5                     | 7.48              | <a href="#">P198</a> | Yes         |
| HTR1B       | 5-HT1B      | LSD                 | 2.8                     | 8.21              | <a href="#">P198</a> | Yes         |
| HTR1D       | 5-HT1D      | 5-HT                | 3.5                     | 7.16              | <a href="#">P198</a> | Yes         |
| HTR1D       | 5-HT1D      | LSD                 | 2.9                     | 8.45              | <a href="#">P198</a> | Yes         |
| HTR1E       | 5-HT1E      | 5-HT                | 3.9                     | 8.42              | <a href="#">P199</a> | Yes         |
| HTR1E       | 5-HT1E      | LSD                 | 3.6                     | 8.07              | <a href="#">P199</a> | Yes         |
| HTR1F       | 5-HT1F      | 5-HT                | 100                     | 6.96              | <a href="#">P199</a> | Yes         |
| HTR1F       | 5-HT1F      | LSD                 | 83                      | 7.42              | <a href="#">P199</a> | Yes         |
| HTR2A       | 5-HT2A      | 5-HT                | 20                      | 6.60              | <a href="#">P199</a> | Yes         |
| HTR2A       | 5-HT2A      | LSD                 | 10                      | 8.72              | <a href="#">P199</a> | Yes         |
| HTR2B       | 5-HT2B      | 5-HT                | 3.3                     | 8.75              | <a href="#">P200</a> | Yes         |
| HTR2B       | 5-HT2B      | LSD                 | 2.9                     | 8.93              | <a href="#">P200</a> | Yes         |
| HTR2C-INI   | 5-HT2C      | 5-HT                | 2.0                     | 7.42              | <a href="#">P200</a> | Yes         |
| HTR2C-INI   | 5-HT2C      | LSD                 | 3.1                     | 7.41              | <a href="#">P200</a> | Yes         |
| HTR2C-VGV   |             | LSD                 |                         |                   |                      | NO          |
| HTR2C-VNV   |             | LSD                 | 1.8                     | 8.30              | <a href="#">P200</a> | TBO         |
| HTR2C-VSV   |             | LSD                 | 1.6                     | 8.82              | <a href="#">P200</a> | TBO         |
| HTR4        | 5-HT4       | 5-HT                | 47                      | 7.69              | <a href="#">P201</a> | Yes         |
| HTR4        | 5-HT4       | LSD                 | 3.8                     | <5                | <a href="#">P201</a> | Yes         |
| HTR5A       | 5-HT5A      | 5-HT                | 12                      | 8.07              | <a href="#">P201</a> | Yes         |
| HTR5A       | 5-HT5A      | LSD                 | 16                      | 9.83              | <a href="#">P201</a> | Yes         |
| HTR6        | 5-HT6       | 5-HT                | 136                     | 6.53              | <a href="#">P201</a> | Yes         |
| HTR6        | 5-HT6       | LSD                 | 23                      | 7.85              | <a href="#">P201</a> | Yes         |
| HTR7D       | 5-HT7D      | 5-HT                |                         |                   |                      | TBO         |
| HTR7        | 5-HT7D      | LSD                 |                         |                   |                      | NO          |
| HTR7A       | 5-HT7A      | 5-HT, LSD           |                         |                   |                      | In Progress |

| Target Gene | IUPHAR Name | Reference Ligand(s)             | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page                               | Status* |
|-------------|-------------|---------------------------------|-------------------------|-------------------|--|---------|
| ADORA1      | A1          | NECA                            | 4.9                     | 7.98              | <a href="#">P202</a>                         | Yes     |
| ADORA2A     | A2A         | NECA                            |                         |                   |  | No      |
| ADORA2B     | A2B         | NECA                            |                         |                   |  | No      |
| ADORA3      | A3          | NECA                            |                         |                   |  | No      |
| ADRA1A      | α1A         | epinephrine                     |                         |                   |  | No      |
| ADRA1B      | α1B         | epinephrine                     | 10.4                    | 6.64              | <a href="#">P202</a>                         | Yes     |
| ADRA1D      | α1D         | norepinephrine                  | 2.4                     | 5.82              | <a href="#">P202</a>                         | Yes     |
| ADRA2A      | α2A         | clonidine                       | 5.4                     | 8.36              | <a href="#">P202</a>                         | No      |
| ADRA2B      | α2B         | norepinephrine                  | 8.2                     | 7.28              | <a href="#">P202</a>                         | Yes     |
| ADRA2C      | α2C         | clonidine                       | 152.9                   | 7.18              | <a href="#">P202</a>                         | Yes     |
| ADRB1       | β1          | epinephrine                     | 38.7                    | 5.49              | <a href="#">P203</a>                         | Yes     |
| ADRB2       | β2          | epinephrine                     | 21.4                    | 6.24              | <a href="#">P203</a>                         | Yes     |
| ADRB3       | β3          | epinephrine                     |                         |                   |  | No      |
| AGTR1       | AT1         | angiotensin II                  | 7.8                     | 7.97              | <a href="#">P203</a>                         | Yes     |
| AGTR2       | AT2         | angiotensin II                  |                         |                   |  | No      |
| APJ         | Apelin      | apelin-17                       | 5.3                     | 8.80              | <a href="#">P203</a>                         | Yes     |
| AVPR1A      | V1a         | vasopressin                     | 14.2                    | 8.58              | <a href="#">P222</a>                         | Yes     |
| AVPR1B      | V1b         | vasopressin                     | 6.3                     | 8.63              | <a href="#">P222</a>                         | Yes     |
| AVPR2       | V2          | vasopressin                     | 9.7                     | 8.18              | <a href="#">P222</a>                         | Yes     |
| BDKRB1      | B1          | bradykinin                      | 8.0                     | 5.41              | <a href="#">P203</a>                         | Yes     |
| BDKRB2      | B2          | bradykinin                      | 4.8                     | 8.84              | <a href="#">P203</a>                         | Yes     |
| NMBR        | BB1         | Bombesin                        | 3.3                     | 8.16              | <a href="#">P204</a>                         | Yes     |
| NMBR        | BB1         | Gastrin-releasing peptide (GRP) | 4.6                     | 6.98              | <a href="#">P204</a><br><a href="#">P215</a> | Yes     |
| NMBR        | BB1         | Neuromedin B                    | 6.0                     | 8.68              | <a href="#">P204</a>                         | Yes     |
| GRPR        | BB2         | Bombesin                        | 6.2                     | 8.50              | <a href="#">P204</a>                         | Yes     |
| GRPR        | BB2         | Gastrin-releasing peptide (GRP) | 6.1                     | 8.32              | <a href="#">P204</a>                         | Yes     |
| GPRP        | BB2         | Neuromedin B                    | 5.3                     | 7.46              | <a href="#">P204</a>                         | Yes     |
| BRS3        | BB3         | PDSP compound                   | 13.2                    | 6.15              | <a href="#">P204</a>                         | Yes     |
| C3AR1       | C3a         | C3a (70-77)                     | 38.3                    | 6.15              | <a href="#">P204</a>                         | Yes     |
| CALCRb      | CT          | calcitonin                      |                         |                   |  | No      |
| CALCRL      | CT-like     | calcitonin                      |                         |                   |  | No      |
| CASR        | CaS         | spermine                        |                         |                   |  | No      |
| CCKAR       | CCK1        | [Thr28, Nle31]CCK(25-33)        | 3.7                     | 7.45              | <a href="#">P204</a>                         | Yes     |
| CCR4        | CCR4        | CCL22                           |                         |                   | <a href="#">P206</a>                         | Yes     |
| CCR6        | CCR6        | CCL20 (Exodus-1)                | 2.4                     | 7.79              | <a href="#">P206</a>                         | Yes     |
| CMKLR1      | GPCR27      | Chemerin                        | 15.7                    | 7.96              | <a href="#">P205</a>                         | Yes     |

| Target Gene | IUPHAR Name     | Reference Ligand(s)      | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page       | Status* |
|-------------|-----------------|--------------------------|-------------------------|-------------------|----------------------|---------|
| CNR1        | CB1             | WIN 55212-2              | 3.1                     | 6.25              | <a href="#">P207</a> | Yes     |
| CNR2        | CB2             | CP55940                  | 10.1                    | 6.57              | <a href="#">P207</a> | Yes     |
| CXCR1       | CXCR1           | IL-8                     | 38.1                    | <5.0              | <a href="#">P205</a> | Yes     |
| CXCR1       | CXCR1           | CXCL6                    | 3.5                     | 6.60              | <a href="#">P205</a> | Yes     |
| CXCR1       | CXCR1           | CXCL8                    | 4.1                     | 6.95              | <a href="#">P205</a> | Yes     |
| CXCR2       | CXCR2           | IL-8                     | 2.0                     | 8.53              | <a href="#">P205</a> | Yes     |
| CXCR2       | CXCR2           | CXCL6                    | 1.8                     | 7.30              | <a href="#">P205</a> | Yes     |
| CXCR2       | CXCR2           | CXCL8                    | 1.9                     | 6.88              | <a href="#">P205</a> | Yes     |
| CXCR4       | CXCR4           | SDF1- $\alpha$           | 2.0                     | 8.62              | <a href="#">P206</a> | Yes     |
| CXCR6       | CXCR6           | CXCL16 (SRPOX)           | 10.3                    | <5.0              | <a href="#">P206</a> | Yes     |
| CX3CR1      | CX3CR1          | SDF-1 $\alpha$           | 5.8                     | 9.52              | <a href="#">P206</a> | Yes     |
| CYSLTR1     | CysLT1          | leukotriene D4           | 4.5                     | 7.42              | <a href="#">P206</a> | Yes     |
| DRD1        | D1              | cabergoline              | 12.7                    | 5.09              | <a href="#">P207</a> | Yes     |
| DRD2        | D2              | LSD                      | 149.9                   | 9.30              | <a href="#">P207</a> | Yes     |
| DRD3        | D3              | quinpirole               | 23.4                    | 7.44              | <a href="#">P207</a> | Yes     |
| DRD4        | D4              | LSD                      |                         |                   |                      | No      |
| DRD5        | D5              | LSD                      | 40.3                    | 6.59              | <a href="#">P207</a> | Yes     |
| EDNRA       | ET <sub>A</sub> | endothelin-1             | 77.5                    | 9.21              | <a href="#">P208</a> | Yes     |
| EDNRB       | ET <sub>B</sub> | endothelin-1             |                         |                   |                      | No      |
| F2R         | PAR1            | TRAP                     |                         |                   |                      | No      |
| F2RL1       | PAR2            | TRAP                     |                         |                   |                      | No      |
| F2RL2       | PAR3            | TRAP                     |                         |                   |                      | No      |
| F2RL3       | PAR4            | TRAP                     |                         |                   |                      | No      |
| FPR1        | FPR1            | fMLP                     | 150.5                   | 8.64              | <a href="#">P208</a> | Yes     |
| FPR2        | FPR2/ALX        | fMLP                     | 18.3                    | <5.0              | <a href="#">P208</a> | Yes     |
| FPR3        | FPR3            | fMLP                     | 1.2                     | 6.64              | <a href="#">P208</a> | Yes     |
| FSHR        | FSH             | FSH                      |                         |                   |                      | No      |
| GAL1        | GAL1            | galanin                  | 8.9                     | 6.36              | <a href="#">P208</a> | Yes     |
| GAL2        | GAL2            | galnon                   | 3.6                     | 8.29              | <a href="#">P205</a> | Yes     |
| GAL3        | GAL3            | galanin                  | 3.6                     | 9.99              | <a href="#">P208</a> | Yes     |
| GHRH        | GHRH            | GRF (1-29)               |                         |                   |                      | No      |
| GHSR        | Ghrelin         | ghrelin                  | 1.7                     | 7.30              | <a href="#">P210</a> | Yes     |
| GIPR        | GIP             | GIP (1-42)               |                         |                   |                      | No      |
| GLP1R       | GLP-1           | glucagon                 | 6.7                     | <5.0              | <a href="#">P209</a> | Yes     |
| GLP2R       | GLP-2           | glucagon                 |                         |                   |                      | No      |
| GNRHR       | GnRH            | leuprolide               | 2.6                     | 8.47              | <a href="#">P209</a> | Yes     |
|             |                 |                          |                         |                   |                      |         |
| GPBA        | GPBA            | sodium taurodeoxycholate | 2.3                     | <5.0              | <a href="#">P209</a> | Yes     |

| Target Gene | IUPHAR Name  | Reference Ligand(s)       | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page                               | Status* |
|-------------|--------------|---------------------------|-------------------------|-------------------|--|---------|
| GPR44       | GPR44        | prostaglandin D2          | 4.0                     | 7.22              | <a href="#">P209</a>                         | Yes     |
| GPR75       | GPR75        | RANTES                    |                         |                   |  | No      |
| GRM1        | mGluR1       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM2        | mGluR2       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM3        | mGluR3       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM4        | mGluR4       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM5        | mGluR5       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM6        | mGluR6       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM7        | mGluR7       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| GRM8        | mGluR8       | D-homocysteic acid, L-Glu |                         |                   |  | No      |
| HRH1        | Histamine H1 | N-methylhistaprodifen     | 1.3                     | 7.21              | <a href="#">P210</a>                         | Yes     |
| HRH2        | H2           | Histamine                 | 6.9                     | 5.40              | <a href="#">P210</a>                         | Yes     |
| HRH3        | H3           | n-methylhistamine         | 15.8                    | 7.75              | <a href="#">P210</a>                         | Yes     |
| HRH4        | H4           | histamine                 | 3.1                     | 6.88              | <a href="#">P210</a>                         | Yes     |
| HCRTR1      | OX1          | orexin-A                  | 109.3                   | 7.78              | <a href="#">P209</a>                         | Yes     |
| HCRTR2      | OX2          | orexin-A                  | 92.8                    | 8.16              | <a href="#">P209</a>                         | Yes     |
| LHCGR       | LH           | LH luteinizing hormone    |                         |                   |  | No      |
| LPAR1       | LPA1         | 1-oleoyl LPA              | 6.4                     | 5.29              | <a href="#">P223</a>                         | Yes     |
| LPAR2       | LPA2         | 1-oleoyl LPA              | 291.0                   | <5.0              | <a href="#">P223</a>                         | Yes     |
| LPAR3       | LPA3         | LPA                       |                         |                   |  | No      |
| LPAR4       | LPA4         | LPA                       |                         |                   |  | No      |
| LPAR5       | LPA5         | 1-oleoyl LPA              | 14.2                    | <5.0              | <a href="#">P223</a>                         | Yes     |
| LPAR5       | LPA5         | LPA                       |                         |                   |  | No      |
| LTB4R       | BLT1         | leukotriene D4            | 33.2                    | 7.31              | <a href="#">P210</a>                         | Yes     |
| LTB4R       | BLT1         | Leukotriene B4            | 27.7                    | 7.55              | <a href="#">P211</a>                         | Yes     |
| LTB4R2b     | BLT2         | leukotriene D4            |                         |                   |  | No      |
| CHRM1       | M1           | Acetylcholine             | 2.1                     | 4.58              | <a href="#">P213</a>                         | Yes     |
| CHRM1       | M1           | Carbachol                 | 1.5                     | 6.74              | <a href="#">P212</a><br><a href="#">P213</a> | Yes     |

| Target Gene | IUPHAR Name | Reference Ligand(s)    | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page                               | Status* |
|-------------|-------------|------------------------|-------------------------|-------------------|--|---------|
| CHRM1       | M1          | Arecoline              | 3.7                     | <5.0              | <a href="#">P213</a>                         | Yes     |
| CHRM2       | M2          | Acetylcholine          | 8.1                     | 6.41              | <a href="#">P213</a>                         | Yes     |
| CHRM2       | M2          | Carbachol              | 3.6                     | 6.60              | <a href="#">P212</a><br><a href="#">P213</a> | Yes     |
| CHRM2       | M2          | Arecoline              | 22.7                    | 5.67              | <a href="#">P213</a>                         | Yes     |
| CHRM3       | M3          | Acetylcholine          | >10                     | <5.0              | <a href="#">P213</a>                         | Yes     |
| CHRM3       | M3          | Carbachol              | 12.2                    | <5.0              | <a href="#">P212</a><br><a href="#">P213</a> | Yes     |
| CHRM3       | M3          | Arecoline              | >10                     | <5.0              | <a href="#">P213</a>                         | Yes     |
| M3D         | M3DREADD    | CNO                    | 6.0                     | 6.57              | <a href="#">P213</a>                         | Yes     |
| CHRM4       | M4          | Acetylcholine          | 11.5                    | 5.61              | <a href="#">P213</a>                         | Yes     |
| CHRM4       | M4          | Carbachol              | 34.4                    | 5.67              | <a href="#">P212</a><br><a href="#">P213</a> | Yes     |
| CHRM4       | M4          | Arecoline              | 33.8                    | 5.62              | <a href="#">P213</a>                         | Yes     |
| CHRM5       | M5          | Acetylcholine          | <2.0                    | <5.0              | <a href="#">P213</a>                         | Yes     |
| CHRM5       | M5          | Carbachol              | 1.6                     | 6.91              | <a href="#">P212</a><br><a href="#">P213</a> | Yes     |
| CHRM5       | M5          | Arecoline              | 8.9                     | <5.0              | <a href="#">P213</a>                         | Yes     |
| MC1R        | MC1         | α-MSH                  | 2.7                     | 7.39              | <a href="#">P214</a>                         | Yes     |
| MC2R        | MC2         | ACTH                   |                         |                   |  | No      |
| MC3Rb       | MC3         | α-MSH                  |                         |                   |  | No      |
| MC4R        | MC4         | ACTH                   | 7.0                     | 6.53              | <a href="#">P214</a>                         | Yes     |
| MC5R        | MC5         | α-MSH                  | 2.0                     | 6.38              | <a href="#">P214</a>                         | Yes     |
| MCHR1       | MCH1        | [Ala-17]-MCH           | 7.0                     | 7.66              | <a href="#">P214</a>                         | Yes     |
| MCHR2       | MCH2        | [Ala-17]-MCH           | 17.2                    | 8.42              | <a href="#">P214</a>                         | Yes     |
| MLNR        | Motilin     | motilin                | 158.7                   | 8.08              | <a href="#">P214</a>                         | Yes     |
| MRGRPX1     | MRGPRX1     | BAM-22                 |                         |                   |  | No      |
| MRGRPX2     | MRGPRX2     | BAM-22                 |                         |                   |  | No      |
| MRGRPX3     | MRGPRX3     | BAM-22                 | 4.4                     | 10.14             | <a href="#">P215</a>                         | Yes     |
| MRGRPX4     | MRGPRX4     | BAM-22                 | 4.4                     | 8.62              | <a href="#">P215</a>                         | Yes     |
| MTNR1A      | MT1         | melatonin              | 3.8                     | 10.07             | <a href="#">P215</a>                         | Yes     |
| MTNR1B      | MT2         | melatonin              | 93.1                    | 8.18              | <a href="#">P215</a>                         | Yes     |
| NMUR1       | NMU1        | neuromedin S           | 48.1                    | 8.02              | <a href="#">P216</a>                         | Yes     |
| NMUR2       | NMU2        | neuromedin S           | 14.2                    | 7.92              | <a href="#">P215</a>                         | Yes     |
| NPS         | NPS         | neuropeptide S         | 58.2                    | 6.07              | <a href="#">P216</a>                         | Yes     |
| NPY1R       | NPY1        | pancreatic polypeptide | 3.8                     | 9.27              | <a href="#">P216</a>                         | Yes     |
| NPY2R       | NPY2        | pancreatic polypeptide | 17.0                    | 8.75              | <a href="#">P216</a>                         | Yes     |


| Target Gene | IUPHAR Name | Reference Ligand(s)     | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page                               | Status* |
|-------------|-------------|-------------------------|-------------------------|-------------------|--|---------|
| NPY4R       | NPY4        | neuropeptide Y          | 2.3                     | 6.89              | <a href="#">P216</a>                         | Yes     |
| NPY5R       | NPY5        | neuropeptide Y          |                         |                   |  | No      |
| NTSR1       | NTS1        | neurotensin             |                         |                   |  | No      |
| NTSR2       | NTS2        | SR48692                 | 10.6                    | 6.23              | <a href="#">P216</a>                         | Yes     |
| OPRD1       | δ           | DADLE                   | 7.1                     | 9.16              | <a href="#">P211</a>                         | Yes     |
| KOR         | κ           | salvinorin A            | 8.7                     | 7.74              | <a href="#">P211</a>                         | Yes     |
| OPRM1       | μ           | methadone               | 10.3                    | <5.0              | <a href="#">P211</a>                         | Yes     |
| OPRL1       | NOP         | orphanin                | 151.5                   | 7.92              | <a href="#">P211</a>                         | Yes     |
| OXER1       | OXE         | 5-oxo-ETE               |                         |                   |  | No      |
| OXTR        | Oxytocin    | oxytocin                | 26.9                    | 7.87              | <a href="#">P211</a><br><a href="#">P222</a> | Yes     |
| P2RY1       | P2Y1        | 2-MeSADP                |                         |                   |  | Yes     |
| P2RY12      | P2Y12       | 2-MeSADP                | 46.3                    | 7.34              | <a href="#">P217</a>                         | Yes     |
| P2RY13      | P2Y13       | 2-MeSADP                | 1.9                     | 7.54              | <a href="#">P217</a>                         | Yes     |
| P2RY14      | P2Y14       | UDP-glucose             | 1.5                     | 6.21              | <a href="#">P218</a>                         | Yes     |
| P2RY2       | P2Y2        | UTP                     | 67.1                    | <5.0              | <a href="#">P217</a>                         | Yes     |
| P2RY4       | P2Y4        | UTP                     | 6.0                     | 5.12              | <a href="#">P217</a>                         | Yes     |
| P2RY6       | P2Y6        | UDP                     | 1.8                     | 5.52              | <a href="#">P217</a>                         | Yes     |
| P2YR11      | P2Y11       | ATP                     | 4.7                     | <5.0              | <a href="#">P217</a>                         | Yes     |
| PTAFR       | PAF         | PAF (C16)               | 1.4                     | 6.30              | <a href="#">P218</a>                         | Yes     |
| PTGDR       | DP1         | prostaglandin D2        |                         |                   |  | No      |
| PTGER1      | EP1         | prostaglandin E2        | 8.0                     | 7.35              | <a href="#">P218</a>                         | Yes     |
| PTGER2      | EP2         | prostaglandin E2        | 66.0                    | 5.69              | <a href="#">P218</a>                         | Yes     |
| PTGER3      | EP3         | prostaglandin E2        | 3.5                     | 8.14              | <a href="#">P218</a>                         | Yes     |
| PTGER4      | EP4         | prostaglandin E2        | 38.7                    | 8.91              | <a href="#">P218</a>                         | Yes     |
| PTGFR       | FP          | prostaglandin F2α       | 21.8                    | 7.02              | <a href="#">P219</a>                         | Yes     |
| PTGIR       | IP1         | iloprost                | 10.9                    | <5.0              | <a href="#">P219</a>                         | Yes     |
| PTHR1       | PTH1        | PTH (1-42)              | 9.1                     | 6.41              | <a href="#">P219</a>                         | Yes     |
| PTHR2       | PTH2        | PTH (1-42)              |                         |                   |  | No      |
| S1PR1       | S1P1        | sphingosine-1-phosphate | 9.0                     | 5.97              | <a href="#">P219</a>                         | Yes     |
| S1PR2       | S1P2        | sphingosine-1-phosphate | 15.5                    | 5.77              | <a href="#">P219</a>                         | Yes     |
| S1PR3       | S1P3        | sphingosine-1-phosphate | 3.8                     | 5.68              | <a href="#">P219</a>                         | Yes     |
| S1PR4       | S1P4        | sphingosine-1-phosphate |                         |                   |  | No      |
| S1PR5       | S1P5        | sphingosine-1-phosphate |                         |                   |  | No      |







| Target Gene | IUPHAR Name | Reference Ligand(s) | E <sub>max</sub> (fold) | pEC <sub>50</sub> | Figure on page       | Status* |
|-------------|-------------|---------------------|-------------------------|-------------------|----------------------|---------|
| SCTR        | Secretin    | secretin            | 5.9                     | 5.29              | <a href="#">P220</a> | Yes     |
| SSTR1       | SST1        | somatostatin        | 6.7                     | 5.89              | <a href="#">P220</a> | Yes     |
| SSTR2       | SST2        | somatostatin        | 92.1                    | 5.69              | <a href="#">P220</a> | Yes     |
| SSTR3       | SST3        | somatostatin        | 119.2                   | 5.32              | <a href="#">P220</a> | Yes     |
| SSTR4       | SST4        | somatostatin        | 54.1                    | 6.06              | <a href="#">P220</a> | Yes     |
| SSTR5       | SST5        | somatostatin        | 254.2                   | 5.81              | <a href="#">P220</a> | Yes     |
| TA1         | TA1         | β-phenylethylamine  | 2.6                     | <5.0              | <a href="#">P221</a> | Yes     |
| TACR1       | NK1         | substance P         | 3.3                     | 6.71              | <a href="#">P221</a> | Yes     |
| TACR2       | NK2         | substance P         | 31.4                    | 8.80              | <a href="#">P221</a> | Yes     |
| TACR3       | NK3         | substance P         | 3.1                     | 5.76              | <a href="#">P221</a> | Yes     |
| TBXA2R      | TP          | thromboxane B2      |                         |                   |                      | No      |
| TRHR        | TRH1        | TRH                 |                         |                   |                      | No      |
| TSHR        | TSH         | thyrotropic hormone |                         |                   |                      | No      |
| UTSR        | UT          | urotensin II        | 9.5                     | 8.72              | <a href="#">P221</a> | Yes     |
| VIPR1       | VPAC1       | VIP                 | 7.7                     | <5.0              | <a href="#">P222</a> | Yes     |
| VIPR2       | VPAC2       | VIP                 | 12.9                    | <5.0              | <a href="#">P222</a> | Yes     |







**Table 29a.** List of available Tango constructs for peptide and orphan GPCRs at NIMH-PDSP. Some of them have been verified and are also listed in Table 29. Many others are in the process of being verified and/or optimized. IUPHAR target gene names are used.







|              |              |               |               |
|--------------|--------------|---------------|---------------|
| APJ-Tango    | GPR115-Tango | GPR21-Tango   | GPRC6A-Tango  |
| BB3-Tango    | GPR116-Tango | GPR22-Tango   | HCA1-Tango    |
| C5A-Tango    | GPR119-Tango | GPR25-Tango   | HCA2-Tango    |
| CCR1-Tango   | GPR120-Tango | GPR26-Tango   | HCA3-Tango    |
| CCR2-Tango   | GRP123-Tango | GPR27-Tango   | KISS-Tango    |
| CCR3-Tango   | GPR124-Tango | GPR31-Tango   | LPA4-Tango    |
| CCR4-Tango   | GPR125-Tango | GPR32-Tango   | MAS1L-Tango   |
| CCR5-Tango   | GPR126-Tango | GPR34-Tango   | MAS1-Tango    |
| CCR6-Tango   | GPR12-Tango  | GPR35-Tango   | MRGPRD-Tango  |
| CCR7-Tango   | GPR132-Tango | GPR37L1-Tango | MRGPRED-Tango |
| CCR8-Tango   | GPR133-Tango | GPR37-Tango   | MRGPRF-Tango  |
| CCR10-Tango  | GPR135-Tango | GPR39-Tango   | MRGPRG-Tango  |
| CCRL2-Tango  | GPR141-Tango | GPR3-Tango    | MRGPRX1-Tango |
| CD97-Tango   | GPR142-Tango | GPR45-Tango   | MRGPRX2-Tango |
| CMKLR1-Tango | GPR143-Tango | GPR4-Tango    | MRGPRX3-Tango |
| CX3C1-Tango  | GPR144-Tango | GPR50-Tango   | MRGPRX4-Tango |
| CXCR1-Tango  | GPR146-Tango | GPR52-Tango   | NPBW1-Tango   |
| CXCR2-Tango  | GPR148-Tango | GPR55-Tango   | NPBW2-Tango   |
| CXCR3-Tango  | GPR148-Tango | GPR56-Tango   | NPFF1-Tango   |
| CXCR4-Tango  | GPR150-Tango | GPR61-Tango   | NPFF2-Tango   |
| CXCR5-Tango  | GPR151-Tango | GPR62-Tango   | NPS-Tango     |
| CXCR6-Tango  | GPR152-Tango | GPR63-Tango   | OPN3-Tango    |
| CXCR7-Tango  | GPR153-Tango | GPR64-Tango   | OPN5-Tango    |
| ELTD1-Tango  | GPR156-Tango | GPR65-Tango   | OXGR1-Tango   |
| FFA1-Tango   | GPR157-Tango | GPR68-Tango   | PK1-Tango     |
| FFA2-Tango   | GPR158-Tango | GPR6-Tango    | PK2-Tango     |
| FFA3-Tango   | GPR15-Tango  | GPR75-Tango   | PRRP-Tango    |
| FPR1-Tango   | GPR160-Tango | GPR78-Tango   | QRFP-Tango    |
| FPR2-Tango   | GPR161-Tango | GPR82-Tango   | RXFP1-Tango   |
| FPR3-Tango   | GPR162-Tango | GPR83-Tango   | RXFP2-Tango   |
| GAL1-Tango   | GPR171-Tango | GPR84-Tango   | RXFP3-Tango   |
| GAL2-Tango   | GPR173-Tango | GPR85-Tango   | RXFP4-Tango   |
| GAL3-Tango   | GPR174-Tango | GPR87-Tango   | SUCNR1-Tango  |
| GPBA-Tango   | GPR17-Tango  | GPR88-Tango   | TA1-Tango     |
| GPBR-Tango   | GPR182-Tango | GPR92-Tango   | TAAR2-Tango   |
| GPR101-Tango | GPR183-Tango | GPR97-Tango   | TAAR5-Tango   |
| GPR110-Tango | GPR18-Tango  | GPRC5A-Tango  | TAAR6-Tango   |
| GPR111-Tango | GPR19-Tango  | GPRC5B-Tango  | TAAR8-Tango   |
| GPR113-Tango | GPR1-Tango   | GPRC5C-Tango  | TAAR9-Tango   |
| GPR114-Tango | GPR20-Tango  | GPRC5D-Tango  |               |







**Figure 45.** Representative concentration-response curves for GPCR mediated  $\beta$ -arrestin translocation. Curves are in essentially the same order as in **Table 29**.













































































































































































## 2.7. Parallel receptorome screening using $\beta$ -arrestin recruitment assay

**Main equipment:** Liquid handling workstation for 96- and 384-well plates, luminescence counter

**Main reagent:** BrightGlo<sup>®</sup> from Promega

**Assay buffer:** 20 mM HEPES, 1x HBSS, pH 7.40

**2.7.1. Cell culture.** HTLA cells (a gift from Dr. Richard Axel), stably expressing a tTA-dependent luciferase reporter and a  $\beta$ -arrestin-TEV protease fusion gene, are maintained in DMEM supplemented with 10% FBS and 2  $\mu$ g/ml Puromycin and 100  $\mu$ g/ml Hygromycin. To set up the cells for transfection, HTLA cells are plated in DMEM supplemented with 10% dialyzed FBS in Poly-L-Lys (PLL) coated 384-well white clear bottom cell culture plates at a density of 10,000 to 20,000 cells in 50  $\mu$ l per well and incubated overnight.

**2.7.2. DNA plate.** Each single DNA plasmid is plated using the liquid handling workstation into one well of a 96-well plate at 0.5  $\mu$ g/well (enough for 8x 384-well plates). Each plate includes 80 receptor DNA samples, positive controls in wells A12 and B12 (DNA plasmids for receptors with cognate ligands), transfection controls in wells A1, B1, G12, and H12 (DNA plasmid for YFP), and negative controls with buffer only. We use D<sub>2</sub>-Tango (quinpirole as agonist) and V<sub>2</sub>-Tango (Vasopressin as agonist) constructs as positive assay controls. These plates are kept in the freezer until ready for assay. Immediately before transfection (see below), two DNA plates are manually combined in a cell culture hood into one 384-well plate with duplicate wells for every DNA (see **Figures 46 and 47** for DNA maps in 96- and 384-well plate formats); immediately followed by calcium phosphate transfection (see below).

**2.7.3. Transfection using Calcium phosphate precipitation protocol.** HTLA cells are set up as indicated above overnight before transfection. Plated DNA in 384-well plate (**Figure 47**) is first diluted to final volume of 100  $\mu$ l with 0.25 M CaCl<sub>2</sub> in TE buffer (1 mM Tris-HCl, 0.1 mM EDTA, pH 7.6). An equal volume of 2x HBS solution (50 mM HEPES, 280 mM NaCl, 10 mM KCl, 1.5 mM Na<sub>2</sub>HPO<sub>4</sub>, pH 7.00) is added to the DNA/CaCl<sub>2</sub> solution and dispensed onto the plated cells in duplicate (adjacent wells) as shown in **Figure 47** (i.e. A1 and A2, A3 and A4, ..., P23 and P24). Thus, the DNAs from two 96-well plates (p1 and p2) are transfected in one 384-well plate (P1),

and therefore each DNA plasmid is transfected in duplicate. DNA plate 1 (p1) is transfected in rows A to H and DNA plate 2 (p2) is transfected in rows I to P. Plates are incubated overnight at 37°C.


**2.7.4. Assay procedure.** Each 384-well plate is designed to test one drug at multiple (up to 160) target receptors simultaneously. Odd-numbered columns contain Tango assay buffer to serve as a basal control and even-numbered columns contain drug stimulation solution prepared in filtered Tango assay buffer at 6  $\mu\text{M}$  (final concentration is 1  $\mu\text{M}$ ). The drug plate design is shown in **Figure 48**. On the day of the assay, growth medium is replaced by serum-free medium supplemented with penicillin/streptomycin 4 hours before stimulation. Cells are then stimulated by addition of drugs (10  $\mu\text{l}$  per well) and incubated overnight at 37°C. The following day, medium and drug solutions are removed and 20  $\mu\text{l}$  per well of BrightGlo reagent (diluted by 20-fold with Tango assay buffer) are added. The plate is incubated for at least 20 minutes at room temperature in the dark before luminescence is measured.

**2.7.5. Data processing and analysis.** The luminescence counter records relative luminescence units (RLU) and save files in Excel spreadsheets for easy processing. For receptors that have positive controls (non-orphan receptors), Activation relative to positive control (%) is calculated according to the following formula:

$$\begin{aligned} & \textit{Activation (relative to positive control, \%)} \\ &= \frac{(\text{Test compound RLU}) - (\text{ave. negative control RLU})}{(\text{ave. positive control RLU}) - (\text{ave. negative control RLU})} * 100 \end{aligned}$$

For receptors that lack positive controls (orphan receptors), Activation relative to baseline (%) is calculated according to the formula below:

$$\textit{Activation (relative to baseline, \%)} = \frac{(\text{Test compound RLU})}{(\text{ave. baseline RLU})} * 100$$




**Figure 46.** 96-well DNA map. Green wells contain a YFP plasmid for transfection control. Red wells contain non-orphan receptors for use as positive assay controls (D2-Tango and V2-tango constructs). Grey wells have DNA plasmids. White wells contain assay buffer only for use as a negative control.

|           |   | 1     | 2     | 3    | 4    | 5    | 6    | 7    | 8    | 9    | 10   | 11   | 12   | 13   | 14   | 15   | 16   | 17   | 18   | 19   | 20   | 21   | 22   | 23   | 24   |       |
|-----------|---|-------|-------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|-------|
| DNA<br>p1 | A | Green | Green | Blue | Red  | Red   |
|           | B | Green | Green | Blue | Red  | Red   |
|           | C | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | D | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | E | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | F | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | G | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Green |
|           | H | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Green |
| DNA<br>p2 | I | Green | Green | Blue  |
|           | J | Green | Green | Blue  |
|           | K | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | L | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | M | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | N | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue  |
|           | O | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Green |
|           | P | Blue  | Blue  | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Blue | Green |

**Figure 47.** 384-well DNA plate for calcium precipitation and transfection (P1). DNA constructs are transfected in neighboring wells. Green wells are YFP-transfected wells for transfection controls. Red wells are assay controls with cells transfected with non-orphan receptors and stimulated with their cognate ligand. DNA from the 96-well plate (p1) is used to transfect rows A to H. DNA from a second plate (p2) is used to transfect rows I to P.

|   | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 |  |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|--|
| A |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| B |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| C |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| D |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| E |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| F |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| G |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| H |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| I |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| J |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| K |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| L |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| M |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| N |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| O |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |
| P |   |   |   |   |   |   |   |   |   |    |    |    |    |    |    |    |    |    |    |    |    |    |    |    |  |


**Figure 48.** 384-well drug plate design for stimulation. White wells are baseline control with assay buffer; purple wells contain drug solution (6  $\mu\text{M}$ ). Columns 1 and 2 serve as negative controls to determine drug effects on non-transfected cells. Red wells are assay controls and are thus stimulated with cognate ligands (Quinpirole for  $\text{D}_2$  and Vasopressin for  $\text{V}_2$ ).



**Figure 49.** A screen capture of the Navigator interface. This software calculates Activation (%) relative to a positive control as well as relative to baseline. Compound structures, numbers of plates and compounds tested per assay, and calculated values for every receptor/compound pair can be viewed on the main screen.

To calculate and report these values, we use the Navigator software (custom made and developed in-house by the Molecular Modeling Laboratory, College of Pharmacy, UNC at Chapel Hill). The Navigator takes raw output files in Excel spreadsheets and calculates relative activation as indicated above. A screen shot of the Navigator interface is shown in **Figure 49**.

**2.7.6. Representative figures.** As a proof of concept, we screened three compounds (LSD, Rimonabant, and Fluoxetine) against 143 GPCRs (non-orphan, non-olfactory GPCRs) at a final concentration of 1  $\mu$ M. Results (fold of baseline) are reported in **Figure 50**.



**Figure 50.** Sample parallel primary screening at non-orphan non-olfactory GPCRs using the  $\beta$ -arrestin recruitment assay. Activity is reported as Response (fold stimulation over baseline) for LSD (blue), rimonabant (red), and fluoxetine (green).

## 2.8. Thallium flux (FluxOR) assays for hERG activity

**Main equipment:** FLIPR<sup>TETRA</sup> from Molecular Devices (Sunnyvale, CA)

**Main reagents:** FluxOR kit from Invitrogen (Carlsbad, CA)

**FluxOR assay buffer:** 20 mM HEPES, 1x HBSS, 2.5 mM Probenecid, pH 7.40

The following protocol was adapted from Huang et al., (2010) *ASSAY and Drug Development Technologies*, 6:727-742.

**2.8.1. Cell culture:** HEK293 cells stably expressing hERG channels are purchased from ChanTest (Cleveland, OH) and maintained accordingly in DMEM supplemented with 10% FBS and 500 µg/ml G418. The hERG HEK293 cells are subcultured when they reached 80-90% confluency.


**2.8.2. FluxOR assays for hERG modulators:** Thallium (Tl<sup>+</sup>) Flux assays are carried out using the FluxOR Potassium Ion Channel Assay kit (Invitrogen) according to the manufacturer's instructions with modifications. In detail, hERG HEK293 cells are plated into PLL coated 384-well black clear bottom cell culture plates with DMEM supplemented with 1% dialyzed FBS at a density of 15,000 cells in a final volume of 40 µl per well. The plated cells are incubated overnight before use. On the day of the assay, FluxOR dye reagents are reconstituted by mixing Component A (1/1000 dilution) and PowerLoad (1/100 dilution) using the FluxOR assay buffer and loaded onto cells with 20 µl of FluxOR reagent for 90 min in the dark. During the incubation, drug and stimulation solutions are prepared. The stimulation solution is composed as follows: 2.5ml deionized water, 1ml FluxOR chloride-free buffer (Component E), 1 ml K<sub>2</sub>SO<sub>4</sub> (125 mM, Component F), and 0.5 ml Tl<sub>2</sub>SO<sub>4</sub> (50 mM, Component G). At the end of the dye loading period, the dye is removed, and the FLIPR<sup>TETRA</sup> is programmed to transfer drugs from drug plates into cell plates (25 µl per well). Following transfer, cells are incubated with the drugs for 15 min at room temperature in the dark. Then, stimulation solution (6.3 µl per well) is added with FLIPR<sup>TETRA</sup>. Before the drug addition, the fluorescence intensity (excitation at 490 nm and emission at 525 nm) in each well is measured every second for 10 seconds to establish a baseline, and after addition, for additional 90 seconds. Alternatively, the FluxOR assays can be performed with

cryopreserved hERG HEK293 cells. In brief, frozen cells are washed with growth media once to remove the DMSO in the freezing medium, and then plated at 20,000 cells per well as above and assayed 5 hours later. Results using cryopreserved cells are not different from those obtained using fresh cells.


**2.8.3. Assay procedure for hERG trafficking modulators:** The above protocol is designed to measure the acute effect of drugs on hERG channel activity (such as hERG channel inhibitors or activators). For drugs without acute inhibitory effects in the  $TI^+$  flux assays, we have modified the standard protocol to conduct a longer “chronic” study to identify compounds that might be acting through indirect mechanism(s) (e.g., hERG trafficking inhibition, hERG internalization, and  $Na^+-K^+$  ATPase inhibition). Briefly, cells are first treated with drugs for the desired time (up to 16 hours), and are washed once with assay buffer before dye loading. The dye solution is then replaced with assay buffer and the fluorescence intensity is measured upon addition of stimulation solution as in the standard protocol.

**2.8.4. Data processing and analysis:** The fluorescence intensity time course of each well in the 384-well assay plate is processed using ScreenWorks 2.0 to export the slope of curve for the first 15 seconds after addition of the stimulation buffer. The initial slopes are then plotted against drug concentrations and analyzed in Prism 5.0 as outlined in **Section 2.3**. To normalize, the average initial slope with buffer is set as 0% inhibition and the average initial slope with 10  $\mu$ M Cisapride is set as 100% inhibition.


**Figure 51.** Raw fluorescence signal traces recorded on the FLIPR with HEK293 cells stably expressing hERG channels. Recording started 10 seconds before addition of stimulation solution and continued for another 90 seconds. Cisapride reduced signals in a concentration-dependent manner (error bars not shown). Not all concentrations are shown.




**Figure 52.** The Initial slope of curves from **Figure 51** was exported from FLIPR ScreenWorks, plotted against corresponding concentrations of cisapride, and fitted to a four-parameter logistic function using GraphPad Prism 5.0.



**Figure 53.** Acute and chronic effects of the hERG trafficking inhibitor ouabain on hERG channel activity. The activity of ouabain was determined with acute  $TI^+$  flux assay (upper) and chronic  $TI^+$  assay (lower). HEK293 cells stably expressing hERG channels were first incubated with ouabain for desired times at 37°C; the remaining hERG activity was determined.



**Figure 54.** Activation of hERG channels in  $TI^+$  flux assays. Terfenadine completely inhibited hERG channel activity, whereas PD-118057 (upper panel) and Ciclopirox (lower panel) activated the hERG channel. For assays done in the presence of terfenadine, terfenadine was added 5 min before PD-118057 or ciclopirox. For hERG activator activity, data were normalized to their corresponding basal activities in percentage: basal activity as 100% in the absence of terfenadine and basal activity as 0% in the presence of terfenadine.



## 2.9. PatchXpress assays for hERG activity

**Main equipment:** PatchXpress 7000A (MDS Analytical Technologies, Sunnyvale, CA)

**PatchXpress external solution:** 137mM NaCl, 4mM KCl, 1.8mM CaCl<sub>2</sub>, 1mM MgCl<sub>2</sub>, 10mM HEPES, 10mM Glucose, pH7.4 by NaOH.

**PatchXpress internal buffer:** 15mM NaCl, 70mM KF, 60mM KCl, 1mM MgCl<sub>2</sub>, 5mM HEPES, 5mM EGTA, 4mM ATP, and 0.4mM GTP, pH 7.2 by KOH.

**Automated planar patch clamp (APPC):** PatchXpress electrophysiology


The following protocol is adapted from Huang et al., (2010) ASSAY and Drug Development Technologies, 6:727-742.

**2.9.1. External and Internal solutions:** Fresh external and internal solutions are prepared at room temperature on the day of patch clamp assays. The KF-based internal buffer was adopted from a report by Zeng et al., (2008), in which they reported a greater overall success rate with a KF-based internal buffer than with a traditional KCl-internal solution. We added ATP, GTP, and Mg<sup>2+</sup> to the KF-based internal solution to prevent potential current run-down. The osmolality of the buffers was determined with a VAPRO 5520 Vapor Pressure Osmometer (Wescor, INC, Logan UT). The osmolality is usually at 285 ± 10 mmol/Kg for the external buffer and 295 ± 10 mmol/Kg for the internal buffer. The buffers are vacuum filtered to remove any air bubbles or small particles.

**2.9.2. Cell preparation:** For patch clamp assays, hERG-expressing HEK 293 cells are maintained as described above and subcultured into 10-cm dishes two days before the scheduled assays at 1 – 2 million cells per dish in growth media without G418. Cells are not used if they reach more than 90% confluency. To prepare the cells for patch clamping, we follow ChanTest's recommendations with minor modifications. The goal is to prepare a clean, fresh cell suspension immediately before loading cells onto the PatchXpress. Cells are briefly washed with PBS, treated with Accutase (Sigma, 2.5ml per 10-cm dish) for 4 min at room temperature to detach

them, gently transferred and suspended in 20 ml growth media without G418 in a 50-ml centrifuge tube, and allowed to recover from detachment for 30 min at 37°C. At the end of the incubation period, 1 million cells are transferred into another 50ml centrifuge tube and pelleted by centrifugation at 250 x g for 2.5 min at room temperature. The cell pellet is then gently re-suspended into 170 µl of the external buffer, transferred into a 1.5ml microcentrifuge tube, and loaded into the PatchXpress 7000A (MDS Analytical Technologies, Sunnyvale, CA ) in 'waiting mode' for cells. To minimize delay and cell clumping, the APPC system is started well in advance of use and primed with fresh external and internal solutions before preparing cells. The patch clamping procedure is started with a new SealChip immediately after the end of the 30-min incubation and recovery period. By the time the cell suspension is ready to load into the APPC system, the system has then reached the 'waiting mode' for receiving the cells.

**2.9.3. Drug plate preparation:** While cells are in the incubation and recovery period, drug solutions are prepared in 0.5ml polypropylene round-bottom 96-well drug plates (Fisher Scientific). Drugs in 10 mM DMSO stocks are diluted in external buffer at a final volume of 360 µl each well, a sufficient quantity for the drug to be tested in two cells (triple addition of 50 µl for one cell at each concentration). The final concentration of DMSO is 0.3% (v/v) for all dilutions (except for occasional drugs that require 3% DMSO as indicated in Results). For initial assays, 8-point (ranging from 30 nM to 30 µM with 0.3% DMSO in external buffer as a negative control) concentration-response curves are generated. For subsequent assays, the concentration range is adjusted as necessary to give full concentration-response curves (0 to 100% inhibition). The drug plate setup information is manually entered into the APPC system before starting the procedure. The assay is usually started within 15 minutes after drug plate is prepared, and the highest concentration is tested within 90 minutes. A new drug plate is prepared for each new SealChip.



**Figure 55.** A typical APPC drug plate map (minimum of 360  $\mu$ l per well). Each drug has eight serial dilutions 0.5 log units apart as indicated in the plate map, starting with a buffer control on the left and ending at 30  $\mu$ M, followed by a positive control (**red**) on the right, usually Cisapride at 3  $\mu$ M.



**2.9.4. APPC procedures:** The APPC procedure starts with manual loading of a SealChip<sub>16</sub><sup>TM</sup> (AVIVA Biosciences, San Diego, CA) and is executed automatically at room temperature. The system requires approximately 7 min to (1) dry and load the SealChip onto the recording station; (2) add external and internal solutions; (3) check the quality of each of the 16 chambers; (4) enter waiting mode and prompt for cell loading, at which point the cell suspension prepared above is immediately loaded. After loading, the single cell suspension is triturated briefly and aliquoted (35  $\mu$ l,  $\sim$  20,000 cells) by the on-board robotic Cavro pipette to each chamber on the SealChip. First, positive pressure (6 mmHg), then brief suction with negative pressure (-45

mmHg), are applied to help the cells descend quickly, and a single cell is drawn onto the top of the electrode hole of each chamber. A negative pressure ramp to 75mmHg (at 5mmHg per second) is then applied repeatedly until the seal resistance reaches 1G $\Omega$  (Gigaohm) or greater for Giga-seals and over 300 M $\Omega$  (Megaohm) for 2<sup>nd</sup> seals. After obtaining a seal, another negative pressure ramp from -40mmHg to -250mmHg is applied repeatedly to rupture the patched membrane and achieve whole cell configuration. Chambers with low seal resistances (i.e. no cell detection) are terminated within 90 seconds; chambers that do not form seals within 5 minutes are terminated by either a built-in script or by user intervention. After achieving whole cell configuration, the voltage protocol is activated, beginning with a 2-minute wash-out and a 5-minute stabilization period. After stabilization, cells with less than 0.2 nA tail current amplitude are terminated. The voltage protocol consists of the following steps: (1) depolarization from a holding potential of -80mV to -50mV for 50 milliseconds to measure leak current without activation of the hERG channel; (2) further depolarization to +20mV for 5 seconds to activate the hERG channel (hERG channels are activated and quickly inactivated); (3) repolarization to -50mV for 1.7 seconds to remove the inactivation and elicit the outward hERG tail current; (4) repolarization to the holding potential of -80mV to keep hERG channels closed. The pulse pattern is applied repeatedly every 10 seconds (0.1 Hz). The instantaneous current at -50mV before stepping to +20mV is designated as the leak current and is subtracted from the corresponding peak tail current for data processing. Each concentration-response experiment starts with a buffer control to determine the maximal hERG tail current (0% inhibition) in the absence of drugs. Each drug concentration is applied three times (50  $\mu$ l each at 25 $\mu$ l per second) with 11 seconds between additions, and each addition is preceded by aspiration of the buffer or previous drug solution from the chamber down to a 5  $\mu$ l dead volume. This triple addition and aspiration protocol fully exchanges the solution in the chamber with negligible dilution within one minute. The built-in Database script is activated between each concentration. When the drug effect reaches a steady state (less than 0.1% difference from last measurement) or after a maximum of 5 minutes, the next concentration is queued for delivery. Immediately after each concentration-response trial, a 5-minute wash-out is applied to monitor recovery of the hERG

channels from inhibition. At the end of the wash-out, a positive control (usually Cisapride at 3  $\mu\text{M}$ ) is applied to inhibit any recovered or remaining hERG channel activity.

**2.9.5. Data processing and analysis:** APPC data are collected and automatically deposited into the database program DataXpress V2.0 (MDS Analytic Technologies, Sunnyvale, CA). If a cell has a leakcurrent that is as much as 1/3 of the total tail current, the data for that cell is excluded. If a cell shows current run-down of more than 25% of the initial total tail current at the first (lowest) drug concentration, the data for that cell is excluded. The drug concentration range is usually adjusted after an initial assay so that, in subsequent assays, the lowest drug concentration inhibits hERG by less than 10%. The hERG tail currents are transformed and normalized to percent inhibition with a built-in script in the DataXpress program (total initial tail current = 0% inhibition, no current = 100% inhibition). Normalized results from multiple assays are pooled and analyzed with Prism's built-in four-parameter logistic functions as outlined in **Section 2.3**. A built-in statistical comparison function in Prism is then used to determine if a model with a variable Hill slope fit the data better than a model with a standard Hill slope (slope of 1); a p value less than 0.05 is considered significant. Finally, the potency value and corresponding Hill slope value from the best-fit model are reported.

**Figure 56.** Voltage protocol (upper) and corresponding hERG tail current recording (lower) for the automated planar patch clamp assays.



**Figure 57.** A representative whole-cell patch clamp recording showing that the hERG channel current (nA) was severely inhibited by a PDSP compound in a concentration-dependent manner, which was partially recovered during a wash-out cycle, and then completely inhibited again by a positive control (3  $\mu$ M Cisapride). The figure is a captured screen in the DataXpress program showing hERG tail currents in the absence and presence of a PDSP compound. The red vertical line (left) indicates addition of buffer control, the green vertical lines indicate triple additions of drugs with concentrations ( $\mu$ M) above the green lines, the dark blue vertical line indicates wash-out with buffer, and the pink vertical lines indicate the positive control, Cisapride. The red crosses indicate the time when the hERG tail current was measured by the system.







Figure 58. A captured screen in the DataXpress program showing a representative whole-cell patch clamp recording of hERG tail currents in the absence and presence of cisapride. The red vertical line indicates addition of buffer control and the green vertical triple lines indicate triple additions of cisapride with concentrations ( $\mu\text{M}$ ) listed above the green lines. The red crosses indicate the times when the hERG tail current was measured by the system.



**Figure 59.** The hERG tail currents in the **Figure 58** were extracted from the program DataXpress, plotted against cisapride concentrations, and fitted to a four-parameter logistic function using GraphPad Prism 5.0.



**Figure 60.** Comparison of concentration-dependent activity of cisapride in  $^3\text{H}$ -dofetilide competition binding assay;  $\text{TI}^+$  flux assay, and APPC assay. Results from multiple assays ( $n \geq 2$ ) are normalized to percentage inhibition and pooled for curve-fitting in GraphPad Prism 5.0. The Results indicate that the PatchXpress assay is the most sensitive and the  $\text{TI}^+$  flux assay was the least sensitive.



## **2.10. Neurotransmitter transporter assays for DAT (Dopamine transporter), NET (norepinephrine transporter), and SERT (Serotonin transporter).**

**Main equipment:** FlexStation II (Molecular Devices, Sunnyvale, CA)

**Main reagent:** Neurotransmitter transporter uptake assay kit (R8174) (Molecular Devices)

**Assay buffer:** 20 mM HEPES, 1x HBSS, pH 7.40, RT



**2.10.1. Background and principle:** The assay system is designed to use the same fluorophore to measure norepinephrine (NET), dopamine (DAT), and serotonin (SERT) transporter activity. The proprietary fluorophore mimics biogenic neurotransmitters and is actively transported into the cell through the NET, DAT, or SERT transporters. After incubation with test compounds, the dye solution is added to the cells and the fluorescent dye is transported into the cell. External fluorescence is extinguished with a masking dye, which cannot enter cells. Therefore, the fluorophore fluoresces when it enters the cell, and the fluorescence intensity is proportional to transporter activity. The assay can be performed without a washing step and fluorescence intensity can be monitored in kinetic mode or end-point modes.


**2.10.2. Assay procedure:** Neurotransmitter transporter assays are conducted using Molecular Devices' Neurotransmitter Transporter Uptake Assay Kit (R8174) with HEK293 cells stably expressing human DAT, NET, or SERT. In brief, cells are plated in Poly-L-Lys (PLL) coated 384-well black clear bottom cell culture plates in DMEM + 1% dialyzed FBS, at a density of 15,000 cells per well in a total volume of 40  $\mu$ l. The cells are incubated for a minimum of 6 hours before being used for assays. Medium is removed, and 20  $\mu$ l of assay buffer (20 mM HEPES, 1x HBSS, pH 7.40) is added, followed by 5  $\mu$ l of 5x drug solutions (384-well drug plate map #1 or #2, **Figure 22**, for primary assays and see 384-well drug map #4, **Figure 25**, in **Section 2.1**). The plate is incubated at 37°C for 30 min. After incubation, 25  $\mu$ l of dye solution is added and the fluorescence intensity is measured after 30 min at 37°C, using the FlexStation II (bottom read mode, Excitation at 440 nm, Emission at 520 nm with 510 nm cut-off). Results (in Relative Fluorescence Units, RLU) are exported and plotted against drug concentrations in Prism 5.0 for nonlinear regression

to obtain inhibitory potency. Nisoxetine and Cocaine serve as positive controls for DAT, NET, and SERT.

**2.10.3. Data processing and analysis.** Fluorescence intensity values are exported and analyzed in Prism 5.0 to obtain  $IC_{50}$  values using non-linear least-squares curve fitting.

**2.10.4. Representative figures**

**Figure 61.** Representative curves of neurotransmitter reuptake inhibition at DAT, NET, and SERT.



### 2.11. Multidrug Resistance Transporter (MDR-1) assay.

**Main equipment:** FlexStation II (Molecular Devices, Sunnyvale, CA)

**Assay buffer:** Dulbecco's PBS, 10 mM Glucose

**Protocol:** The MDR assay protocol is adapted from PubChem BioAssay ID 377 (<http://pubchem.ncbi.nlm.nih.gov/assay/assay.cgi?aid=377>)

**2.11.1. MDR assay Background:** Assays for modulation of MDR activity are performed using Caco-2 cells, a cultured line derived from human colonic epithelium, or HEK human kidney cells that express MDR. The assay relies on calcein acetoxymethyl ester (calcein-AM), a lipophilic compound that enters cells by passive diffusion across the plasma membrane. Once inside cells, esterases can hydrolyze calcein-AM to calcein, which is trapped in the cytoplasm because it is negatively charged. Importantly, calcein is highly fluorescent, but calcein-AM is not.

The assay is based upon the principle that calcein-AM, but not free calcein, is transported out of cells by MDR. In cells that express MDR, the lipophilic calcein-AM that crosses the plasma membrane will be pumped out of the cells before it can be converted to calcein, and thus the accumulation of fluorescence will be low. Compounds that compete with calcein-AM for MDR will prevent calcein-AM transport, resulting in increased conversion to calcein and increased fluorescence. Therefore, the effect of a compound on fluorescence accumulation can be used to measure its interaction with the MDR protein.

**2.11.2. MDR assay procedure:** The assay monitors the time-dependent increase in calcein fluorescence in live cells in 96 well plates. This is carried out using a FlexStation II fluorimeter (Molecular Devices). Cells are seeded into 96-well clear bottom cell culture plates one day before assay (80,000 cells per well). On the day of the assay, the medium is removed and replaced with 50  $\mu$ l of D-PBS, 10 mM glucose containing no additional compound (negative control), test compound (25  $\mu$ M), or reference compound (cyclosporin A) (25  $\mu$ M). The cells are incubated for 30 min at 37°C, and then the instrument adds calcein-AM to the cells (500 nM

final concentration). The instrument monitors fluorescence over a 4-min period and calculates the slope of the increase in fluorescence. All compounds are assayed in quadruplicate and each assay contains wells with no test compound (negative control) and wells with 25  $\mu$ M cyclosporin A (positive control), an efficient MDR inhibitor. Results for test compounds are calculated from the slope of the fluorescence increase and are normalized so that the value for untreated cells is 0% and the value for cyclosporin A is 100%.


This assay has several features that make it ideal for initial screening of compounds for interaction with MDR. Because the assay is carried out in live cells, compounds must diffuse across lipid bilayers to interact with MDR sites on the cytoplasmic face of the protein. This is similar to the situation in vivo, where compounds must diffuse into the cytoplasm where they interact with MDR. Similarly, the assay provides a means for assessing not only interactions with MDR but also partitioning across cell membranes and thus hydrophobicity.

Although this assay is excellent for initial screening, users should be aware that the assay has several drawbacks. i) The assay does not distinguish between MDR substrates and inhibitors. Both will give similar signals in the assay because they prevent the transport of calcein-AM. ii) Some compounds may give spurious results by inhibiting the esterases that convert calcein-AM to calcein. Finally (iii) activity depends on the cytoplasmic concentration of the compounds. For a compound that is an MDR substrate, this concentration depends on the rate of diffusion across the plasma membrane and the rate at which MDR pumps the compound from the cells. At steady state, the cytoplasmic concentration will be lower than the extracellular concentration, but it cannot be measured easily. Consequently this assay is not the best choice for determining half-maximal concentrations for interacting compounds.

**2.11.3. Data process.** Fluorescence intensity is exported and analyzed in Prism 5.0 using non-linear least-squares curve fitting as above.

#### **2.11.4. Representative figures**

Figure 62. Representative data from an MDR inhibition assay.



## 2.12. Histone Deacetylase (HDAC) inhibition assay.

**Main equipment:** FlexStation II (Fluorescence plate reader) (Molecular Devices, CA)

**Assay buffer:** 50 mM Tris HCl, 137 mM NaCl; 2.7 mM KCl, 1 mM MgCl<sub>2</sub>, pH 8.0



**Main reagent:** The HDAC assay protocol was adapted from BioMol Fluor de Lys assay system (Plymouth Meeting, PA)

**2.12.1. General Assay Procedure.** To identify potential inhibitors of HDAC, we utilize the fluorimetric Fluor de Lys HDAC Assay Kit (Biomol) as instructed by the manufacturer. Briefly, 4X dilutions of test compound or reference compound (trichostatin A) are prepared (final assay concentrations span from 0.1 nM to 10  $\mu$ M) in Assay Buffer and 12.5  $\mu$ l are added to appropriate wells of a proprietary 96-well plate (Biomol). Nuclear extracts containing HDAC activity (Biomol) are diluted to 4X and 12.5  $\mu$ l are added to the wells containing test or reference compound (each concentration is assayed in triplicate). The samples are incubated at room temperature for 10 min to equilibrate the temperature. Then, 25  $\mu$ l of 2X Fluor de Lys HDAC substrate (final HDAC substrate concentration is typically a value between one half its apparent  $K_M$  and the apparent  $K_M$ ; for HDAC1 a concentration of 50  $\mu$ M is used, for HDAC6 a concentration between 10 and 30  $\mu$ M is used) are added to each well. Deacetylation of the substrate, which generates a product that can be made fluorescent, is allowed to proceed for 30 min. Next, the reactions are stopped and the fluorescence of the deacetylated product is developed by adding 50  $\mu$ l of 2X Assay Developer and incubating at room temperature for 15 min. Finally, fluorescence is read on a FLEXStation II plate reader (Molecular Devices) (excitation wavelength 350-380 nm, emission wavelength 440-460 nm).

**2.12.2. Data processing.** Raw data (RFUs) representing fluorescence of the deacetylated substrate are plotted as a function of the logarithm of the molar concentration of the test or reference compound. Non-linear regression of the raw data (normalized to the fluorescence measured in the absence of HDAC inhibitor and test compound) is performed in Prism 5.0 using the built-in three parameter logistic model describing competitive inhibition (one-site).

### 2.12.3. Representative figures.

**Figure 63.** Representative HDAC 1 inhibition curves (upper panel) and HDAC 6 inhibition in the absence and presence of compounds at 10  $\mu$ M (lower panel).



### 2.13. Monoamine Oxidase (MAO) A and B assays.

**Main equipment:** FlexStation II (Fluorescence plate reader) (Molecular Devices, CA)

**Main reagent:** Monoamine A and B Detection Kit from Cell Technology (Mountain View, CA)

**2.13.1. Brief background.** Monoamine oxidases (MAO) are flavin-containing enzymes that catalyse the oxidation of amine-containing neurotransmitters such as serotonin, norepinephrine, epinephrine, and dopamine to yield the corresponding aldehydes. MAO has two isoforms, MAO-A and MAO-B. They exhibit different specificities to substrates and inhibitor selectivity. MAO-A acts preferentially on serotonin and norepinephrine and is inhibited by clorgyline. MAO-B acts preferentially on 2-phenylethylamine and benzylamine and is inhibited by deprenyl and pargyline.



**2.13.2. Assay principle.** We use the MAO-A and -B Detection Kit (Cell Technology) for Monoamine oxidase-A and -B enzyme assays. The assay system utilizes a non-fluorescent proprietary substrate to detect H<sub>2</sub>O<sub>2</sub> released by the conversion of a MAO substrate (Benzylamine for MAO A and Tyramine for MAO B) to its aldehyde. The reaction of H<sub>2</sub>O<sub>2</sub> with the non-fluorescent substrate is catalyzed by Peroxidase at a 1:1 stoichiometry to produce a fluorescent product with an emission wavelength of 590 – 600 nm and an excitation wavelength of 530 – 571 nm.

**2.13.3. Assay procedure.** MAO-A and -B assays are performed in 96-well plates according to the manufacturer's instructions. In brief, samples (test drugs and controls) are added to designated wells, followed by a reaction cocktail containing reaction buffer, non-fluorescent substrate for H<sub>2</sub>O<sub>2</sub>, HRP, and a substrate for MAO-A or -B. The reactions are allowed to proceed for 30 minutes at room temperature in the dark. After 30 min incubation, plates are read in the FlexStation II (Molecular Devices, CA) using an excitation wavelength of 570 nm and an emission wavelength of 590 nm.

**2.13.4. Data processing.** Fluorescence intensity values are exported and analyzed in Prism v 5.0 using non-linear least-squares curve fitting as above.

#### 2.13.5. Representative figures

**Figure 64.** Representative curves for MAO-A activity in the presence of Clorgyline and Pargyline (upper panel) and bar graph showing MAO-A activity in the presence of various test compounds at 10  $\mu$ M (lower panel).



## 2.14. Protein Kinase C activity assay

**Main equipment:** FlexStation II (Fluorescence plate reader) (Molecular Devices, CA)

**Main reagent:** Omnia Ser/Thr Recombinant Kit 8 (#KNZ2081) from Invitrogen (Carlsbad, CA)

**PKC isoforms:** purchased from Invitrogen or Sigma (St. Louis, MO)

**2.14.1. Brief background.** Protein Kinase C (PKC) has 12 isoforms and is classified into three groups (conventional PKCs, Novel PKCs, and Atypical PKCs) based on their requirement for activators (Calcium ions and/or lipids). Conventional PKC isoforms ( $\alpha$ ,  $\beta$ I,  $\beta$ II, and  $\gamma$ ) require Calcium, DAG, and phospholipid as activators; novel PKC isoforms ( $\delta$ ,  $\epsilon$ ,  $\eta$ ,  $\theta$ , and  $\mu$ ) require DAG but not Calcium; atypical PKC isoforms ( $\zeta$  and  $\iota$ ) require neither DAG nor Calcium. We use Invitrogen's PKC assay kit (Kinase Activity Assay Kit, KNZ2081, aka The Omnia(R) Ser/Thr Recombinant Kit 8) for PKC inhibitor screening assays.


**2.14.2. Principle of the PKC assay.** The assay uses a Ser/Thr containing peptide substrate conjugated with the chelation-enhanced fluorophore (CHEF) 8-hydroxy-5-(N,N-dimethylsulfonamido)-2-methylquinoline (Sox). Phosphorylation of the peptide substrate results in  $Mg^{2+}$  chelation and formation of an ion-ion interaction bridge between the Sox moiety and the phosphate group, leading to an increase in fluorescence at 485 nm when being excited at 360 nm.



**2.14.3. Assay procedure.** The PKC kinase activity assays are performed in 96-well plates (1/2 area and low protein binding surface) according to the manufacturer's suggested procedures. In brief, a master reaction mix is made containing the following components: assay buffer, peptide substrate, ATP, DTT, Calcium and/or lipid activator (use buffer for atypical PKC isoforms). This mixture is aliquoted to appropriate wells, drug working solutions (samples) or buffer (as negative control) or known inhibitors (as positive controls) are added, and the mixture is warmed in the FlexStation to 30°C for about 10 min. The reactions start when the PKC isoform is added, and plates are read every minute for 60 min, with an excitation wavelength of 360 nm and an emission wavelength of 485 nm.



**2.14.4. Data analysis.** Fluorescence intensity increases over time. The intensity values are exported when the total activity reaches a plateau (or at 30 min) and are analyzed in Prism v 5.0 using non-linear least-squares curve fitting as above.



**2.14.5. Representative figures.**


**Figure 65.** Time course of PKC activity in the absence and presence of various compounds (10  $\mu$ M).











## 2.15. Checkpoint Kinase 2 (CHK2) assay

**Main equipment:** FlexStation II (Fluorescence plate reader) (Molecular Devices, Sunnyvale, CA)

**Main reagent:** Omnia Ser/Thr Recombinant Kit 3 (#KNZ1031) from Invitrogen (Carlsbad, CA)

CHK2 enzyme (#PV3367) is also purchased from Invitrogen.

**2.15.1. Brief background.** CHK2 is one of the Ser/Thr kinases phosphorylated and activated by upstream signaling apparatus (ATM and ATR) in response to DNA damage. Along with CHK1, these kinases play a critical role in determining cellular responses to DNA damage. Inhibitors for kinases like CHK2 may represent novel anticancer therapies.



**2.15.2. Principle of the CHK2 kinase assay.** We use Invitrogen's Kinase Activity Assay Kit (#KNZ1031, Omnia Ser/Thr Recombinant Kit 3) to measure CHK2 activity. The Omnia<sup>®</sup> Kinase assay uses a Ser/Thr-containing peptide substrate conjugated with the chelation-enhanced fluorophore (CHEF) 8-hydro-5-(N,N-dimethylsulfonamido)-2-methylquinoline (Sox). Phosphorylation of the peptide substrate results in Mg<sup>2+</sup> chelation and formation of an ion-ion interaction bridge between the Sox moiety and phosphate group, leading to an increase in fluorescence at 485 nm when being excited at 360 nm.

**2.15.3. Assay procedure.** CHK2 kinase activity assays are performed in 96-well plates according to the manufacturer's suggested procedure. In brief, a master reaction mix is made with following components: kinase assay buffer, peptide substrate, ATP, DTT, and samples or controls; the mixture and the empty assay plate are warmed for 5 min at 30°C. To start the reaction, kinase solution is added first to the warmed plate, followed by the master mix to appropriate wells. The plate is incubated at 30°C. During the incubation, the plate is read every minute for 60 min, with excitation wavelength of 360 nm and emission wavelength of 485 nm.

**2.15.4. Data analysis.** Fluorescence intensity increases over time. The intensity values are exported when the total activity reaches a plateau (or at 30 min) and are analyzed in Prism v 5.0 using non-linear least-squares curve fitting as above.

### 2.15.5. Representative figures

**Figure 66.** Representative figures for CHK2 kinase activity kinetics (upper panel) and percentage inhibition (lower panel).



## 2.16. Functional assays with human $\alpha 3\beta 4$ and $\alpha 4\beta 2$ nAChRs – $^{86}\text{Rb}^+$ efflux assay

**2.16.1. Brief background.** Nicotinic acetylcholine receptors (nAChRs) are ligand-gated ion channels, that differ from G protein-coupled muscarinic acetylcholine receptors (mAChRs). Each nAChR contains five subunits symmetrically forming a pore; each subunit has four transmembrane domains with both intracellular N- and C-termini.

**2.16.2. Cell culture.** Cells for nAChR functional assays are maintained as in [Section 1.5.1](#).

**2.16.3. General assay procedure.** Agonist and antagonist activities of PDSP compounds on nAChRs are assessed by measuring  $^{86}\text{Rb}^+$  efflux in HEK293 cells stably expressing nAChRs, as described previously (Xiao et al. 1998; Xiao et al., 2006). In brief, aliquots of cells in selection growth medium are plated into Poly-D-Lysine coated 24-well plates. The plated cells are grown at 37°C for 18 to 24 h to reach 70-95% confluence. The cells are then incubated in growth medium (0.5 ml/well) containing  $^{86}\text{RbCl}$  (2  $\mu\text{Ci/ml}$ ) for 4 h at 37°C. The loading mixture is aspirated, and the cells are washed four times with HEPES buffer (15 mM HEPES, 140 mM NaCl, 2 mM KCl, 1 mM  $\text{MgSO}_4$ , 1.8 mM  $\text{CaCl}_2$ , 11 mM Glucose, pH 7.4; 1 ml/well). One ml of buffer, with or without agonists, is then added to each well. After incubation for 2 min, the assay buffer is collected and the amount of  $^{86}\text{Rb}^+$  in the buffer is determined. Cells are lysed by adding 1 ml of 100 mM NaOH to each well, and the lysate is then collected for determination of the amount of  $^{86}\text{Rb}^+$  in the cells at the end of the efflux assay. Radioactivity of assay samples and lysates is measured by liquid scintillation counting. The total amount of  $^{86}\text{Rb}^+$  loaded (cpm) is calculated as the sum of the assay sample and the lysate of each well. The amount of  $^{86}\text{Rb}^+$  efflux is expressed as a percentage of the  $^{86}\text{Rb}^+$  loaded. “Stimulated  $^{86}\text{Rb}^+$  efflux” is defined as the difference between efflux in presence of nicotinic agonists and basal efflux measured in the absence of agonists.

**2.16.3.1. Primary functional assays –  $^{86}\text{Rb}^+$  efflux experiments.** For assessing agonist activity, 4 concentrations of a PDSP compound to be tested, 0.1, 1, 10 and 100  $\mu\text{M}$ , are applied. Agonist activity is scaled as a percentage of the stimulation by 100  $\mu\text{M}$  nicotine (100%). If a PDSP compound shows a concentration-dependent activation, and shows 25% stimulation at any concentration, it is subject to a secondary functional assay for agonist activity. For assessing


antagonist activity, 4 concentrations of a PDSP compound to be tested, 0.1, 1, 10 and 100  $\mu\text{M}$ , are applied in the presence of 100  $\mu\text{M}$  nicotine. Antagonist activity is scaled as a percentage of the inhibition of  $^{86}\text{Rb}^+$  efflux stimulated by 100  $\mu\text{M}$  nicotine. If a PDSP compound shows concentration-dependent inhibition, and shows more than 50% inhibition at 100  $\mu\text{M}$ , it is subject to a secondary functional assay for antagonist activity. All efflux assays are performed in quadruplicate. Nicotine is included in assays to define 100% agonist activity, as well as to function as a control.


**2.16.3.2. Secondary functional assays –  $^{86}\text{Rb}^+$  efflux experiments.** For assessing agonist activity, 8 concentrations of a PDSP compound to be tested are used. Agonist activity is scaled as % of the stimulation by 100  $\mu\text{M}$  nicotine (100%). For assessing antagonist activity, 8 concentrations of a PDSP compound to be tested are used in the presence of 100  $\mu\text{M}$  nicotine. Antagonist activity is scaled as a percentage of the inhibition of  $^{86}\text{Rb}^+$  efflux stimulated by 100  $\mu\text{M}$  nicotine. All efflux assays are performed in quadruplicate. Nicotine is included in the assays to define 100% agonist activity, as well as to function as a control.

**2.16.4. Data analysis and representative figures.** Raw counts of radioactivity (in cpm/well) are exported and analyzed in Prism 5.0 by nonlinear least-squares regression to estimate  $\text{EC}_{50}$  or  $\text{IC}_{50}$  values.

## References

- Xiao Y, Meyer EL, Thompson JM, Surin A, Wroblewski J, Kellar KJ (1998). Rat alpha3/beta4 subtype of neuronal nicotinic acetylcholine receptor stably expressed in a transfected cell line: pharmacology of ligand binding and function. *Mol Pharmacol* 54(2): 322-333.
- Xiao Y, Fang H, Musachio JL, Wei ZL, Chellappan SK, Kozikowski AP, Kellar KJ (2006). Sazetidide-A, a novel ligand that desensitizes alpha4beta2 nicotinic acetylcholine receptors without activating them. *Mol Pharmacol* 70(4): 1454-1460.

**Figure 67.** Representative figures of agonist activity at nAChRs

**Figure 68.** Representative figures of antagonist activity at nAChRs.

### Section 3. Master table of PDSP targets.

**Table 30. Complete list of targets and corresponding functional and radioligand binding assays available at the PDSP.**  $^3\text{H}$  or  $^{125}\text{I}$  radioligands are available for radioligand binding assays. A checkmark (√) indicates that the functional assay is available: G<sub>q</sub> - Calcium mobilization or Inositol phosphate accumulation assay; G<sub>i</sub> or G<sub>s</sub> - split luciferase cAMP biosensor assay; β-arrestin for GPCR mediated arrestin translocation assay; “Other Assay” indicates target-specific assays as detailed in the main text and indicated in the specific section; “?” indicates assays being developed, verified or optimized. The rightmost column (PDSP V1) indicates assays available under the previous PDSP contract, in which ‘B’ represents radioligand binding assays and ‘F’ represents functional assays. See individual sections for detailed information regarding binding and functional assays.

| Target<br>(IUPHAR name) | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|-------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                         | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| 5-HT receptors          |                                       |                |                                  |            |             |                    |            |
| 5-HT1A                  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT1B                  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT1D                  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT1E                  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT1F                  |                                       |                | √                                | √          |             |                    | B, F       |
| 5-HT2A                  | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| 5-HT2B                  | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| 5-HT2C (INI)            | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| 5-HT2C (VGI)            | 3H                                    | √              |                                  |            |             |                    | B, F       |
| 5-HT2C (VGV)            | 3H                                    | √              |                                  |            |             |                    |            |
| 5-HT2C (VNV)            | 3H                                    | √              |                                  | √          |             |                    |            |
| 5-HT2C (VSV)            | 3H                                    | √              |                                  | √          |             |                    |            |
| 5-HT3                   | 3H                                    |                |                                  |            |             |                    | B          |
| 5-HT4                   | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT5A                  | 3H                                    |                | ?                                | √          |             |                    | B, F       |
| 5-HT6                   | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT7A                  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| 5-HT7B                  |                                       |                |                                  |            |             |                    | B, F       |
| 5-HT7D                  |                                       |                |                                  |            |             |                    |            |
| Muscarinic receptors    |                                       |                |                                  |            |             |                    |            |
| M1                      | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| M2                      | 3H                                    |                | √                                | √          |             |                    | B, F       |
| M3                      | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| M4                      | 3H                                    |                | √                                | √          |             |                    | B, F       |

| Target<br>(IUPHAR name)  | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|--------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                          | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| M5                       | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| Muscarinic DREADDs       |                                       |                |                                  |            |             |                    |            |
| M1D                      | 3H                                    | √              |                                  |            |             |                    |            |
| M2D                      | 3H                                    |                | √                                |            |             |                    |            |
| M3D                      | 3H                                    | √              |                                  | √          |             |                    |            |
| M4D                      | 3H                                    |                | √                                |            |             |                    |            |
| M5D                      | 3H                                    | √              |                                  |            |             |                    |            |
| Nicotinic receptors      |                                       |                |                                  |            |             |                    |            |
| α2β2                     | 3H                                    |                |                                  |            | √           | Section 2.16       | B          |
| α2β4                     | 3H                                    |                |                                  |            |             |                    | B          |
| α3β2                     | 3H                                    |                |                                  |            |             |                    | B          |
| α3β4                     | 3H                                    |                |                                  |            | √           | Section 2.16       | B          |
| α4β2                     | 3H                                    |                |                                  |            | √           | Section 2.16       | B          |
| α4β4                     | 3H                                    |                |                                  |            |             |                    | B          |
| α4β2 (Rat Brain)         | 3H                                    |                |                                  |            |             |                    | B          |
| α2β2                     | 3H                                    |                |                                  |            |             |                    | B          |
| α7                       | 125I                                  |                |                                  |            |             |                    | B          |
| Adrenergic receptors     |                                       |                |                                  |            |             |                    |            |
| α1A                      | 3H                                    | √              |                                  |            |             |                    | B, F       |
| α1B                      | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| α1D                      | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| α2A                      | 3H                                    |                | √                                |            |             |                    | B, F       |
| α2B                      | 3H                                    |                | √                                | √          |             |                    | B, F       |
| α2C                      | 3H                                    |                | √                                | √          |             |                    | B, F       |
| β1                       | 125I                                  |                | √                                | √          |             |                    | B, F       |
| β2                       | 3H                                    |                | √                                | √          |             |                    | B, F       |
| β3                       | 3H                                    |                | √                                |            |             |                    | B, F       |
| Adenosine receptors      |                                       |                |                                  |            |             |                    |            |
| A1                       | 3H                                    |                | √                                | √          |             |                    | B, F       |
| A2A                      | 3H                                    |                | √                                |            |             |                    | B, F       |
| A2B                      | 3H                                    |                |                                  |            |             |                    | B, F       |
| A3                       | 3H                                    |                |                                  |            |             |                    | B, F       |
| Angiotensin II receptors |                                       |                |                                  |            |             |                    |            |
| AT1                      | 3H                                    | √              |                                  | √          |             |                    | B          |
| AT2                      | 3H                                    |                |                                  |            |             |                    | B          |
| Apelin receptor          |                                       |                |                                  |            |             |                    |            |
| Apelin                   |                                       |                |                                  | √          |             |                    |            |
| Bile Acid receptors      |                                       |                |                                  |            |             |                    |            |
| GPBA                     |                                       |                |                                  | √          |             |                    |            |
| Bombesin receptors       |                                       |                |                                  |            |             |                    |            |
| BB1                      |                                       | √              |                                  | √          |             |                    |            |
| BB2                      |                                       | √              |                                  | √          |             |                    |            |
| BB2, mouse               |                                       | √              |                                  |            |             |                    |            |

| Target<br>(IUPHAR name)                  | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|--|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|  | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| BB3                                      |                                       | √              |                                  | √          |             |                    |            |
| Bradykinin receptors                     |                                       |                |                                  |            |             |                    |            |
| B1                                       |                                       |                |                                  | √          |             |                    | F          |
| B2                                       |                                       | √              |                                  | √          |             |                    | F          |
| Calcitonin receptors                     |                                       |                |                                  |            |             |                    |            |
| CT                                       |                                       |                |                                  |            |             |                    |            |
| CT-like                                  |                                       |                |                                  |            |             | CALCR              | F          |
| Calcium sensing receptors                |                                       |                |                                  |            |             |                    |            |
| CaS                                      |                                       |                |                                  | √          |             |                    | F          |
| GPRC6                                    |                                       |                |                                  |            |             |                    |            |
| Cannabinoid receptors                    |                                       |                |                                  |            |             |                    |            |
| CB1 (rat brain)                          | 3H                                    |                |                                  | √          |             |                    | B, F       |
| CB2                                      | 3H                                    |                |                                  | √          |             |                    | B, F       |
| Chemokine receptors                      |                                       |                |                                  |            |             |                    |            |
| CCR1                                     |                                       |                |                                  |            |             |                    | F          |
| CCR2                                     |                                       |                |                                  |            |             |                    | F          |
| CCR3                                     |                                       |                |                                  |            |             |                    | F          |
| CCR4                                     |                                       |                |                                  | √          |             |                    | F          |
| CCR5                                     |                                       |                |                                  |            |             |                    | F          |
| CCR6                                     |                                       |                |                                  | √          |             |                    | F          |
| CCR7                                     |                                       |                |                                  |            |             |                    | F          |
| CCR8                                     |                                       |                |                                  |            |             |                    | F          |
| CCR9                                     |                                       |                |                                  |            |             |                    | F          |
| CCR10                                    |                                       |                |                                  |            |             |                    | F          |
| CXCR1                                    |                                       |                |                                  | √          |             |                    | F          |
| CXCR2                                    |                                       |                |                                  | √          |             |                    | F          |
| CXCR3                                    |                                       |                |                                  |            |             |                    | F          |
| CXCR4                                    |                                       |                |                                  | √          |             |                    | F          |
| CXCR5                                    |                                       |                |                                  |            |             |                    | F          |
| CXCR6                                    |                                       |                |                                  | √          |             |                    | F          |
| CXCR7                                    |                                       |                |                                  |            |             |                    | F          |
| CX3CR1                                   |                                       |                |                                  | √          |             |                    | F          |
| XCR1                                     |                                       |                |                                  |            |             |                    | F          |
| CCRL1                                    |                                       |                |                                  |            |             | CCBP2              | F          |
| CCRL2                                    |                                       |                |                                  |            |             |                    | F          |
| Cholecystokinin receptors                |                                       |                |                                  |            |             |                    |            |
| CCK1                                     |                                       |                |                                  | √          |             |                    |            |
| CCK2                                     |                                       | √              |                                  |            |             |                    |            |
| Complement peptide receptors             |                                       |                |                                  |            |             |                    |            |
| C3a                                      |                                       |                |                                  |            |             |                    | F          |
| C5a1                                     |                                       |                |                                  |            |             |                    | F          |
| C5a2                                     |                                       |                |                                  |            |             | C5L2               | F          |
| Corticotropin-releasing factor receptors |                                       |                |                                  |            |             |                    |            |
| CRF1                                     | 3H                                    |                | √                                |            |             |                    | B, F       |

| Target<br>(IUPHAR name)                | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|--|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|  | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| CRF2                                   | 3H                                    |                | √                                |            |             |                    | B, F       |
| Dopamine receptors                     |                                       |                |                                  |            |             |                    |            |
| D1                                     | 3H                                    |                | √                                | √          |             |                    |            |
| D2                                     | 3H                                    |                | √                                | √          |             |                    |            |
| D3                                     | 3H                                    |                |                                  | √          |             |                    |            |
| D4                                     | 3H                                    |                | √                                |            |             |                    |            |
| D5                                     | 3H                                    |                | √                                | √          |             |                    |            |
| Endothelin receptors                   |                                       |                |                                  |            |             |                    |            |
| ETA                                    |                                       |                |                                  | √          |             |                    | F          |
| ETB                                    |                                       |                |                                  |            |             |                    | F          |
| Estrogen receptor                      |                                       |                |                                  |            |             |                    |            |
| GPBR                                   |                                       |                |                                  |            |             |                    |            |
| Formylpeptide receptors                |                                       |                |                                  |            |             |                    |            |
| FPR1                                   |                                       |                |                                  | √          |             |                    |            |
| FPR2/ALX                               |                                       |                |                                  | √          |             | FPRL1              | F          |
| FPR3                                   |                                       |                |                                  | √          |             | FPRL2              | F          |
| Free fatty acid receptors              |                                       |                |                                  |            |             |                    |            |
| FFA1                                   |                                       | √              |                                  |            |             | GPR40              |            |
| FFA2                                   |                                       |                |                                  |            |             |                    |            |
| FFA3                                   |                                       |                |                                  |            |             |                    |            |
| Frizzled receptors                     |                                       |                |                                  |            |             |                    |            |
| FZD1                                   |                                       |                |                                  |            |             |                    |            |
| FZD2                                   |                                       |                |                                  |            |             |                    |            |
| FZD3                                   |                                       |                |                                  |            |             |                    |            |
| FZD4                                   |                                       |                |                                  |            |             |                    |            |
| FZD5                                   |                                       |                |                                  |            |             |                    |            |
| FZD6                                   |                                       |                |                                  |            |             |                    |            |
| FZD7                                   |                                       |                |                                  |            |             |                    |            |
| FZD8                                   |                                       |                |                                  |            |             |                    |            |
| FZD9                                   |                                       |                |                                  |            |             |                    |            |
| FZD10                                  |                                       |                |                                  |            |             |                    |            |
| SMO                                    | 3H                                    |                | √?                               | ?          | Gli-Luc (?) |                    |            |
| GABAA receptors (BZP = benzodiazepine) |                                       |                |                                  |            |             |                    |            |
| GABAA a1                               | 3H, BZP                               |                |                                  |            |             |                    | B          |
| GABAA a2                               | 3H, BZP                               |                |                                  |            |             |                    | B          |

| Target<br>(IUPHAR name)                   | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|---|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|   | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| GABAA a3                                  | 3H, BZP                               |                |                                  |            |             |                    | B          |
| GABAA a5                                  | 3H, BZP                               |                |                                  |            |             |                    | B          |
| GABAA a6                                  | 3H, BZP                               |                |                                  |            |             |                    | B          |
| GABAA (rat brain)                         | 3H, BZP                               |                |                                  |            |             |                    | B          |
| GABAB (rat brain)                         | 3H                                    |                |                                  |            |             |                    | B          |
| GABAB receptors                           |                                       |                |                                  |            |             |                    |            |
| GABAB1                                    |                                       |                |                                  |            |             |                    | B          |
| GABAB2                                    |                                       |                |                                  |            |             |                    |            |
| Galanin receptors                         |                                       |                |                                  |            |             |                    |            |
| GAL1                                      |                                       |                |                                  | √          |             |                    |            |
| GAL2                                      |                                       |                |                                  | √          |             |                    |            |
| GAL3                                      |                                       |                |                                  | √          |             |                    |            |
| Ghrelin receptor                          |                                       |                |                                  |            |             |                    |            |
| Ghrelin                                   | 3H                                    | √              |                                  | √          |             |                    | F          |
| Glucagon receptors                        |                                       |                |                                  |            |             |                    |            |
| GHRH                                      |                                       |                |                                  |            |             |                    | F          |
| GIP                                       |                                       |                |                                  |            |             |                    | F          |
| GLP-1                                     |                                       |                |                                  | √          |             |                    | F          |
| GLP-2                                     |                                       |                |                                  |            |             |                    | F          |
| Glucagon                                  |                                       |                |                                  |            |             | GCGR               | F          |
| Secretin                                  |                                       |                |                                  | √          |             | SCTR               | F          |
| Glycoprotein receptors                    |                                       |                |                                  |            |             |                    |            |
| FSH                                       |                                       |                |                                  |            |             |                    |            |
| LH  |                                       |                |                                  |            |             |                    |            |
| TSH                                       |                                       |                |                                  |            |             |                    |            |
| Gonadotrophin-releasing hormone receptors |                                       |                |                                  |            |             |                    |            |
| GnRH                                      |                                       |                |                                  | √          |             |                    |            |
| GnRH2                                     |                                       |                |                                  |            |             |                    |            |
| Histamine receptors                       |                                       |                |                                  |            |             |                    |            |
| H1  | 3H                                    | √              |                                  | √          |             |                    | B, F       |
| H2  | 3H                                    | √              | √                                | √          |             |                    | B, F       |
| H3  | 3H                                    |                | √                                | √          |             |                    | B, F       |
| H4  |                                       |                |                                  | √          |             |                    | B, F       |
| Hydroxycarboxylic acid receptors          |                                       |                |                                  |            |             |                    |            |
| HCA1                                      |                                       |                |                                  |            |             |                    |            |
| HCA2                                      |                                       |                |                                  |            |             |                    |            |
| HCA3                                      |                                       |                |                                  |            |             |                    |            |
| Imidazoline receptors                     |                                       |                |                                  |            |             |                    |            |
| I1 (rat brain)                            | 3H                                    |                |                                  |            |             |                    | B          |
| I2 (rat brain)                            | 3H                                    |                |                                  |            |             |                    | B          |
| Kisspeptin receptor                       |                                       |                |                                  |            |             |                    |            |
| KISS (GPR54)                              |                                       |                |                                  |            |             |                    |            |
| Leukotriene receptors                     |                                       |                |                                  |            |             |                    |            |

| Target<br>(IUPHAR name)                 | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|---|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|   | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| BLT1                                    |                                       |                |                                  | √          |             |                    | F          |
| BLT2                                    |                                       |                |                                  |            |             |                    | F          |
| CysLT1                                  |                                       |                |                                  | √          |             |                    | F          |
| CysLT2                                  |                                       |                |                                  |            |             |                    | F          |
| OXE                                     |                                       |                |                                  |            |             |                    |            |
| FRP2/ALX                                |                                       |                |                                  | √          |             |                    |            |
| Lysophospholipid (LPA) receptors        |                                       |                |                                  |            |             |                    |            |
| LPA1                                    |                                       |                |                                  | √          |             | EDG2               | F          |
| LPA2                                    |                                       |                |                                  | √          |             | EDG4               | F          |
| LPA3                                    |                                       |                |                                  |            |             | EDG7               | F          |
| LPA4                                    |                                       |                |                                  |            |             |                    | F          |
| LPA5                                    |                                       |                |                                  | √          |             |                    | F          |
| Lysophospholipid (S1P) receptors        |                                       |                |                                  |            |             |                    |            |
| S1P1                                    |                                       |                |                                  | √          |             |                    | F          |
| S1P2                                    |                                       |                |                                  | √          |             | EDG5               | F          |
| S1P3                                    |                                       |                |                                  | √          |             | EDG3               | F          |
| S1P4                                    |                                       |                |                                  |            |             | EDG6               | F          |
| S1P5                                    |                                       |                |                                  |            |             | EDG8               | F          |
| Melanin-concentrating hormone receptors |                                       |                |                                  |            |             |                    |            |
| MCH1                                    |                                       |                |                                  | √          |             |                    | F          |
| MCH2                                    |                                       |                |                                  | √          |             |                    | F          |
| Melanocortin receptors                  |                                       |                |                                  |            |             |                    |            |
| MC1                                     |                                       | √              |                                  | √          |             |                    | F          |
| MC2                                     |                                       |                |                                  |            |             |                    | F          |
| MC3                                     |                                       | √              |                                  |            |             |                    | F          |
| MC4                                     |                                       |                |                                  | √          |             |                    | F          |
| MC5                                     |                                       |                |                                  | √          |             |                    | F          |
| Melatonin receptors                     |                                       |                |                                  |            |             |                    |            |
| MT1                                     |                                       |                |                                  | √          |             |                    |            |
| MT2                                     |                                       |                |                                  | √          |             |                    |            |
| Glutamate receptors                     |                                       |                |                                  |            |             |                    |            |
| NMDA (rat brain)                        | 3H                                    |                |                                  |            |             |                    |            |
| AMPA (rat brain)                        | 3H                                    |                |                                  |            |             |                    |            |
| PCP (rat brain)                         | 3H                                    |                |                                  |            |             |                    |            |
| Kainate<br>(rat brain)                  | 3H                                    |                |                                  |            |             |                    |            |
| NR1 (rat brain)                         | 3H                                    |                |                                  |            |             |                    |            |
| NR2B (rat brain)                        | 3H                                    |                |                                  |            |             |                    |            |
| Metabotropic Glutamate receptors        |                                       |                |                                  |            |             |                    |            |
| mGluR1                                  | 3H?                                   |                |                                  |            |             |                    | F          |
| mGluR2                                  | 3H?                                   |                |                                  |            |             |                    | F          |
| mGluR3                                  | 3H?                                   |                |                                  |            |             |                    | F          |
| mGluR4                                  | 3H?                                   |                |                                  |            |             |                    | F          |
| mGluR5                                  | 3H?                                   |                |                                  |            |             |                    | F          |

| Target<br>(IUPHAR name)                | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|--|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|  | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| mGluR5<br>(rat brain)                  | 3H                                    |                |                                  |            |             |                    |            |
| mGluR6                                 | 3H?                                   |                |                                  |            |             |                    |            |
| mGluR7                                 | 3H?                                   |                |                                  |            |             |                    | F          |
| mGluR8                                 | 3H?                                   |                |                                  |            |             |                    | F          |
| Motilin receptor                       |                                       |                |                                  |            |             |                    |            |
| Motilin                                |                                       |                |                                  | √          |             |                    | F          |
| Neuromedin U receptors                 |                                       |                |                                  |            |             |                    |            |
| NMU1                                   |                                       |                |                                  | √          |             |                    | F          |
| NMU2                                   |                                       |                |                                  | √          |             |                    | F          |
| Neuropeptide FF receptors              |                                       |                |                                  |            |             |                    |            |
| NPFF1                                  |                                       |                |                                  |            |             |                    |            |
| NPFF2                                  |                                       |                |                                  |            |             |                    |            |
| Neuropeptide S receptor                |                                       |                |                                  |            |             |                    |            |
| NPS                                    |                                       |                |                                  | √          |             |                    |            |
| Neuropeptide W/Neuropeptide B receptor |                                       |                |                                  |            |             |                    |            |
| NPBW1                                  |                                       |                |                                  |            |             |                    |            |
| NPBW2                                  |                                       |                |                                  |            |             |                    |            |
| Neuropeptide Y receptors               |                                       |                |                                  |            |             |                    |            |
| Y1                                     |                                       |                |                                  | √          |             |                    | B, F       |
| Y2                                     |                                       |                |                                  | √          |             |                    | B, F       |
| Y4                                     |                                       |                |                                  | √          |             |                    | B, F       |
| Y5                                     |                                       |                |                                  |            |             |                    | B, F       |
| Y6                                     |                                       |                |                                  |            |             |                    | B, F       |
| Neurotensin receptors                  |                                       |                |                                  |            |             |                    |            |
| NTS1                                   | 3H                                    | √              |                                  |            |             |                    | F          |
| NTS2                                   | 3H                                    |                |                                  | √          |             |                    | F          |
| Opioid receptors                       |                                       |                |                                  |            |             |                    |            |
| δ (DOR)                                | 3H                                    |                | √                                | √          |             |                    | B,F        |
| κ (KOR)                                | 3H                                    |                | √                                | √          |             |                    | B,F        |
| μ (MOR)                                | 3H                                    |                | √                                | √?         |             |                    | B,F        |
| NOP                                    | 3H                                    |                | √                                | √          |             |                    | B,F        |
| Orexin receptors                       |                                       |                |                                  |            |             |                    |            |
| OX1                                    |                                       |                |                                  | √          |             |                    |            |
| OX2                                    |                                       |                |                                  | √          |             |                    | F          |
| P2Y receptors                          |                                       |                |                                  |            |             |                    |            |
| P2Y1                                   |                                       | √              |                                  | √          |             |                    | F          |
| P2Y2                                   |                                       | √              |                                  | √          |             |                    | F          |
| P2Y4                                   |                                       | √              |                                  | √          |             |                    | F          |
| P2Y6                                   |                                       | √              |                                  | √          |             |                    | F          |
| P2Y11                                  |                                       | √              |                                  | √          |             |                    | F          |
| P2Y12                                  |                                       |                |                                  | √          |             |                    |            |
| P2Y13                                  |                                       |                |                                  | √          |             |                    |            |
| P2Y14                                  |                                       |                |                                  | √          |             |                    |            |

| Target<br>(IUPHAR name)                     | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|---|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|   | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| Parathyroid hormone receptors               |                                       |                |                                  |            |             |                    |            |
| PTH1  |                                       |                |                                  | √          |             |                    | F          |
| PTH2  |                                       |                |                                  |            |             |                    | F          |
| Peptide P518 receptor                       |                                       |                |                                  |            |             |                    |            |
| QRFP  |                                       |                |                                  |            |             |                    |            |
| Platelet-activating factor receptor         |                                       |                |                                  |            |             |                    |            |
| PAF   | 3H                                    | √              |                                  | √          |             |                    |            |
| Prokineticin receptors                      |                                       |                |                                  |            |             |                    |            |
| PKR1  |                                       |                |                                  |            |             |                    |            |
| PKR2  |                                       |                |                                  |            |             |                    |            |
| Prolactin-releasing peptide receptor        |                                       |                |                                  |            |             |                    |            |
| PRRP  |                                       |                |                                  |            |             |                    |            |
| Prostanoid receptors                        |                                       |                |                                  |            |             |                    |            |
| DP1   |                                       |                |                                  | √          |             |                    |            |
| DP2   |                                       |                |                                  |            |             |                    |            |
| EP1   | 3H                                    |                |                                  | √          |             |                    | B, F       |
| EP2   | 3H                                    |                |                                  | √          |             |                    | B, F       |
| EP3   | 3H                                    |                |                                  | √          |             |                    | B, F       |
| EP4   | 3H                                    |                |                                  | √          |             |                    | B, F       |
| FP  |                                       |                |                                  | √          |             |                    |            |
| IP1   |                                       |                |                                  | √          |             |                    | B, F       |
| TP  |                                       |                |                                  |            |             |                    |            |
| Protease-activated receptors                |                                       |                |                                  |            |             |                    |            |
| PAR1  |                                       | √              |                                  |            |             |                    |            |
| PAR2  |                                       |                |                                  |            |             |                    |            |
| PAR3  |                                       |                |                                  |            |             |                    |            |
| PAR4  |                                       |                |                                  |            |             |                    |            |
| Relaxin family peptide receptors            |                                       |                |                                  |            |             |                    |            |
| RXFP1                                       |                                       |                |                                  |            |             |                    |            |
| RXFP2                                       |                                       |                |                                  |            |             |                    |            |
| RXFP3                                       |                                       |                |                                  |            |             | SALPR              | F          |
| RXFP4                                       |                                       |                |                                  |            |             |                    |            |
| Sigma receptors                             |                                       |                |                                  |            |             |                    |            |
| Sigma 1<br>(Guinea Pig)                     | 3H                                    |                |                                  |            |             |                    | B          |
| Sigma 2 (PC12)                              | 3H                                    |                |                                  |            |             |                    | B          |
| Somatostatin receptors                      |                                       |                |                                  |            |             |                    |            |
| sst1  |                                       |                |                                  | √          |             |                    | F          |
| sst2  |                                       |                |                                  | √          |             |                    | F          |
| sst3  |                                       |                |                                  | √          |             |                    | F          |
| sst4  |                                       |                |                                  | √          |             |                    | F          |
| sst5  |                                       |                | √                                | √          |             |                    | F          |
| Tachykinin receptors (Neurokinin receptors) |                                       |                |                                  |            |             |                    |            |
| NK1   |                                       | √              |                                  | √          |             |                    | F          |

| Target<br>(IUPHAR name)                 | Available assays at PDSP (new system) |                |                                  |            |                                       | Synonym<br>Or Note         | PDSP<br>V1 |
|---|---------------------------------------|----------------|----------------------------------|------------|---------------------------------------|----------------------------|------------|
|   | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay                           |                            |            |
| NK2                                     |                                       | √              |                                  | √          |                                       |                            | F          |
| NK3                                     |                                       | √              |                                  | √          |                                       |                            | F          |
| Thyrotropin-releasing hormone receptors |                                       |                |                                  |            |                                       |                            |            |
| TRH1                                    |                                       |                |                                  |            |                                       |                            |            |
| TRH2                                    |                                       |                |                                  |            |                                       |                            |            |
| Trace amine receptor                    |                                       |                |                                  |            |                                       |                            |            |
| TA1                                     |                                       |                |                                  | √          |                                       |                            |            |
| Urotensin receptor                      |                                       |                |                                  |            |                                       |                            |            |
| UT                                      |                                       |                |                                  | √          |                                       |                            | F          |
| Vasopressin and Oxytocin receptors      |                                       |                |                                  |            |                                       |                            |            |
| V1a                                     | 3H                                    | √              |                                  | √          |                                       |                            | B, F       |
| V1b                                     | 3H                                    | √              |                                  | √          |                                       |                            | B, F       |
| V2                                      | 3H                                    | √              |                                  | √          |                                       |                            | B, F       |
| OT                                      | 3H                                    | √              |                                  | √          |                                       |                            | B, F       |
| VIP and PACAP receptors                 |                                       |                |                                  |            |                                       |                            |            |
| PAC1                                    |                                       |                |                                  |            |                                       | ADCYAP1R                   | F          |
| VPAC1                                   |                                       |                |                                  | √          |                                       |                            |            |
| VPAC2                                   |                                       |                |                                  | √          |                                       |                            | F          |
| Neurotransmitter transporters           |                                       |                |                                  |            |                                       |                            |            |
| DAT                                     | 3H                                    |                |                                  |            | √                                     | Section 2.10               | B,F        |
| NET                                     | 3H                                    |                |                                  |            | √                                     | Section 2.10               | B,F        |
| SERT                                    | 3H                                    |                |                                  |            | √                                     | Section 2.10               | B,F        |
| VMAT1                                   |                                       |                |                                  |            |                                       |                            |            |
| VMAT2                                   | 3H                                    |                |                                  |            |                                       |                            |            |
| Ion channels                            |                                       |                |                                  |            |                                       |                            |            |
| hERG K+                                 | 3H                                    |                |                                  |            | √ (TI <sup>+</sup> Flux)<br>√ (Patch) | Section 2.8<br>Section 2.9 | B,F        |
| Ca <sup>2+</sup> (rat brain)            | 3H                                    |                |                                  |            |                                       |                            | B          |
| Na <sup>+</sup> site II                 | 3H                                    |                |                                  |            |                                       |                            | B          |
| Other targets                           |                                       |                |                                  |            |                                       |                            |            |
| MDR1                                    |                                       |                |                                  |            | √                                     | Section 2.11               |            |
| PKCα                                    | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCβ I                                  | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCβ II                                 | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCγ                                    | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCδ                                    | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCε                                    | 3H                                    |                |                                  |            | √                                     | Section 2.14               | B, F       |
| PKCη                                    |                                       |                |                                  |            | √                                     | Section 2.14               |            |
| PKCζ                                    |                                       |                |                                  |            | √                                     | Section 2.14               |            |
| PKCτ                                    |                                       |                |                                  |            | √                                     | Section 2.14               |            |
| PBR (rat brain)                         | 3H                                    |                |                                  |            |                                       |                            |            |
| PNR                                     |                                       |                |                                  |            |                                       |                            | F          |
| Vanilloid receptor                      |                                       |                |                                  |            |                                       | TRPs                       | F          |

| Target<br>(IUPHAR name) | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|-------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                         | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| CHK2                    |                                       |                |                                  |            |             |                    | F          |
| HDAC                    |                                       |                |                                  |            |             |                    | F          |
| MAO A                   |                                       |                |                                  |            |             |                    | F          |
| MAO B                   |                                       |                |                                  |            |             |                    | F          |
|                         | Class A orphan GPCRs                  |                |                                  |            |             |                    |            |
| CMKLR1                  |                                       |                |                                  | √          |             | CHEMERINR          | F          |
| GPR1                    |                                       |                |                                  |            |             |                    |            |
| GPR3                    |                                       |                |                                  |            |             |                    |            |
| GPR4                    |                                       |                |                                  |            |             |                    |            |
| GPR6                    |                                       |                |                                  |            |             |                    |            |
| GPR12                   |                                       |                |                                  |            |             |                    |            |
| GPR15                   |                                       |                |                                  |            |             |                    |            |
| GPR17                   |                                       |                |                                  |            |             |                    |            |
| GPR18                   |                                       |                |                                  |            |             |                    |            |
| GPR19                   |                                       |                |                                  |            |             |                    |            |
| GPR20                   |                                       |                |                                  |            |             |                    |            |
| GPR21                   |                                       |                |                                  |            |             |                    | F          |
| GPR22                   |                                       |                |                                  |            |             |                    |            |
| GPR25                   |                                       |                |                                  |            |             |                    |            |
| GPR26                   |                                       |                |                                  |            |             |                    |            |
| GPR27                   |                                       |                |                                  |            |             |                    | F          |
| GPR31                   |                                       |                |                                  |            |             |                    |            |
| GPR32                   |                                       |                |                                  |            |             |                    |            |
| GPR33                   |                                       |                |                                  |            |             |                    |            |
| GPR34                   |                                       |                |                                  |            |             |                    |            |
| GPR35                   |                                       |                |                                  |            |             |                    |            |
| GPR37                   |                                       |                |                                  |            |             | EDNRBL             | F          |
| GPR37L1                 |                                       |                |                                  |            |             |                    |            |
| GPR39                   |                                       |                |                                  |            |             |                    |            |
|                         |                                       |                |                                  |            |             |                    |            |
| GPR45                   |                                       |                |                                  |            |             |                    | F          |
| GPR50                   |                                       |                |                                  |            |             |                    |            |
| GPR52                   |                                       |                |                                  |            |             |                    |            |
| GPR55                   |                                       |                |                                  |            |             |                    |            |
| GPR61                   |                                       |                |                                  |            |             |                    | F          |
| GPR62                   |                                       |                |                                  |            |             |                    |            |
| GPR63                   |                                       |                |                                  |            |             |                    |            |
| GPR65                   |                                       |                |                                  |            |             |                    |            |
| GPR68                   |                                       | ?              | √                                |            |             |                    |            |
| GPR75                   |                                       |                |                                  |            |             |                    |            |
| GPR78                   |                                       |                |                                  |            |             |                    |            |
| GPR79                   |                                       |                |                                  |            |             |                    |            |
| GPR82                   |                                       |                |                                  |            |             |                    |            |
| GPR83                   |                                       |                |                                  |            |             |                    |            |

| Target<br>(IUPHAR name) | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|-------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                         | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| GPR84                   |                                       |                |                                  |            |             |                    |            |
| GPR85                   |                                       |                |                                  |            |             |                    |            |
| GPR87                   |                                       |                |                                  |            |             |                    |            |
| GPR88                   |                                       |                | √                                | ?          |             |                    |            |
| GPR101                  |                                       |                |                                  |            |             |                    |            |
| GPR119                  |                                       |                |                                  |            |             |                    |            |
| GPR120                  |                                       |                |                                  |            |             |                    |            |
| GPR132                  |                                       |                |                                  |            |             |                    |            |
| GPR135                  |                                       |                |                                  |            |             |                    |            |
| GPR139                  |                                       |                |                                  |            |             |                    |            |
| GPR141                  |                                       |                |                                  |            |             |                    |            |
| GPR142                  |                                       |                |                                  |            |             |                    |            |
| GPR146                  |                                       |                |                                  |            |             |                    |            |
| GPR148                  |                                       |                |                                  |            |             |                    |            |
| GPR150                  |                                       |                |                                  |            |             |                    |            |
| GPR151                  |                                       |                |                                  |            |             |                    |            |
| GPR152                  |                                       |                |                                  |            |             |                    |            |
| GPR153                  |                                       |                |                                  |            |             |                    |            |
| GPR160                  |                                       |                |                                  |            |             |                    |            |
| GPR161                  |                                       |                |                                  |            |             |                    |            |
| GPR162                  |                                       |                |                                  |            |             |                    |            |
| GPR171                  |                                       |                |                                  |            |             |                    |            |
| GPR173                  |                                       |                |                                  |            |             |                    |            |
| GPR174                  |                                       |                |                                  |            |             |                    |            |
| GPR176                  |                                       |                |                                  |            |             |                    |            |
| GPR182                  |                                       |                |                                  |            |             |                    |            |
| GPR183                  |                                       |                |                                  |            |             |                    |            |
| LGR4                    |                                       |                |                                  |            |             |                    |            |
| LGR5                    |                                       |                |                                  |            |             |                    |            |
| LGR6                    |                                       |                |                                  |            |             |                    |            |
| LPAR6                   |                                       |                |                                  |            |             |                    |            |
| MAS1                    |                                       |                |                                  |            |             |                    |            |
| MAS1L                   |                                       |                |                                  |            |             |                    |            |
| MRGPRD                  |                                       |                |                                  |            |             |                    |            |
| MRGPRE                  |                                       |                |                                  |            |             |                    |            |
| MRGPRF                  |                                       |                |                                  |            |             |                    |            |
| MRGPRG                  |                                       |                |                                  |            |             |                    |            |
| MRGPRX1                 |                                       |                |                                  |            |             |                    |            |
| MRGPRX2                 |                                       |                |                                  |            |             |                    |            |
| MRGPRX3                 |                                       |                |                                  | √          |             |                    |            |
| MRGPRX4                 |                                       |                |                                  | √          |             |                    |            |
| OPN3                    |                                       |                |                                  |            |             |                    |            |
| OPN5                    |                                       |                |                                  |            |             |                    |            |
| OXGR1                   |                                       |                |                                  |            |             |                    |            |

| Target<br>(IUPHAR name) | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|-------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                         | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| P2RY8                   |                                       |                |                                  |            |             |                    |            |
| P2RY10                  |                                       |                |                                  |            |             |                    |            |
| SUCNR1                  |                                       |                |                                  |            |             |                    |            |
| TAAR2                   |                                       |                |                                  |            |             |                    |            |
| TAAR3                   |                                       |                |                                  |            |             |                    | F          |
| TAAR4                   |                                       |                |                                  |            |             |                    | F          |
| TAAR5                   |                                       |                |                                  |            |             |                    | F          |
| TAAR6                   |                                       |                |                                  |            |             |                    |            |
| TAAR8                   |                                       |                |                                  |            |             |                    |            |
| TAAR9                   |                                       |                |                                  |            |             |                    |            |
| Class B orphan GPCRs    |                                       |                |                                  |            |             |                    |            |
| BAI1                    |                                       |                |                                  |            |             |                    |            |
| BAI2                    |                                       |                |                                  |            |             |                    |            |
| BAI3                    |                                       |                |                                  |            |             |                    |            |
| CD97                    |                                       |                |                                  |            |             |                    |            |
| CELSR1                  |                                       |                |                                  |            |             |                    |            |
| CELSR2                  |                                       |                |                                  |            |             |                    |            |
| CELSR3                  |                                       |                |                                  |            |             |                    |            |
| ELTD1                   |                                       |                |                                  |            |             |                    |            |
| EMR1                    |                                       |                |                                  |            |             |                    | F          |
| EMR2                    |                                       |                |                                  |            |             |                    | F          |
| EMR3                    |                                       |                |                                  |            |             |                    | F          |
| EMR4P                   |                                       |                |                                  |            |             |                    |            |
| GPR56                   |                                       |                |                                  |            |             |                    |            |
| GPR64                   |                                       |                |                                  |            |             |                    |            |
| GPR97                   |                                       |                |                                  |            |             |                    |            |
| GPR98                   |                                       |                |                                  |            |             | MASS1              | F          |
| GPR110                  |                                       |                |                                  |            |             |                    |            |
| GPR112                  |                                       |                |                                  |            |             |                    |            |
| GPR113                  |                                       |                |                                  |            |             |                    |            |
| GPR114                  |                                       |                |                                  |            |             |                    |            |
| GPR115                  |                                       |                |                                  |            |             |                    | F          |
| GPR116                  |                                       |                |                                  |            |             |                    |            |
| GPR123                  |                                       |                |                                  |            |             |                    | F          |
| GPR124                  |                                       |                |                                  |            |             |                    |            |
| GPR125                  |                                       |                |                                  |            |             |                    | F          |
| GRP126                  |                                       |                |                                  |            |             |                    |            |
| GPR127                  |                                       |                |                                  |            |             |                    |            |
| GPR128                  |                                       |                |                                  |            |             |                    | F          |
| GPR133                  |                                       |                |                                  |            |             |                    | F          |
| GPR144                  |                                       |                |                                  |            |             |                    |            |
| GPR157                  |                                       |                |                                  |            |             |                    | F          |
| LPHN1                   |                                       |                |                                  |            |             |                    |            |
| LPHN2                   |                                       |                |                                  |            |             |                    |            |

| Target<br>(IUPHAR name) | Available assays at PDSP (new system) |                |                                  |            |             | Synonym<br>Or Note | PDSP<br>V1 |
|-------------------------|---------------------------------------|----------------|----------------------------------|------------|-------------|--------------------|------------|
|                         | Binding                               | G <sub>q</sub> | G <sub>i</sub> or G <sub>s</sub> | β-arrestin | Other Assay |                    |            |
| LPHN3                   |                                       |                |                                  |            |             |                    |            |
| Class C orphan GPCRs    |                                       |                |                                  |            |             |                    |            |
| GPR156                  |                                       |                |                                  |            |             |                    |            |
| GPR158                  |                                       |                |                                  |            |             |                    |            |
| GPR179                  |                                       |                |                                  |            |             |                    |            |
| GPRC5A                  |                                       |                |                                  |            |             |                    |            |
| GPRC5B                  |                                       |                |                                  |            |             |                    |            |
| GPRC5C                  |                                       |                |                                  |            |             |                    |            |
| GPRC5D                  |                                       |                |                                  |            |             |                    |            |
| Other 7TM receptors     |                                       |                |                                  |            |             |                    |            |
| GPR107                  |                                       |                |                                  |            |             |                    |            |
| GPR137                  |                                       |                |                                  |            |             | TM7SF1             | F          |
| GPR143                  |                                       |                |                                  |            |             |                    |            |

## Section 4: List of publications supported by NIMH-PDSP

Complete list of research publications co-authored or otherwise supported by the NIMH PDSP. References are listed in the order to publication data, from the most recent to the oldest. Detailed validated protocols and modified protocols for all assays can also be found in these peer-reviewed publications.

1. Wang C, Jiang Y, Ma J, Wu H, Wacker D, Katritch V, Han GW, Liu W, Huang XP, Vardy E, McCorvy JD, Gao X, Zhou E, Melcher K, Zhang C, Bai F, Yang H, Yang L, Jiang H, **Roth BL**, Cherezov V, Stevens RC and Xu H Structural Basis for Molecular Recognition at Serotonin Receptors Science 2013. NIH Manuscript ID (NIHMSID) 350378.
2. Wacker D, Wang C, Katritch V, Han GW, Huang XP, Vardy E, McCorvy JD, Jiang Y, Chu M, Siu F, Liu W, Xu HE, Cherezov V, **Roth BL** and Stevens RC Structural Features for Functional Selectivity at Serotonin Receptors Science 2013. NIMH Manuscript ID (NIHMSID) 450382 (\*BRL and RCS = Co-corresponding authors)
3. Vardy E, **Roth BL**. Conformational Ensembles in GPCR Activation. Cell 152(3): 385-6, 2013
4. Lin H, Sassano MF, **Roth BL**, Shoichet BK. A pharmacological organization of G protein-coupled receptors. Nat Methods 10(2): 140-6. (\*BLR and BKS = Co-corresponding authors)
5. Iversen L, Gibbons S, Treble R, Setola V, Huang XP, **Roth BL**. Neurochemical profiles of some novel psychoactive substances. Eur J Pharmacol 700(1-3): 147-51, 2013
6. Besnard J, Ruda GF, Setola V, Abecassis K, Rodriguiz RM, Huang XP, Norval S, Sassano MF, Shin AI, Webster LA, Simeons FR, Stojanovski L, Prat A, Seidah NG, Constam DB, Bickerton GR, Read KD, Wetsel WC, Gilbert IH, **Roth BL**, Hopkins AL. Automated design of ligands to polypharmacological profiles. Nature 492(7428): 215-20, 2012. (\*BRL and ALH = Co-corresponding authors)
7. **Roth BL**, Marshall FH. NOBEL 2012 Chemistry: Studies of a ubiquitous receptor family. Nature 492(7427): 57, 2012.
8. Namjoshi OA, Wang ZJ, Rallapalli SK, Johnson EM, Jr., Johnson YT, Ng H, Ramerstorfer J, Varagic Z, Sieghart W, Majumder S, **Roth BL**, Rowlett JK, Cook JM. Search for

- alpha3beta(2)/(3)gamma2 subtype selective ligands that are stable on human liver microsomes. *Bioorg Med Chem* 21(1): 93-101, 2012
9. Pigott A, Frescas S, McCorvy JD, Huang XP, **Roth BL**, Nichols DE. trans-2-(2,5-Dimethoxy-4-iodophenyl)cyclopropylamine and trans-2-(2,5-dimethoxy-4-bromophenyl)cyclopropylamine as potent agonists for the 5-HT(2) receptor family. *Beilstein J Org Chem* 8(1705-9), 2012.
  10. White KL, **Roth BL**. Psychotomimetic effects of kappa opioid receptor agonists. *Biol Psychiatry* 72(10): 797-8, 2012.
  11. Houle S, Mamo DC, Mulsant BH, Pollock BG, van Oosten EM, Vasdev N, Wilson AA. Towards the development of new subtype-specific muscarinic receptor radiopharmaceuticals-- radiosynthesis and ex vivo biodistribution of [<sup>18</sup>F]3-(4-(2-(2-(2-fluoroethoxy)ethoxy)ethylthio)-1,2,5-thiadiazol-3-yl)-1-methyl-1, 2,5,6-tetrahydropyridine. *Canadian Journal of Chemistry* 88:1222, 2012.
  12. Hulubei V, Meikrantz SB, Quincy DA, Houle T, McKenna JI, Rogers ME, Steiger S, Natale NR. 4-Isoxazolyl-1,4-dihydropyridines exhibit binding at the multidrug-resistance transporter. *Bioorg Med Chem* 20(22):6613-6620, 2012.
  13. Farrell MS, **Roth BL**. Pharmacosynthetics: Reimagining the pharmacogenetic approach. *Brain Res* 2012, doi: 10.1016/j.brainres.2012.09.043. [Epub ahead of print]
  14. Sidique S, Dhanya RP, Sheffler DJ, Nickols HH, Yang L, Dahl R, Mangravita-Novo A, Smith LH, D'Souza MS, Semenova S, Conn PJ, Markou A, Cosford ND. Orally Active Metabotropic Glutamate Subtype 2 Receptor Positive Allosteric Modulators: Structure-Activity Relationships and Assessment in a Rat Model of Nicotine Dependence. *J Med Chem* 55(22): 9434-9445, 2012.
  15. Liu S, Patel RY, Daga PR, Liu H, Fu G, Doerksen RJ, Chen Y, Wilkins DE. Combined rule extraction and feature elimination in supervised classification. *IEEE Trans Nanobioscience* 11(3): 228-36, 2012.
  16. Banister SD, Manoli M, Doddareddy MR, Hibbs DE, Kassiou M. A sigma(1) receptor pharmacophore derived from a series of N-substituted 4-azahexacyclo[5.4.1.0(2,6).0(3,10).0(5,9).0(8,11)]dodecan-3-ols (AHDs). *Bioorg Med Chem Lett* 22(19): 6053-8, 2012.

17. Zhu H, Urban DJ, Blashka J, McPheeters MT, Kroeze WK, Mieczkowski P, Overholser JC, Jurjus GJ, Dieter L, Mahajan GJ, Rajkowska G, Wang Z, Sullivan PF, Stockmeier CA, **Roth BL**. Quantitative analysis of focused a-to-I RNA editing sites by ultra-high-throughput sequencing in psychiatric disorders. *PLoS One* 7(8): e43227, 2012.
18. Little JW, Chen Z, Doyle T, Porreca F, Ghaffari M, Bryant L, Neumann WL, Salvemini D. Supraspinal peroxynitrite modulates pain signaling by suppressing the endogenous opioid pathway. *J Neurosci* 32(32): 10797-808, 2012.
19. Twiner MJ, Doucette GJ, Rasky A, Huang XP, **Roth BL**, Sanguinetti MC. Marine algal toxin azaspiracid is an open-state blocker of hERG potassium channels. *Chem Res Toxicol* 25(9): 1975-84, 2012.
20. Twiner MJ, Doucette GJ, Rasky A, Huang XP, **Roth BL**, Sanguinetti MC. Marine algal toxin azaspiracid is an open-state blocker of HERG potassium channels. *Chem Res Toxicol* 25(9): 1975-84, 2012.
21. Chen X, Sassano MF, Zheng L, Setola V, Chen M, Bai X, Frye SV, Wetsel WC, **Roth BL**, Jin J. Structure--functional selectivity relationship studies of beta-arrestin-biased dopamine D(2) receptor agonists. *J Med Chem* 55(16): 7141-53, 2012.
22. Schulz SB, Heidmann KE, Mike A, Klafit ZJ, Heinemann U, Gerevich Z. First and second generation antipsychotics influence hippocampal gamma oscillations by interactions with 5-HT(3) and D(3) receptors. *Br J Pharmacol* 2012.
23. Chen G, Cho SJ, Huang XP, Jensen NH, Svennebring A, Sassano MF, **Roth BL**, Kozikowski AP. Rational Drug Design Leading to the Identification of a Potent 5-HT(2C) Agonist Lacking 5-HT(2B) Activity. *ACS Med Chem Lett* 2(12): 929-932, 2012.
24. Mangin F, Dilly S, Joly B, Scuvee-Moreau J, Evans J, Setola V, **Roth BL**, Liegeois JF. Moderate chemical modifications of WAY-100635 improve the selectivity for 5-HT1A versus D4 receptors. *Bioorg Med Chem Lett* 22(14): 4550-4, 2012.
25. Frankowski KJ, Hedrick MP, Gosalia P, Li K, Shi S, Whipple D, Ghosh P, Prisinzano TE, Schoenen FJ, Su Y, Vasile S, Sergienko E, Gray W, Hariharan S, Milan L, Heynen-Genel S, Mangravita-Novo A, Vicchiarelli M, Smith LH, Streicher JM, Caron MG, Barak LS, Bohn LM, Chung TD, Aube J. Discovery of Small Molecule Kappa Opioid Receptor Agonist and

- Antagonist Chemotypes through a HTS and Hit Refinement Strategy. *ACS Chem Neurosci* 3(3): 221-236, 2012.
26. Gregori-Puigjane E, Setola V, Hert J, Crews BA, Irwin JJ, Lounkine E, Marnett L, **Roth BL**, Shoichet BK. Identifying mechanism-of-action targets for drugs and probes. *Proc Natl Acad Sci U S A* 109(28): 11178-83, 2012.
27. Almeida C, Hemberger Y, Schmitt SM, Bouhired S, Natesan L, Kehraus S, Dimas K, Gutschow M, Bringmann G, König GM. Marilines A-C: novel phthalimidines from the sponge-derived fungus *Stachylidium* sp. *Chemistry* 18(28): 8827-34, 2012.
28. Lin X, Huang XP, Chen G, Whaley R, Peng S, Wang Y, Zhang G, Wang SX, Wang S, **Roth BL**, Huang N. Life beyond kinases: structure-based discovery of sorafenib as nanomolar antagonist of 5-HT receptors. *J Med Chem* 55(12): 5749-59, 2012.
29. Williams AA, Ingram WM, Levine S, Resnik J, Kamel CM, Lish JR, Elizalde DI, Janowski SA, Shoker J, Kozlenkov A, Gonzalez-Maeso J, Gallitano AL. Reduced levels of serotonin 2A receptors underlie resistance of *Egr3*-deficient mice to locomotor suppression by clozapine. *Neuropsychopharmacology* 37(10): 2285-98, 2012.
30. Banister SD, Rendina LM, Kassiou M. 7-Azabicyclo[2.2.1]heptane as a scaffold for the development of selective sigma-2 (sigma2) receptor ligands. *Bioorg Med Chem Lett* 22(12): 4059-63, 2012.
31. Thompson AA, Liu W, Chun E, Katritch V, Wu H, Vardy E, Huang XP, Trapella C, Guerrini R, Calo G, **Roth BL**, Cherezov V, Stevens RC. Structure of the nociceptin/orphanin FQ receptor in complex with a peptide mimetic. *Nature* 485(7398): 395-9, 2012.
32. Hajjo R, Setola V, **Roth BL**, Tropsha A. Chemocentric informatics approach to drug discovery: identification and experimental validation of selective estrogen receptor modulators as ligands of 5-hydroxytryptamine-6 receptors and as potential cognition enhancers. *J Med Chem* 55(12): 5704-19, 2012.
33. Zhu XY, Etukala JR, Eyunni SV, Setola V, **Roth BL**, Ablordeppey SY. Benzothiazoles as probes for the 5HT1A receptor and the serotonin transporter (SERT): a search for new dual-acting agents as potential antidepressants. *Eur J Med Chem* 53(124-32), 2012.

34. Baitz HA, Thornton AE, Procyshyn RM, Smith GN, MacEwan GW, Kopala LC, Barr AM, Lang DJ, Honer WG. Antipsychotic medications: linking receptor antagonism to neuropsychological functioning in first episode psychosis. *J Int Neuropsychol Soc* 18(4): 717-27, 2012.
35. Wu H, Li J, Fronczek FR, Ferreira D, Burandt CL, Jr., Setola V, **Roth BL**, Zjawiony JK. Labdane Diterpenoids from *Leonotis leonurus*. *Phytochemistry* 2012.
36. Garner AR, Rowland DC, Hwang SY, Baumgaertel K, **Roth BL**, Kentros C, Mayford M. Generation of a synthetic memory trace. *Science* 335(6075): 1513-6, 2012.
37. Wu H, Wacker D, Mileni M, Katritch V, Han GW, Vardy E, Liu W, Thompson AA, Huang XP, Carroll FI, Mascarella SW, Westkaemper RB, Mosier PD, **Roth BL**, Cherezov V, Stevens RC. Structure of the human kappa-opioid receptor in complex with JDTic. *Nature* 485(7398): 327-32, 2012.
38. Williams AJ, Ekins S, Tkachenko V. Towards a gold standard: regarding quality in public domain chemistry databases and approaches to improving the situation. *Drug Discov Today* 17(13-14): 685-701, 2012.
39. Smith AL, Freeman SM, Stehouwer JS, Inoue K, Voll RJ, Young LJ, Goodman MM. Synthesis and evaluation of C-11, F-18 and I-125 small molecule radioligands for detecting oxytocin receptors. *Bioorg Med Chem* 20(8): 2721-38, 2012.
40. Norman TR. The effect of agomelatine on 5HT(2C) receptors in humans: a clinically relevant mechanism? *Psychopharmacology (Berl)* 221(1): 177-8; author reply 179, 2012.
41. Peprah K, Zhu XY, Eyunni SV, Etukala JR, Setola V, **Roth BL**, Ablordeppey SY. Structure-activity relationship studies of SYA 013, a homopiperazine analog of haloperidol. *Bioorg Med Chem* 20(5): 1671-8, 2012.
42. Bao X, Lu S, Liow JS, Zoghbi SS, Jenko KJ, Clark DT, Gladding RL, Innis RB, Pike VW. Radiosynthesis and evaluation of an (18)F-labeled positron emission tomography (PET) radioligand for brain histamine subtype-3 receptors based on a nonimidazole 2-aminoethylbenzofuran chemotype. *J Med Chem* 55(5): 2406-15, 2012.
43. Lee HM, **Roth BL**. Hallucinogen actions on human brain revealed. *Proc Natl Acad Sci U S A* 109(6): 1820-1, 2012.

44. Mavel S, Mincheva Z, Meheux N, Carcenac Y, Guilloteau D, Abarbri M, Emond P. QSAR study and synthesis of new phenyltropanes as ligands of the dopamine transporter (DAT). *Bioorg Med Chem* 20(4): 1388-95, 2012.
45. Duncan GE, Koller BH, Moy SS. Effects of the selective kainate receptor antagonist ACET on altered sensorimotor gating in a genetic model of reduced NMDA receptor function. *Brain Res* 1443(98-105), 2012.
46. Toll L, Zaveri NT, Polgar WE, Jiang F, Khroyan TV, Zhou W, Xie XS, Stauber GB, Costello MR, Leslie FM. AT-1001: a high affinity and selective alpha3beta4 nicotinic acetylcholine receptor antagonist blocks nicotine self-administration in rats. *Neuropsychopharmacology* 37(6): 1367-76, 2012.
47. Peprah K, Zhu XY, Eyunni SV, Setola V, **Roth BL**, Ablordeppey SY. Multi-receptor drug design: Haloperidol as a scaffold for the design and synthesis of atypical antipsychotic agents. *Bioorg Med Chem* 20(3): 1291-7, 2012.
48. Almeida C, Kehraus S, Prudencio M, Konig GM. Marilonones A-C, phthalides from the sponge-derived fungus *Stachylidium* sp. *Beilstein J Org Chem* 7(1636-42), 2012.
49. Halberstadt AL, Nichols DE, Geyer MA. Behavioral effects of alpha, alpha, beta, beta-tetra deuterio-5-MeO-DMT in rats: comparison with 5-MeO-DMT administered in combination with a monoamine oxidase inhibitor. *Psychopharmacology (Berl)* 221(4): 709-18, 2012.
50. Kim MK, Lee HS, Kim S, Cho SY, **Roth BL**, Chong Y, Choo H. 4-Aminoethylpiperazinyl aryl ketones with 5-HT(1)A/5-HT(7) selectivity. *Bioorg Med Chem* 20(2): 1139-48, 2012.
51. Carlsson J, Coleman RG, Setola V, Irwin JJ, Fan H, Schlessinger A, Sali A, Roth BL, and Shoichet BK. Ligand discovery from a dopamine D3 receptor homology model and crystal structure. *Nature Chemical Biology* 7: 769-778, 2011. (BRL and BKS = Co-corresponding authors).
52. Huang HS, Allen JA, Mabb AM, King IF, Miriyala J, Taylor-Blake B, Sciaky N, Dutton JW, Jr., Lee HM, Chen X, Jin J, Bridges AS, Zylka MJ, **Roth BL**, Philpot BD. Topoisomerase inhibitors unsilence the dormant allele of *Ube3a* in neurons. *Nature* 481(7380): 185-9, 2011. (BLR, MJZ, and BDP = Co-corresponding authors).

53. Yu LF, Tuckmantel W, Eaton JB, Caldarone B, Fedolak A, Hanania T, Brunner D, Lukas RJ, Kozikowski AP. Identification of novel alpha4beta2-nicotinic acetylcholine receptor (nAChR) agonists based on an isoxazole ether scaffold that demonstrate antidepressant-like activity. *J Med Chem* 55(2): 812-23, 2011.
54. You H, Gadotti VM, Petrov RR, Zamponi GW, Diaz P. Functional characterization and analgesic effects of mixed cannabinoid receptor/T-type channel ligands. *Mol Pain* 7(89), 2011.
55. Hedrick MP, Gosalia P, Li K, Frankowski KJ, Shi S, Prisinzano TE, Schoenen F, Aube J, Su Y, Stonich D, Vasile S, Sergienko E, Gray W, Hariharan S, Milan L, Heynen-Genel S, Vicchiarelli M, Mangravita-Novo A, Streicher JM, Smith LH, Chung TDY, Caron M, Bohn LM, Barak LS. Antagonist for the Kappa Opioid Receptor. 2011.
56. Lapish CC, Belardetti F, Ashby DM, Ahn S, Butts KA, So K, Macrae CM, Hynd JJ, Miller JJ, Phillips AG. A preclinical assessment of d.l-govadine as a potential antipsychotic and cognitive enhancer. *Int J Neuropsychopharmacol* 15(10): 1441-55, 2011.
57. Marder SR, Roth B, Sullivan PF, Scolnick EM, Nestler EJ, Geyer MA, Welnerberger DR, Karayiorgou M, Guidotti A, Gingrich J, Akbarian S, Buchanan RW, Lieberman JA, Conn PJ, Haggarty SJ, Law AJ, Campbell B, Krystal JH, Moghaddam B, Sawa A, Caron MG, George SR, Allen JA, Solis M. Advancing drug discovery for schizophrenia. *Ann N Y Acad Sci* 1236:30-43, 2011.
58. Allen JA, Yost JM, Setola V, Chen X, Sassano MF, Chen M, Peterson S, Yadav PN, Huang XP, Feng B, Jensen NH, Che X, Bai X, Frye SV, Wetsel WC, Caron MG, Javitch JA, **Roth BL**, Jin J. Discovery of beta-arrestin-biased dopamine D2 ligands for probing signal transduction pathways essential for antipsychotic efficacy. *Proc Natl Acad Sci U S A* 108(45): 18488-93, 2011.
59. Al Hussainy R, Verbeek J, van der Born D, Booij J, Herscheid JK. Design, synthesis and in vitro evaluation of bridgehead fluoromethyl analogs of N-{2-[4-(2-methoxyphenyl)piperazin-1-yl]ethyl}-N-(pyridin-2-yl)cyclohexanecarboxamide (WAY-100635) for the 5-HT(1A) receptor. *Eur J Med Chem* 46(12): 5728-35, 2011.

60. Dutta AK, Gopishetty B, Gogoi S, Ali S, Zhen J, Reith M. The novel trisubstituted pyran derivative D-142 has triple monoamine reuptake inhibitory activity and exerts potent antidepressant-like activity in rodents. *Eur J Pharmacol* 671(1-3): 39-44, 2011.
61. Quednow BB, Komater M, Geyer MA, Vollenweider FX. Psilocybin-induced deficits in automatic and controlled inhibition are attenuated by ketanserin in healthy human volunteers. *Neuropsychopharmacology* 37(3): 630-40, 2011.
62. Holtman JR, Dwoskin LP, Dowell C, Wala EP, Zhang Z, Crooks PA, McIntosh JM. The novel small molecule  $\alpha 9 \alpha 10$  nicotinic acetylcholine receptor antagonist ZZ-204G is analgesic. *Eur J Pharmacol* 670(2-3): 500-8, 2011.
63. Liu J, Yu LF, Eaton JB, Caldarone B, Cavino K, Ruiz C, Terry M, Fedolak A, Wang D, Ghavami A, Lowe DA, Brunner D, Lukas RJ, Kozikowski AP. Discovery of isoxazole analogues of sazetidine-A as selective  $\alpha 4 \beta 2$ -nicotinic acetylcholine receptor partial agonists for the treatment of depression. *J Med Chem* 54(20): 7280-8, 2011.
64. Lebar MD, Hahn KN, Mutka T, Maignan P, McClintock JB, Amsler CD, van Olphen A, Kyle DE, Baker BJ. CNS and antimalarial activity of synthetic meridianin and psammopemmin analogs. *Bioorg Med Chem* 19(19): 5756-62, 2011.
65. Moussa IA, Banister SD, Akladios FN, Chua SW, Kassiou M. Effects of linker elongation in a series of N-(2-benzofuranylmethyl)-N'-(methoxyphenylalkyl)piperazine sigma(1) receptor ligands. *Bioorg Med Chem Lett* 21(19): 5707-10, 2011.
66. Petrov RR, Ferrini ME, Jaffar Z, Thompson CM, Roberts K, Diaz P. Design and evaluation of a novel fluorescent CB2 ligand as probe for receptor visualization in immune cells. *Bioorg Med Chem Lett* 21(19): 5859-62, 2011.
67. Ettrup A, Mikkelsen JD, Lehel S, Madsen J, Nielsen EO, Palner M, Timmermann DB, Peters D, Knudsen GM. 11C-NS14492 as a novel PET radioligand for imaging cerebral  $\alpha 7$  nicotinic acetylcholine receptors: in vivo evaluation and drug occupancy measurements. *J Nucl Med* 52(9): 1449-56, 2011.
68. Chu UB, Hajipour AR, Ramachandran S, Ruoho AE. Characterization of interactions of 4-nitrophenylpropyl-N-alkylamine with  $\sigma$  receptors. *Biochemistry* 50(35): 7568-78, 2011.

69. Banister SD, Yoo DT, Chua SW, Cui J, Mach RH, Kassiou M. N-Arylalkyl-2-azaadamantanes as cage-expanded polycarbocyclic sigma (sigma) receptor ligands. *Bioorg Med Chem Lett* 21(18): 5289-92, 2011.
70. Stephenson KA, Wilson AA, Houle S, Vasdev N. Synthesis and in vitro evaluation of derivatives of the beta(1)-adrenergic receptor antagonist HX-CH 44. *Bioorg Med Chem Lett* 21(18): 5506-9, 2011.
71. Vedadi M, Barsyte-Lovejoy D, Liu F, Rival-Gervier S, Allali-Hassani A, Labrie V, Wigle TJ, Dimaggio PA, Wasney GA, Siarheyeva A, Dong A, Tempel W, Wang SC, Chen X, Chau I, Mangano TJ, Huang XP, Simpson CD, Pattenden SG, Norris JL, Kireev DB, Tripathy A, Edwards A, **Roth BL**, Janzen WP, Garcia BA, Petronis A, Ellis J, Brown PJ, Frye SV, Arrowsmith CH, Jin J. A chemical probe selectively inhibits G9a and GLP methyltransferase activity in cells. *Nat Chem Biol* 7(8): 566-74, 2011.
72. Yadav PN, Kroeze WK, Farrell MS, **Roth BL**. Antagonist functional selectivity: 5-HT<sub>2A</sub> serotonin receptor antagonists differentially regulate 5-HT<sub>2A</sub> receptor protein level in vivo. *J Pharmacol Exp Ther* 339(1): 99-105, 2011.
73. Loyd DR, Weiss G, Henry MA, Hargreaves KM. Serotonin increases the functional activity of capsaicin-sensitive rat trigeminal nociceptors via peripheral serotonin receptors. *Pain* 152(10): 2267-76, 2011.
74. Sromek AW, Si YG, Zhang T, George SR, Seeman P, Neumeyer JL. Synthesis and Biological Evaluation of N-Fluoroalkyl and 2-Fluoroalkoxy Substituted Aporphines: Potential PET Ligands for Dopamine D(2) Receptors. *ACS Med Chem Lett* 2(3): 189-194, 2011.
75. Tarselli MA, Raehal KM, Brasher AK, Streicher JM, Groer CE, Cameron MD, Bohn LM, Micalizio GC. Synthesis of conolidine, a potent non-opioid analgesic for tonic and persistent pain. *Nat Chem* 3(6): 449-53, 2011.
76. Banister SD, Moussa IA, Jorgensen WT, Chua SW, Kassiou M. Molecular hybridization of 4-azahexacyclo[5.4.1.0(2,6).0(3,10).0(5,9).0(8,11)]dodecane-3-ol with sigma (sigma) receptor ligands modulates off-target activity and subtype selectivity. *Bioorg Med Chem Lett* 21(12): 3622-6, 2011.

77. Al Hussainy R, Verbeek J, van der Born D, Braker AH, Leysen JE, Knol RJ, Booij J, Herscheid JK. Design, synthesis, radiolabeling, and in vitro and in vivo evaluation of bridgehead iodinated analogues of N-{2-[4-(2-methoxyphenyl)piperazin-1-yl]ethyl}-N-(pyridin-2-yl)cyclohexanecarboxamide (WAY-100635) as potential SPECT ligands for the 5-HT<sub>1A</sub> receptor. *J Med Chem* 54(10): 3480-91, 2011.
78. Ates-Alagoz Z, Sun S, Wallach J, Adejare A. Syntheses and pharmacological evaluations of novel N-substituted bicyclo-heptane-2-amines at N-methyl-D-aspartate receptors. *Chem Biol Drug Des* 78(1): 25-32, 2011.
79. Schurer SC, Vempati U, Smith R, Southern M, Lemmon V. BioAssay ontology annotations facilitate cross-analysis of diverse high-throughput screening data sets. *J Biomol Screen* 16(4): 415-26, 2011.
80. Gopishetty B, Hazeldine S, Santra S, Johnson M, Modi G, Ali S, Zhen J, Reith M, Dutta A. Further structure-activity relationship studies on 4-(((3S,6S)-6-benzhydryltetrahydro-2H-pyran-3-yl)amino)methylphenol: identification of compounds with triple uptake inhibitory activity as potential antidepressant agents. *J Med Chem* 54(8): 2924-32, 2011.
81. Donohue SR, Dannals RF, Halldin C, Pike VW. N-(4-cyanotetrahydro-2H-pyran-4-yl) and N-(1-cyanocyclohexyl) derivatives of 1,5-diarylpyrazole-3-carboxamides showing high affinity for 18 kDa translocator protein and/or cannabinoid receptors. *J Med Chem* 54(8): 2961-70, 2011.
82. Rogan SC, **Roth BL**. Remote control of neuronal signaling. *Pharmacol Rev* 63(2): 291-315, 2011.
83. Yow TT, Pera E, Absalom N, Heblinski M, Johnston GA, Hanrahan JR, Chebib M. Naringin directly activates inwardly rectifying potassium channels at an overlapping binding site to tertiapin-Q. *Br J Pharmacol* 163(5): 1017-33, 2011.
84. Krashes MJ, Koda S, Ye C, Rogan SC, Adams AC, Cusher DS, Maratos-Flier E, **Roth BL**, Lowell BB. Rapid, reversible activation of AgRP neurons drives feeding behavior in mice. *J Clin Invest* 121(4): 1424-8, 2011.
85. Zou MF, Cao J, Rodriguez AL, Conn PJ, Newman AH. Design and synthesis of substituted N-(1,3-diphenyl-1H-pyrazol-5-yl)benzamides as positive allosteric modulators of the metabotropic glutamate receptor subtype 5. *Bioorg Med Chem Lett* 21(9): 2650-4, 2011.

86. **Roth BL.** Irving Page Lecture: 5-HT(2A) serotonin receptor biology: interacting proteins, kinases and paradoxical regulation. *Neuropharmacology* 61(3): 348-54, 2011.
87. McKenna DJ, Ruiz JM, Hoyer TR, **Roth BL**, Shoemaker AT. Receptor screening technologies in the evaluation of Amazonian ethnomedicines with potential applications to cognitive deficits. *J Ethnopharmacol* 134(2): 475-92, 2011.
88. Allen JA, Yadav PN, Setola V, Farrell M, **Roth BL.** Schizophrenia risk gene CAV1 is both pro-psychotic and required for atypical antipsychotic drug actions in vivo. *Transl Psychiatry* 1(e33), 2011.
89. Ettrup A, Hansen M, Santini MA, Paine J, Gillings N, Palner M, Lehel S, Herth MM, Madsen J, Kristensen J, Begtrup M, Knudsen GM. Radiosynthesis and in vivo evaluation of a series of substituted 11C-phenethylamines as 5-HT (2A) agonist PET tracers. *Eur J Nucl Med Mol Imaging* 38(4): 681-93, 2010.
90. Huang XP, Mangano T, Hufeisen S, Setola V, **Roth BL.** Identification of human Ether-a-go-go related gene modulators by three screening platforms in an academic drug-discovery setting. *Assay Drug Dev Technol* 8(6): 727-42, 2010.
91. Halberstadt AL, Koedood L, Powell SB, Geyer MA. Differential contributions of serotonin receptors to the behavioral effects of indoleamine hallucinogens in mice. *J Psychopharmacol* 25(11): 1548-61, 2010.
92. Banister SD, Moussa IA, Beinat C, Reynolds AJ, Schiavini P, Jorgensen WT, Kassiou M. Trishomocubane as a scaffold for the development of selective dopamine transporter (DAT) ligands. *Bioorg Med Chem Lett* 21(1): 38-41, 2010.
93. Brudeli B, Moltzau LR, Andressen KW, Krobert KA, Klaveness J, Levy FO. Synthesis and pharmacological properties of novel hydrophilic 5-HT<sub>4</sub> receptor antagonists. *Bioorg Med Chem* 18(24): 8600-13, 2011.
94. Keiser MJ, Irwin JJ, Shoichet BK. The chemical basis of pharmacology. *Biochemistry* 49(48): 10267-76, 2010.
95. Yadav PN, Abbas AI, Farrell MS, Setola V, Sciaky N, Huang XP, Kroeze WK, Crawford LK, Piel DA, Keiser MJ, Irwin JJ, Shoichet BK, Deneris ES, Gingrich J, Beck SG, **Roth BL.** The presynaptic

component of the serotonergic system is required for clozapine's efficacy.

Neuropsychopharmacology 36(3): 638-51, 2010.

96. Lin Z, Antemano RR, Hughen RW, Tianero MD, Peraud O, Haygood MG, Concepcion GP, Olivera BM, Light A, Schmidt EW. Pulicatins A-E, neuroactive thiazoline metabolites from cone snail-associated bacteria. *J Nat Prod* 73(11): 1922-6, 2010.
97. Hajjo R, Grulke CM, Golbraikh A, Setola V, Huang XP, **Roth BL**, Tropsha A. Development, validation, and use of quantitative structure-activity relationship models of 5-hydroxytryptamine (2B) receptor ligands to identify novel receptor binders and putative valvulopathic compounds among common drugs. *J Med Chem* 53(21): 7573-86, 2010.
98. Allen JA, **Roth BL**. Strategies to discover unexpected targets for drugs active at G protein-coupled receptors. *Annu Rev Pharmacol Toxicol* 51(117-44), 2010.
99. Reynolds A, Hanani R, Hibbs D, Damont A, Da Pozzo E, Selleri S, Dolle F, Martini C, Kassiou M. Pyrazolo[1,5-a]pyrimidine acetamides: 4-Phenyl alkyl ether derivatives as potent ligands for the 18 kDa translocator protein (TSPO). *Bioorg Med Chem Lett* 20(19): 5799-802, 2010.
100. Frankowski KJ, Ghosh P, Setola V, Tran TB, **Roth BL**, Aube J. N-Alkyl-octahydroisoquinolin-1-one-8-carboxamides: a Novel Class of Selective, Nonbasic, Nitrogen-Containing kappa-Opioid Receptor Ligands. *ACS Med Chem Lett* 1(5): 189-193, 2010.
101. Gianni D, Taulet N, Zhang H, DerMardirossian C, Kister J, Martinez L, Roush WR, Brown SJ, Bokoch GM, Rosen H. A novel and specific NADPH oxidase-1 (Nox1) small-molecule inhibitor blocks the formation of functional invadopodia in human colon cancer cells. *ACS Chem Biol* 5(10): 981-93, 2010.
102. Kwon YJ, Saubern S, Macdonald JM, Huang XP, Setola V, **Roth BL**. N-tetrahydrothiochromenoisoxazole-1-carboxamides as selective antagonists of cloned human 5-HT<sub>2B</sub>. *Bioorg Med Chem Lett* 20(18): 5488-90, 2010.
103. Sakurada T, Gill MB, Frausto S, Copits B, Noguchi K, Shimamoto K, Swanson GT, Sakai R. Novel N-methylated 8-oxoisoguanines from Pacific sponges with diverse neuroactivities. *J Med Chem* 53(16): 6089-99, 2010.
104. Green BR, White KL, McDougle DR, Zhang L, Klein B, Scholl EA, Pruess TH, White HS, Bulaj G. Introduction of lipidization-cationization motifs affords systemically bioavailable

- neuropeptide Y and neurotensin analogs with anticonvulsant activities. *J Pept Sci* 16(9): 486-95, 2010.
105. Hedlund PB, Leopoldo M, Caccia S, Sarkisyan G, Fracasso C, Martelli G, Lacivita E, Berardi F, Perrone R. LP-211 is a brain penetrant selective agonist for the serotonin 5-HT(7) receptor. *Neurosci Lett* 481(1): 12-6, 2010.
106. Milak MS, Severance AJ, Prabhakaran J, Kumar JS, Majo VJ, Ogden RT, Mann JJ, Parsey RV. In vivo serotonin-sensitive binding of [<sup>11</sup>C]CUMI-101: a serotonin 1A receptor agonist positron emission tomography radiotracer. *J Cereb Blood Flow Metab* 31(1): 243-9, 2010.
107. Kozikowski AP, Cho SJ, Jensen NH, Allen JA, Svennebring AM, **Roth BL**. HTS and rational drug design to generate a class of 5-HT(2C)-selective ligands for possible use in schizophrenia. *ChemMedChem* 5(8): 1221-5, 2010.
108. Dong S, Allen JA, Farrell M, **Roth BL**. A chemical-genetic approach for precise spatio-temporal control of cellular signaling. *Mol Biosyst* 6(8): 1376-80, 2010.
109. Peddibhotla S, Shi R, Khan P, Smith LH, Mangravita-Novo A, Vicchiarelli M, Su Y, Okolotowicz KJ, Cashman JR, Reed JC, Roth GP. Inhibition of protein kinase C-driven nuclear factor-kappaB activation: synthesis, structure-activity relationship, and pharmacological profiling of pathway specific benzimidazole probe molecules. *J Med Chem* 53(12): 4793-7, 2010.
110. Hooker JM, Kim SW, Reibel AT, Alexoff D, Xu Y, Shea C. Evaluation of [(11)C]metergoline as a PET radiotracer for 5HTR in nonhuman primates. *Bioorg Med Chem* 18(22): 7739-45, 2011.
111. Scorza MC, Castane A, Bortolozzi A, Artigas F. Clozapine does not require 5-HT1A receptors to block the locomotor hyperactivity induced by MK-801 Clz and MK-801 in KO1A mice. *Neuropharmacology* 59(1-2): 112-20, 2010.
112. McCormick PN, Kapur S, Graff-Guerrero A, Raymond R, Nobrega JN, Wilson AA. The antipsychotics olanzapine, risperidone, clozapine, and haloperidol are D2-selective ex vivo but not in vitro. *Neuropsychopharmacology* 35(8): 1826-35, 2010.

113. Zhang P, Zou MF, Rodriguez AL, Conn PJ, Newman AH. Structure-activity relationships in a novel series of 7-substituted-aryl quinolines and 5-substituted-aryl benzothiazoles at the metabotropic glutamate receptor subtype 5. *Bioorg Med Chem* 18(9): 3026-35, 2010.
114. Davies MA, Conley Y, Puskar K. Incorporating evidence from pharmacologic and pharmacogenetic studies of atypical antipsychotic drugs into advanced psychiatric nursing practice. *Perspect Psychiatr Care* 46(2): 98-107, 2010.
115. Nemeth CL, Paine TA, Rittiner JE, Beguin C, Carroll FI, **Roth BL**, Cohen BM, Carlezon WA, Jr. Role of kappa-opioid receptors in the effects of salvinorin A and ketamine on attention in rats. *Psychopharmacology (Berl)* 210(2): 263-74, 2010.
116. Ogawa M, Nishiyama S, Tsukada H, Hatano K, Fuchigami T, Yamaguchi H, Matsushima Y, Ito K, Magata Y. Synthesis and evaluation of new imaging agent for central nicotinic acetylcholine receptor alpha7 subtype. *Nucl Med Biol* 37(3): 347-55, 2010.
117. Dong S, Rogan SC, **Roth BL**. Directed molecular evolution of DREADDs: a generic approach to creating next-generation RASSLs. *Nat Protoc* 5(3): 561-73, 2010.
118. Gao M, Wang M, Miller KD, Hutchins GD, Zheng QH. Synthesis and in vitro biological evaluation of carbon-11-labeled quinoline derivatives as new candidate PET radioligands for cannabinoid CB2 receptor imaging. *Bioorg Med Chem* 18(6): 2099-106, 2010.
119. Canal CE, Olaghere da Silva UB, Gresch PJ, Watt EE, Sanders-Bush E, Airey DC. The serotonin 2C receptor potently modulates the head-twitch response in mice induced by a phenethylamine hallucinogen. *Psychopharmacology (Berl)* 209(2): 163-74, 2010.
120. Bortolozzi A, Masana M, Diaz-Mataix L, Cortes R, Scorza MC, Gingrich JA, Toth M, Artigas F. Dopamine release induced by atypical antipsychotics in prefrontal cortex requires 5-HT(1A) receptors but not 5-HT(2A) receptors. *Int J Neuropsychopharmacol* 13(10): 1299-314, 2010.
121. Farrell M, **Roth BL**. Allosteric antipsychotics: m4 muscarinic potentiators as novel treatments for schizophrenia. *Neuropsychopharmacology* 35(4): 851-2, 2010.
122. Strachan RT, Allen JA, Sheffler DJ, **Roth BL**. p90 Ribosomal S6 kinase 2, a novel GPCR kinase, is required for growth factor-mediated attenuation of GPCR signaling. *Biochemistry* 49(12): 2657-71, 2010.
123. Ray TS. Psychedelics and the human receptorome. *PLoS One* 5(2): e9019, 2010.

124. Palner M, McCormick P, Gillings N, Begtrup M, Wilson AA, Knudsen GM. Radiosynthesis and ex vivo evaluation of (R)-(-)-2-chloro-N-[1-11C-propyl]n-propylnorapomorphine. *Nucl Med Biol* 37(1): 35-40, 2010.
125. Eng KT, Berdis AJ. A novel non-natural nucleoside that influences P-glycoprotein activity and mediates drug resistance. *Biochemistry* 49(8): 1640-8, 2010.
126. Pei Y, Dong S, **Roth BL**. Generation of designer receptors exclusively activated by designer drugs (DREADDs) using directed molecular evolution. *Curr Protoc Neurosci Chapter 4*(Unit 4 33, 2010.
127. Setola V. Receptorome screening: a powerful, facile approach to better understand and engineer drugs. *Drug News Perspect* 22(8): 459-66, 2009.
128. Legendre O, Pecic S, Chaudhary S, Zimmerman SM, Fantegrossi WE, Harding WW. Synthetic studies and pharmacological evaluations on the MDMA ('Ecstasy') antagonist nantenine. *Bioorg Med Chem Lett* 20(2): 628-31, 2009.
129. Banister SD, Moussa IA, Jordan MJ, Coster MJ, Kassiou M. Oxo-bridged isomers of aza-trishomocubane sigma (sigma) receptor ligands: Synthesis, in vitro binding, and molecular modeling. *Bioorg Med Chem Lett* 20(1): 145-8, 2009.
130. Mavel S, Meheux N, Guilloteau D, Emond P. Synthesis and in vitro evaluation of fluorinated diphenyloxide derivatives and sulfur analogs as serotonin transporter ligands. *Bioorg Med Chem* 18(1): 236-41, 2009.
131. Gillman PK. Triptans, serotonin agonists, and serotonin syndrome (serotonin toxicity): a review. *Headache* 50(2): 264-72, 2009.
132. Jones KA, Srivastava DP, Allen JA, Strachan RT, **Roth BL**, Penzes P. Rapid modulation of spine morphology by the 5-HT<sub>2A</sub> serotonin receptor through kalirin-7 signaling. *Proc Natl Acad Sci U S A* 106(46): 19575-80, 2009.
133. Keiser MJ, Setola V, Irwin JJ, Laggner C, Abbas AI, Hufeisen SJ, Jensen NH, Kuijter MB, Matos RC, Tran TB, Whaley R, Glennon RA, Hert J, Thomas KL, Edwards DD, Shoichet BK, **Roth BL**. Predicting new molecular targets for known drugs. *Nature* 462(7270): 175-81, 2009. (BLR and BKS = Co-corresponding authors).

134. Brothers SP, Saldanha SA, Spicer TP, Cameron M, Mercer BA, Chase P, McDonald P, Wahlestedt C, Hodder PS. Selective and brain penetrant neuropeptide  $\gamma$   $\gamma$ 2 receptor antagonists discovered by whole-cell high-throughput screening. *Mol Pharmacol* 77(1): 46-57, 2009.
135. Donohue SR, Varnas K, Jia Z, Gulyas B, Pike VW, Halldin C. Synthesis and in vitro autoradiographic evaluation of a novel high-affinity radioiodinated ligand for imaging brain cannabinoid subtype-1 receptors. *Bioorg Med Chem Lett* 19(21): 6209-12, 2009.
136. Mestres J, Gregori-Puigjane E. Conciliating binding efficiency and polypharmacology. *Trends Pharmacol Sci* 30(9): 470-4, 2009.
137. Allen JA, Yu JZ, Dave RH, Bhatnagar A, **Roth BL**, Rasenick MM. Caveolin-1 and lipid microdomains regulate Gs trafficking and attenuate Gs/adenylyl cyclase signaling. *Mol Pharmacol* 76(5): 1082-93, 2009.
138. Peng Y, Zhang Q, Arora S, Keenan SM, Kortagere S, Wannemacher KM, Howells RD, Welsh WJ. Novel delta opioid receptor agonists exhibit differential stimulation of signaling pathways. *Bioorg Med Chem* 17(17): 6442-50, 2009.
139. Huang XP, Setola V, Yadav PN, Allen JA, Rogan SC, Hanson BJ, Revankar C, Robers M, Doucette C, **Roth BL**. Parallel functional activity profiling reveals valvulopathogens are potent 5-hydroxytryptamine(2B) receptor agonists: implications for drug safety assessment. *Mol Pharmacol* 76(4): 710-22, 2009.
140. Yan F, Bikbulatov RV, Mocanu V, Dicheva N, Parker CE, Wetsel WC, Mosier PD, Westkaemper RB, Allen JA, Zjawiony JK, **Roth BL**. Structure-based design, synthesis, and biochemical and pharmacological characterization of novel salvinorin A analogues as active state probes of the kappa-opioid receptor. *Biochemistry* 48(29): 6898-908, 2009.
141. Hanson BJ, Wetter J, Bercher MR, Kopp L, Fuerstenau-Sharp M, Vedvik KL, Zielinski T, Doucette C, Whitney PJ, Revankar C. A homogeneous fluorescent live-cell assay for measuring 7-transmembrane receptor activity and agonist functional selectivity through beta-arrestin recruitment. *J Biomol Screen* 14(7): 798-810, 2009.
142. McIntosh JM, Absalom N, Chebib M, Elgoyhen AB, Vincler M. Alpha9 nicotinic acetylcholine receptors and the treatment of pain. *Biochem Pharmacol* 78(7): 693-702, 2009.

143. Kimple AJ, Soundararajan M, Hutsell SQ, Roos AK, Urban DJ, Setola V, Temple BR, **Roth BL**, Knapp S, Willard FS, Siderovski DP. Structural determinants of G-protein alpha subunit selectivity by regulator of G-protein signaling 2 (RGS2). *J Biol Chem* 284(29): 19402-11, 2009.
144. Kulkarni SS, Zou MF, Cao J, Deschamps JR, Rodriguez AL, Conn PJ, Newman AH. Structure-activity relationships comparing N-(6-methylpyridin-yl)-substituted aryl amides to 2-methyl-6-(substituted-arylethynyl)pyridines or 2-methyl-4-(substituted-arylethynyl)thiazoles as novel metabotropic glutamate receptor subtype 5 antagonists. *J Med Chem* 52(11): 3563-75, 2009.
145. Herth MM, Piel M, Debus F, Schmitt U, Luddens H, Rosch F. Preliminary in vivo and ex vivo evaluation of the 5-HT<sub>2A</sub> imaging probe [(18F)MH.MZ. *Nucl Med Biol* 36(4): 447-54, 2009.
146. Fantegrossi WE, Bauzo RM, Manvich DM, Morales JC, Votaw JR, Goodman MM, Howell LL. Role of dopamine transporters in the behavioral effects of 3,4-methylenedioxymethamphetamine (MDMA) in nonhuman primates. *Psychopharmacology (Berl)* 205(2): 337-47, 2009.
147. Booth RG, Fang L, Huang Y, Wilczynski A, Sivendran S. (1R, 3S)-(-)-trans-PAT: a novel full-efficacy serotonin 5-HT<sub>2C</sub> receptor agonist with 5-HT<sub>2A</sub> and 5-HT<sub>2B</sub> receptor inverse agonist/antagonist activity. *Eur J Pharmacol* 615(1-3): 1-9, 2009.
148. Seeman P, Guan HC, Hirbec H. Dopamine D<sub>2</sub>High receptors stimulated by phencyclidines, lysergic acid diethylamide, salvinorin A, and modafinil. *Synapse* 63(8): 698-704, 2009.
149. Abbas AI, Hedlund PB, Huang XP, Tran TB, Meltzer HY, **Roth BL**. Amisulpride is a potent 5-HT<sub>7</sub> antagonist: relevance for antidepressant actions in vivo. *Psychopharmacology (Berl)* 205(1): 119-28, 2009.
150. Zolkowska D, Jain R, Rothman RB, Partilla JS, **Roth BL**, Setola V, Prisinzano TE, Baumann MH. Evidence for the involvement of dopamine transporters in behavioral stimulant effects of modafinil. *J Pharmacol Exp Ther* 329(2): 738-46, 2009.
151. Tang H, Wang XS, Huang XP, **Roth BL**, Butler KV, Kozikowski AP, Jung M, Tropsha A. Novel inhibitors of human histone deacetylase (HDAC) identified by QSAR modeling of known

- inhibitors, virtual screening, and experimental validation. *J Chem Inf Model* 49(2): 461-76, 2009.
152. Sikazwe DM, Nkansah NT, Altundas R, Zhu XY, **Roth BL**, Setola V, Ablordeppey SY. Synthesis and evaluation of ligands for D2-like receptors: the role of common pharmacophoric groups. *Bioorg Med Chem* 17(4): 1716-23, 2009.
153. Pei Y, Rogan SC, Yan F, **Roth BL**. Engineered GPCRs as tools to modulate signal transduction. *Physiology (Bethesda)* 23(313-21), 2008.
154. Przybyla JA, Cueva JP, Chemel BR, Hsu KJ, Riese DJ, 2nd, McCorvy JD, Chester JA, Nichols DE, Watts VJ. Comparison of the enantiomers of (+/-)-doxanthrine, a high efficacy full dopamine D(1) receptor agonist, and a reversal of enantioselectivity at D(1) versus alpha(2C) adrenergic receptors. *Eur Neuropsychopharmacol* 19(2): 138-46, 2009.
155. Buccafusco JJ, Beach JW, Terry AV, Jr. Desensitization of nicotinic acetylcholine receptors as a strategy for drug development. *J Pharmacol Exp Ther* 328(2): 364-70, 2009.
156. Bug WJ, Ascoli GA, Grethe JS, Gupta A, Fennema-Notestine C, Laird AR, Larson SD, Rubin D, Shepherd GM, Turner JA, Martone ME. The NIFSTD and BIRNLex vocabularies: building comprehensive ontologies for neuroscience. *Neuroinformatics* 6(3): 175-94, 2008.
157. Harmar AJ, Hills RA, Rosser EM, Jones M, Buneman OP, Dunbar DR, Greenhill SD, Hale VA, Sharman JL, Bonner TI, Catterall WA, Davenport AP, Delagrang P, Dollery CT, Foord SM, Gutman GA, Laudet V, Neubig RR, Ohlstein EH, Olsen RW, Peters J, Pin JP, Ruffolo RR, Searls DB, Wright MW, Spedding M. IUPHAR-DB: the IUPHAR database of G protein-coupled receptors and ion channels. *Nucleic Acids Res* 37(Database issue): D680-5, 2009.
158. Mavel S, Vercouillie J, Garreau L, Raguza T, Ravna AW, Chalon S, Guilloteau D, Emond P. Docking study, synthesis, and in vitro evaluation of fluoro-MADAM derivatives as SERT ligands for PET imaging. *Bioorg Med Chem* 16(19): 9050-5, 2008.
159. Lazarova N, Zoghbi SS, Hong J, Seneca N, Tuan E, Gladding RL, Liow JS, Taku A, Innis RB, Pike VW. Synthesis and evaluation of [N-methyl-11C]N-desmethyl-loperamide as a new and improved PET radiotracer for imaging P-gp function. *J Med Chem* 51(19): 6034-43, 2008.

160. Richtand NM, Welge JA, Logue AD, Keck PE, Jr., Strakowski SM, McNamara RK. Role of serotonin and dopamine receptor binding in antipsychotic efficacy. *Prog Brain Res* 172(155-75, 2008.
161. Donohue SR, Pike VW, Finnema SJ, Truong P, Andersson J, Gulyas B, Halldin C. Discovery and labeling of high-affinity 3,4-diarylpyrazolines as candidate radioligands for in vivo imaging of cannabinoid subtype-1 (CB1) receptors. *J Med Chem* 51(18): 5608-16, 2008.
162. Basura GJ, Abbas AI, O'Donohue H, Lauder JM, **Roth BL**, Walker PD, Manis PB. Ontogeny of serotonin and serotonin2A receptors in rat auditory cortex. *Hear Res* 244(1-2): 45-50, 2008.
163. Liang Y, Shaw AM, Boules M, Briody S, Robinson J, Oliveros A, Blazar E, Williams K, Zhang Y, Carlier PR, Richelson E. Antidepressant-like pharmacological profile of a novel triple reuptake inhibitor, (1S,2S)-3-(methylamino)-2-(naphthalen-2-yl)-1-phenylpropan-1-ol (PRC200-SS). *J Pharmacol Exp Ther* 327(2): 573-83, 2008.
164. Stephenson KA, van Oosten EM, Wilson AA, Meyer JH, Houle S, Vasdev N. Synthesis and preliminary evaluation of [(18)F]-fluoro-(2S)-Exaprolol for imaging cerebral beta-adrenergic receptors with PET. *Neurochem Int* 53(5): 173-9, 2008.
165. Jensen NH, **Roth BL**. Massively parallel screening of the receptorome. *Comb Chem High Throughput Screen* 11(6): 420-6, 2008.
166. Allen JA, Yadav PN, **Roth BL**. Insights into the regulation of 5-HT2A serotonin receptors by scaffolding proteins and kinases. *Neuropharmacology* 55(6): 961-8, 2008.
167. Marona-Lewicka D, Chemel BR, Nichols DE. Dopamine D4 receptor involvement in the discriminative stimulus effects in rats of LSD, but not the phenethylamine hallucinogen DOI. *Psychopharmacology (Berl)* 203(2): 265-77, 2009.
168. Dutta AK, Ghosh B, Biswas S, Reith ME. D-161, a novel pyran-based triple monoamine transporter blocker: behavioral pharmacological evidence for antidepressant-like action. *Eur J Pharmacol* 589(1-3): 73-9, 2008.
169. Gunosewoyo H, Guo JL, Bennett MR, Coster MJ, Kassiou M. Cubyl amides: novel P2X7 receptor antagonists. *Bioorg Med Chem Lett* 18(13): 3720-3, 2008.

170. Bojarski AJ, Kuran B, Kossakowski J, Koziol A, Jagiello-Wojtowicz E, Chodkowska A. Synthesis and serotonin receptor activity of the arylpiperazine alkyl/propoxy derivatives of new azatricycloundecanes. *Eur J Med Chem* 44(1): 152-64, 2009.
171. Wang XS, Tang H, Golbraikh A, Tropsha A. Combinatorial QSAR modeling of specificity and subtype selectivity of ligands binding to serotonin receptors 5HT1E and 5HT1F. *J Chem Inf Model* 48(5): 997-1013, 2008.
172. **Roth BL**. High-dose olanzapine for treatment-resistant schizophrenia. *J Clin Psychiatry* 69(2): 176-7, 2008.
173. Reid AE, Ding YS, Eckelman WC, Logan J, Alexoff D, Shea C, Xu Y, Fowler JS. Comparison of the pharmacokinetics of different analogs of 11C-labeled TZTP for imaging muscarinic M2 receptors with PET. *Nucl Med Biol* 35(3): 287-98, 2008.
174. Giorgioni G, Ruggieri S, Claudi F, Di Stefano A, Ljung E, Carlsson T. Synthesis and pharmacological evaluation of 4-phenoxy-1,2,3,4-tetrahydroisoquinolines and 4,5,6,6a-tetrahydrochromeno[2,3,4-de]isoquinolines. *Med Chem* 4(1): 1-10, 2008.
175. Burns AC, Forsyth CJ. Intramolecular Diels-Alder/Tsuji allylation assembly of the functionalized trans-decalin of salvinorin A. *Org Lett* 10(1): 97-100, 2008.
176. Jensen NH, Rodriguiz RM, Caron MG, Wetsel WC, Rothman RB, **Roth BL**. N-desalkylquetiapine, a potent norepinephrine reuptake inhibitor and partial 5-HT1A agonist, as a putative mediator of quetiapine's antidepressant activity. *Neuropsychopharmacology* 33(10): 2303-12, 2008.
177. Beguin C, Potter DN, Dinieri JA, Munro TA, Richards MR, Paine TA, Berry L, Zhao Z, **Roth BL**, Xu W, Liu-Chen LY, Carlezon WA, Jr., Cohen BM. N-methylacetamide analog of salvinorin A: a highly potent and selective kappa-opioid receptor agonist with oral efficacy. *J Pharmacol Exp Ther* 324(1): 188-95, 2008.
178. Buhler AV, Proudfit HK, Gebhart GF. Neurotensin-produced antinociception in the rostral ventromedial medulla is partially mediated by spinal cord norepinephrine. *Pain* 135(3): 280-90, 2008.