

# [<sup>18</sup>F]FBR For Injection: Quality Control Record

PET Radiopharmaceutical Sciences Section,  
Molecular Imaging Branch,  
National Institute of Mental Health,  
National Institutes of Health,  
Bldg. 10, Rm. B3 C346,  
Bethesda, MD 20892

Date of review: DRAFT 8/17/2009

Approved by: \_\_\_\_\_ Initials: \_\_\_\_\_ Date: \_\_\_\_\_

Victor W Pike, Ph.D., Chief, PET, Radiopharmaceutical Sciences, NIMH

Batch # FBR- \_\_\_\_\_ Date: \_\_\_\_\_

## Quality Control Instrument and Materials Verification

Verification of Ionization Chamber	Test Data	Acceptance Criteria	Acceptance Criteria Met
<sup>137</sup> Cs source reference ID 970-31-4 204 μCi on 3/1/2003 Expected radioactivity _____ μCi	<sup>137</sup> Cs measured _____ μCi	Does measured radioactivity fall within acceptable range (± 5%) for both <sup>137</sup> Cs and <sup>57</sup> Co?	<b>Y / N?</b>
<sup>57</sup> Co source reference ID _____ _____ μCi on _____ Expected radioactivity _____ μCi	<sup>57</sup> Co measured _____ μCi		

## Verification of Analytical HPLC system

Verification of Analytical HPLC system	Test Data	Acceptance Criteria	Acceptance Criteria Met
Curve Date _____ Slope (area/ μmol) _____ Volume Injected _____ μL Expected Peak area _____	Measured Peak Area _____  t <sub>R</sub> _____ min	Peak area within ± 10% of expected peak area?	<b>Y / N?</b>

## Post-Production Measurements and Release Tests

Volume	Data
Subtract the weight of the empty dose vial from the full dose vial. Assume the density of the final product is approximately 1 g/ mL to determine volume. Weight of full dose vial _____ Weight of empty dose vial _____	Volume of [ <sup>18</sup> F]FBR for Injection _____ mL
Yield	
Measure the activity of the whole dose after removal of the QC sample. Calculate the percent yield without correcting for decay.	Measured Activity _____ mCi @ _____  Non-decay corrected _____ percent yield

## Release Tests

Test	Test Data	Acceptance Criteria	Criteria Met
pH			
Dispense one to two drops of [ <sup>18</sup> F]FBR for Injection on pH paper.	Measured pH _____	pH measured between 4.5 to 8.0?	<b>Y / N?</b>
Appearance			
Visually inspect the [ <sup>18</sup> F]FBR for Injection in the dose vial.	Appearance:	Colorless solution free of particulate matter?	<b>Y / N?</b>

# [<sup>18</sup>F]FBR For Injection: Quality Control Record

PET Radiopharmaceutical Sciences Section,  
Molecular Imaging Branch,  
National Institute of Mental Health,  
National Institutes of Health,  
Bldg. 10, Rm. B3 C346,  
Bethesda, MD 20892

Date of review: DRAFT 8/17/2009

<b>Membrane Filter Integrity</b>			
Attach the syringe filter and needle used for sterile filtration of [ <sup>18</sup> F]FBR for Injection to the compressed air source. Submerge needle tip in water. Pressurize to 45 psi.	Observation: Bubbles / No Bubbles	No bubbles observed from the submerged needle tip at 45 psi?	Y / N?
<b>Radionuclidic identity</b>			
Determine the half life experimentally from two time points separated by at least three minutes. Compare to the theoretical	t <sub>1</sub> _____ mCi @ _____ t <sub>2</sub> _____ mCi @ _____ Calc half life _____ min	Experimental half-life 110 ± 10 minutes?	Y / N?
<b>Residual Solvent</b>			
Determine the residual ethanol and acetonitrile by GC	Ethanol in sample _____ ng/ µL Acetonitrile in sample _____ ng/ µL Max. Injectable Volume _____ mL	Ethanol in sample ≤ 1 x 10 <sup>5</sup> ng/ µL? Acetonitrile in injectable volume ≤ 4 x 10 <sup>2</sup> ?	Y / N? Y / N?
<b>Tests Based on HPLC Analysis<sup>1</sup></b>			
<b>Radioconcentration</b>			
Measure the radioactivity of a 100 µL aliquot of [ <sup>18</sup> F]FBR for Injection _____ µCi @ _____ Measure the empty syringe after injection. _____ µCi @ _____	Net activity injected _____ µCi @ _____ Calculated conc. = _____ mCi/mL	Greater than 0.5 mCi/ mL	Y / N?
<b>Radiochemical Identity</b>			
Compare the retention time of [ <sup>18</sup> F]FBR for Injection to the FBR Standard retention time.	STD t <sub>R</sub> _____ min Product t <sub>R</sub> _____ min	Difference in retention time is less than 1.0 minutes?	Y / N?
<b>Radiochemical Purity</b>			
Integrate the peaks in the gamma HPLC trace. Determine the percent area represented by the product peak.	% Area Bioscan _____ %	% area radioactive trace ≥ 95%?	Y / N?
<b>Chemical Purity</b>			
Calculate the concentration of carrier and FBR equivalent impurity in the [ <sup>18</sup> F]FBR for Injection from the peak area at 320 nm.	FBR carrier _____ µg/mL Impurity _____ µg/mL	Maximum volume contains no more than 10 µg of carrier or 1 µg FBR equivalent impurity.	Y / N?
Calculate the maximum injection volume as the lesser volume by maximum allowable carrier or maximum allowable FRB equivalent impurity.	Max Vol by carrier _____ mL Max Vol by imp. _____ mL		
<b>Specific Radioactivity</b>			
Calculate the specific activity of the [ <sup>18</sup> F]FBR for Injection in units of mCi/ µmol	Calculated Specific Activity _____ mCi/ µmol	Specific activity > 500 mCi/ µmol?	Y / N?

Attach a copy of the Calculations worksheet to this form

Chemist Name	Signature	Date

<sup>1</sup> A calculations worksheet is typically used to perform the calculations required. A copy of the worksheet may be found in Appendix C of the document "[<sup>18</sup>F]-FBR for Injection: Standard Operating Procedures"  
Document 4:[<sup>18</sup>F]FBR for Injection: Quality Control Record